



EXCOR[®] Pediatric VAD

Ventricular Assist Device

***with Stationary Driving Unit Ikus Rev. 2.1
for Pediatric Use***

Physician's Manual Rev. 2

For products in USA:

Humanitarian Device. Authorized by Federal law for use in the treatment of pediatric patients with severe isolated left ventricular or biventricular dysfunction who are candidates for cardiac transplant and require circulatory support. The effectiveness of this device for this use has not been demonstrated.

Berlin Heart, Inc.
200 Valleywood, Suite A500
The Woodlands, TX 77380
USA

Tel. 281.863.9700
Fax 281.863.9701
Hotline: 866.249.0128

info@berlinheart.com
www.berlinheart.com

This page was left blank intentionally

This physician's manual (PM) corresponds to the following product versions:

- *Ikus* software: from V 3.41 forward
- Laptop software: from V 3.41 V forward
- Laptop from CF30 forward

Manufacturer

Berlin Heart GmbH
Wiesenweg 10
12247 Berlin
Germany

Distributor

Berlin Heart, Inc.
200 Valleywood, Suite A500
The Woodlands, TX 77380
USA

Office Phone: 1.281.863.9700

Office Fax: 1.281.863.9701

Humanitarian Device. Authorized by Federal law for use in the treatment of patients with severe isolated left ventricular or biventricular dysfunction who are candidates for cardiac transplant and require circulatory support. The effectiveness of this device for this use has not been demonstrated.

Imprint

© 2012-03

Berlin Heart GmbH
Wiesenweg 10
12247 Berlin

All rights reserved.


This PM may not be reproduced, transmitted, transferred onto an electronic medium, stored in an electronic data processing system nor translated into another language, either in parts or as a whole, without the express prior written permission of *Berlin Heart GmbH*. Excluded are the sample copies contained in this instruction for use.

This instruction for use is intended for information purposes only. The contents of the instruction for use may be supplemented, modified or updated at any time without prior notice.

Any previous versions of the PM become obsolete upon release of this version.

All translations of this PM have been prepared and checked to the best of the translator's knowledge. The German version of this instruction for use shall prevail for all matters of interpretation within a legal context.

Licenses, trademarks and trade names

EXCOR® and  **Berlin Heart'** are trademarks of Berlin Heart GmbH. They are protected by the relevant legislation in Germany and other countries.

All other trademarks or brands mentioned in this instruction for use are acknowledged as being the property of the respective owners. These rights are specifically recognized and respected whenever such trademarks or brands are mentioned in this instruction for use.

Text and layout

Berlin Heart GmbH

This page was left blank intentionally

Table of contents

1	Important safety information	11
1.1	Warnings	11
1.1.1	Storage and durability	11
1.1.2	Device configurations.....	11
1.1.3	Procedural techniques - Ikus.....	12
1.1.4	Packaging and sterilization	12
1.1.5	Procedural techniques - pumps, cannulae, accessories	13
1.1.6	System.....	15
1.1.7	Procedures to minimize risk of thrombosis	15
1.1.8	Cleaning the components	16
1.1.9	Errors and corrective measures	16
1.1.10	Replacing the blood pump(s)	17
1.1.11	Driving blood pump(s) with the manual pump	17
1.1.12	Ambient conditions	18
1.1.13	Interaction with other procedures and therapies	18
1.2	Precautions.....	19
1.2.1	VAD placement technique.....	19
1.2.2	Ambient conditions	19
1.2.3	Caution while using the Ikus.....	19
1.3	Obligations of the operator	20
2	General Information	21
2.1	Device description.....	21
2.2	Indications for use	21
2.3	IDE Clinical Study Summary.....	21
2.4	Intended operation environment.....	21
2.5	Contraindications.....	22
2.6	Storage and durability	22
3	IDE Clinical Study Summary	23
3.1	Indications for use	23
3.2	Contraindications.....	23
3.3	Alternative Practices or Procedures.....	23
3.4	Marketing History	23
3.5	Potential Adverse Effects	23
3.6	Summary of Clinical Studies.....	26
3.6.1	IDE Clinical Study Summary	26
3.6.2	Study Cohorts	26
3.6.3	Inclusion/Exclusion Criteria	27
3.6.4	Study Enrollment	28
3.6.5	Subject Demographics	29

3.6.6	Results	32
3.6.7	Conclusion	42
4	Description: blood pump, cannulae and accessories	45
4.1	EXCOR blood pumps	46
4.2	EXCOR cannulae.....	46
4.3	EXCOR accessories	47
5	Implantation: Preparations in the operating room	49
5.1	Preparing the components and materials required	49
5.2	Checking and adjusting the settings of the cable tie gun	49
5.3	Unpacking the sterile components.....	50
5.4	Moving the membrane to the end-of-diastole position	51
5.5	De-airing the blood pump.....	51
5.5.1	Inserting the de-airing needle.....	51
5.5.2	Rinsing and filling the blood pump.....	52
6	Implantation - surgical procedure	53
6.1	Cannula exit sites.....	54
6.2	Use of the cannula tunneling tip.....	54
6.3	Cannulae and connector set	55
6.4	Access	56
6.5	LV apex cannula.....	56
6.5.1	Anastomosis of inflow cannula with LV apex	57
6.5.2	Creating a transcutaneous tunnel for the LV apex cannula	58
6.6	Atrial cannula(e)	59
6.6.1	Creating a transcutaneous tunnel for atrial cannula(e)	59
6.6.2	Anastomosis of atrial cannulae	60
6.7	Arterial cannula(e)	60
6.7.1	Creating a transcutaneous tunnel for arterial cannula	61
6.7.2	Anastomosis of the arterial cannula	61
6.8	Shortening the cannulae if necessary	62
6.9	Connecting the blood pumps to the cannulae.....	63
6.10	Intraoperative drive management.....	64
6.10.1	Connecting the blood pump(s) to the Ikus	64
6.10.2	De-airing the blood pumps in single-step mode.....	65
6.10.3	Starting the blood pump (changing to standard view)	65
6.10.4	Checking the parameters when the pump is started and adjusting them	66
6.10.5	Switching from CPB support to VAD support	68
6.10.6	Possible complications.....	68
6.11	Removing the de-airing needle	69

6.12	Securing the connections	70
6.13	Postoperative drive management.....	71
6.13.1	After transfer to the ward	71
6.13.2	Follow-up treatment	71
7	Implantation - anesthesia	73
8	Wound care and treatment	75
8.1	Removing the old dressings	76
8.2	Cleaning the blood pump	76
8.3	Cleaning of the wound	77
8.4	The new dressing	77
8.4.1	Preparing a new dressing	77
8.4.2	Applying a new dressing.....	78
8.5	Regular checks of blood pump(s) and cannulae	79
8.5.1	Visual inspection: pump filling and ejection	80
8.5.2	Visual inspection: deposits	81
8.5.3	Checks using the monitor program	81
8.5.4	Replacing the blood pump due to growth of the patient	82
9	Anticoagulation therapy	83
9.1	Before Implantation of the EXCOR.....	83
9.1.1	General considerations	83
9.1.2	Pre implantation	83
9.2	During Implantation - Cardiopulmonary Bypass	83
9.2.1	Cardiopulmonary Bypass (CPB).....	83
9.2.2	Post CPB	83
9.3	Postoperative anticoagulation therapy.....	83
9.3.1	General Considerations.....	83
9.3.2	Starting anticoagulation therapy.....	83
9.3.3	Unfractionated heparin therapy (i.v.) Patient < 12 months	84
9.3.4	Unfractionated heparin therapy (i.v.) Patient ≥ 12 months	84
9.3.5	Thrombelastography (TEG®).....	84
9.4	Low Molecular Weight Heparin	85
9.5	Oral Anticoagulation Therapy (only for patients ≥ 12 months of age who are taking a full oral diet)	85
9.6	Monitoring of Blood Count and Anticoagulation Status	85
9.7	Postoperative platelet inhibition therapy	86
9.7.1	Start of therapy	86
9.8	Adjunctive Medication	86
9.9	Anticoagulation Therapy	87
9.9.1	Therapeutic Heparin administration and adjustment	87
9.9.2	Oral Anticoagulation Therapy	88

10	Troubleshooting and correcting faults	89
10.1	Replacing the blood pump(s)	89
10.1.1	Preparing a replacement blood pump	90
10.1.2	Replacing the right blood pump (RVAD/ BVAD)	90
10.1.3	Replacing the left blood pump (LVAD/ BVAD)	91
10.2	Driving blood pump(s) with the manual pump	92
11	Weaning and Explantation for BTR and BTT	95
11.1	Weaning Procedure	95
11.1.1	Introduction	95
11.1.2	Indication	95
11.1.3	Eligibility Criteria	95
11.1.4	Weaning Protocol.....	95
11.1.5	10 ml pump.....	96
11.1.6	25/ 30 ml pump	100
11.1.7	50/ 60 ml pump	104
11.1.8	Explantation Criteria	107
11.2	Explantation for BTR.....	108
11.2.1	Explantation with univentricular support.....	108
11.2.2	Explantation after biventricular support	108
11.3	Explantation for BTT.....	109
12	Appendix	111
12.1	Overview: Product range and possible combinations	111
12.1.1	Blood pumps.....	111
12.1.2	Overview: Relationship: body weight – pump size.....	111
12.1.3	LV apex cannulae	112
12.1.4	Atrial cannulae	112
12.1.5	Arterial cannulae	112
12.1.6	Overview: Which cannulae should be used for which pump?	113
12.1.7	System accessories	113
12.1.8	Driving unit	113
12.1.9	Special components	113
12.1.10	Pump-cannula combinations	114
12.1.11	Blood pump combinations in biventricular mode	114
12.1.12	Relative systolic duration.....	114
12.2	Sample copy: EXCOR Implantation log.....	115
12.3	Sample copy: EXCOR pump log	118
12.3.1	Explanations on the pump log	118
12.3.2	Sample copy: EXCOR pump log.....	119
12.4	Sample copy: EXCOR parameter log	120

Dear readers,

This Physician's Manual (PM) is intended for all medical personnel involved in caring for a patient who is being supported by an EXCOR® *Pediatric VAD* (referred to as EXCOR in this instruction for use).








The PM provides recommendations on treatment and application of the EXCOR in conjunction with the *Stationary Driving Unit Ikus* (referred to as *Ikus* in this PM). To ensure patient safety and comfort, please read this PM carefully.

Always make sure that only professional medical personnel who have been specifically trained in the use of the product are permitted to work with EXCOR.

Note: The recommendations in this PM are based on *Berlin Hearts* experience with the EXCOR. The decisions related to implantation, the components to be used, and patient care remain with the patients physicians.

Note: The technical aspects of *Ikus* are described in the EXCOR® *Pediatric VAD* Instructions for Use from Rev. 4 (IFU 1000721x05) forward. This PM applies exclusively in connection with the IFU.

The following pictograms and symbols are used in this instruction for use:

 DANGER	Indicates a hazardous situation which, if not avoided, will result in death or serious injury to the patient.
 WARNING	Indicates a hazardous situation which, if not avoided, could result in death or serious injury to the patient.
 CAUTION	Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury to the patient and/ or damage to the device.
 NOTICE	Notes are practices not related to personal injury. Possible damage to the device.
 ADVICE	This symbol identifies measures and procedures which have proved useful and successful in conjunction with EXCOR and which we therefore recommend.
 HOTLINE	866.249.0128 This is the telephone number of the emergency hotline. The hotline desk is in operation 24 hours a day. This number is intended for use by medical personnel and should be used in cases of emergency only.
 INSTRUCTION	

1. Individual steps of the instructions are numbered in sequential order.

Definition of the used font formates

Description	Meaning
bold, blue	software texts (messages and menus) except in headings and lists
<i>cursive</i>	proper names (except in headings and in registered trademarks)
"text"	quotation
<key>	key on the laptop keyboard
<<filler text>>	e.g. if texts in error messages are various
[dimension unit]	dimension units in tables; e.g. [mmHg]

1 Important safety information

NOTE: This chapter omits safety instructions, information and procedures that refer to the *Ikus* exclusively. Please refer also to the IFU.

1.1 Warnings



Before using EXCOR, read the PM and the IFU carefully.

Only qualified medical personnel trained specifically in the use of the system are permitted to work with EXCOR. Training courses can be arranged with *Berlin Heart, Inc.* Use by untrained personnel can pose a risk to the patient and the EXCOR.

On the system EXCOR only use components of this system. Never use other components than those delivered by *Berlin Heart GmbH/ Berlin Heart Inc.* Otherwise the warranty is no longer valid.

The system *EXCOR Pediatric* and its components are permitted to be used only by prescription of the attending physician.

Unintended use can pose a risk to the patient and the EXCOR.

Do not use the EXCOR if there is any visible damage of the *Ikus* or any of its components.

If there is any malfunction of the *Ikus* while the driving unit is connected to the patient, the *Ikus* must immediately be replaced.

1.1.1 Storage and durability



The expiration date of each *EXCOR* product is found on the product labels located on both the outer and inner packaging. The pumps, cannulae and accessories must not be used after the expiration date and even not be re-sterilized. Otherwise there is a risk of patient infection.

An *EXCOR* blood pump may not be used on a patient for more than 1 year. After this it shall be replaced with new products.

1.1.2 Device configurations



EXCOR was not designed to be used in combination with other systems, nor do any of the currently granted approvals allow for this. Use by untrained personnel poses a risk to the patient and to the *EXCOR*.

In univentricular operation: Always connect the driving tube of the blood pump to the red marked connector.

The units may only be operated with the disposable products and accessories specified in this document. Also see section 12.1.2: Overview: Relationship: body weight – pump size, page 111. Otherwise there is a risk of functional limitation and/or damage to the *Ikus*. Failure to observe this stipulation will invalidate all warranty agreements by *Berlin Heart Inc.*

The connection between the connector *External alarm* (Nurse call) and the internal alarm system of the clinic is not failsafe. The use of this feature does not release the user from supervising the *Ikus* and the displayed messages and alarms.

1.1.3 Procedural techniques - *Ikus*



Whenever the *Ikus* is running in battery operation, the patient must be accompanied by a person trained to use the manual pump. Thus the patient shall be guaranteed care in an emergency.

1.1.4 Packaging and sterilization



EXCOR blood pumps and cannulae are intended for single-use only. Otherwise there is a risk of infection.

The sterile components are sterilized using ETO and are packed in a double-layer sterile package. Check that the various layers of the sterile packaging are not damaged in any way before they are opened. Do not use the components if either of the sterile packages are damaged. The same applies to sterile components which have exceeded the expiration date as printed on the label. Otherwise there is a risk that the product is no longer sterile.

EXCOR sterile components may not be resterilized by the user. Any opened product must be used or sent back to Berlin Heart. If product expires please contact Berlin Heart for exchange.

An aluminum-coated external packaging protects the *Carmeda*[®] *BioActive Surface (CBAS)* of the blood pump and its sterile packaging against fluctuations in relative humidity. Do not use blood pumps with damaged external packaging. Otherwise there is a risk that the *CBAS* coating may be compromised.

The following items are delivered in sterile condition: blood pumps, cannulae, connecting sets, driving tubes, de-airing set, de-airing hammer, tube connecting set, membrane set .

The external packaging and the outer surface of the outer sterile packaging are not sterile. These 2 packaging layers must be removed before the inner sterile packaging containing the product is handed over to the sterile field. Otherwise there is a risk that the sterile field will be contaminated.

1.1.5 Procedural techniques - pumps, cannulae, accessories



The preparation and use of blood pumps should only be performed by trained personnel. Surgical, nursing and perfusion personnel without experience in the use of EXCOR must complete the EXCOR Training Course which provides theoretical introduction and hands-on practical exercises in the operation of this system. The training program is organized and offered by *Berlin Heart, Inc.*

Only use sterile components which have been delivered in undamaged sterile condition (sterile packaging intact, expiration date not expired).

Only use blood pumps which have an undamaged aluminum-coated outer packaging.

The long-term storage conditions for all sterile products must be observed: temperature +15°C to 25°C, relative humidity: 35 % to 50 %. Store in a dry place! Otherwise there is a risk that the product is no longer sterile.

In order to prevent infection, use strict aseptic techniques during implantation and exercise extreme caution throughout the period of EXCOR cardiac support. Danger of infection!

The distal end of the cannulae can be trimmed. At least 5 cm (2 inches) of material without polyester velour covering should remain to allow visual inspection of the cannula/ titanium-connector junction. Otherwise there is a risk that possible deposits if formed, cannot be visualized.

Ensure proper placement of the cannulae, especially with respect to orientation of the LV apex cannula, to prevent suction of the myocardial wall.

Prior to initial operation of the blood pump(s) minimal initial start parameters have to be set on the laptop to ensure smooth transition from CPB to VAD support.

When connecting the blood pump(s) to the cannulae always observe the arrows on the inflow and outflow stubs. They show the blood flow direction. There is a risk of injury to the patient and severe pump malfunction if the titanium connectors on the end of the inflow and outflow stubs are not connected to the appropriate cannulae.

Do not touch or manipulate the blood pumps and cannulae with pointed or sharp-edged objects (surgical instruments, wire brushes, etc.). Otherwise there is a risk of blood pump and cannula leakage.

Creating a transcutaneous tunnel for the LV apex cannula: Always use cannula tunnelling tip, never use a sharp surgical instrument directly on the cannula.

If an *EXCOR connecting set* is required for implantation and the length of the tube part needs to be reduced, the tube part should be cut but only to achieve the following minimum lengths:

Part Number	Diameter Reduction	Minimum Length
A12-016	16 to 12 mm	90 mm
A09-012	12 to 9 mm	75 mm
A06-009	9 to 6 mm	60 mm

Tab. 1-1 Connector set: minimum length of connector tube

⚠ WARNING

Follow exactly the instructions for using the de-airing set. Otherwise there is a risk of membrane damage.

Ensure that cannulae, blood pump(s) and driving tubes are not subject to external forces, like compression, traction or torsion forces, and are free of knots or sharp bends. Prevent the cannulae and connectors from being exposed to tensile forces. Otherwise there is a risk of obstruction of the air and blood flow.

When positioning the driving tubes follow hospital policies to mitigate the risk of adverse tubing and line incidents by routing the driving tubes in a clear pattern toward the feet and to the side.

Do not initiate cardiac support with the EXCOR blood pumps until the blood pumps have been completely de-aired. After connecting the cannulae, ensure removal of all air that is still in the atria or ventricle by performing single steps (**step left, step right**) with subsequent removal of the bubbles inside the pump via the de-airing needle. Otherwise there is a risk of embolism.

When removing the de-airing needle, never pull on the de-airing tube, but rather only on the de-airing needle.

Once the de-airing needle has been removed it cannot be re-inserted.

Rates < 60 bpm are intended to be used only for implantation and explantation. Never use the *Ikus* with a rate < 60 bpm without constant supervision.,

Under circumstances, the messages **Please check left pump and driving tube!** or **Please check right pump and driving tube!** are not generated with the 10 ml EXCOR blood pump due to the low volume of air which is moved in the pump. Therefore in pumps of this size, pay special attention to the movement of the membrane and ensure that each pump fills and empties completely.

Secure each connection between blood pump and cannula with at least one cable tie as soon as the proper function of the EXCOR is established (see section 6.12: Securing the connections, page 70). Otherwise there is a risk of loose connections and inadequate blood supply to the patient..

At least every 4 hours, visually check that the blood pump(s) is (are) filling and ejecting completely over a period of several pump cycles. If a pump is not filling and/ or ejecting completely, institute the appropriate corrective action.

Do not kink the drivelines. Otherwise there might not be sufficient pump output.

In no case should the cannulae either be kinked directly at the connector to the blood pump or at the transition area between velour and silicone.

Do not kink the cannulae needlessly. Otherwise there might not be sufficient pump output. Moreover, cannulae might be damaged.

Wound care and treatment: Before cleaning the wound (see 8.3: Cleaning of the wound, page 77), put on sterile disposable gloves, cap and mask.

Weaning: If the patient does not meet the eligibility criteria at any time during the weaning process: Resume pumping at rate prior to any weaning (initial rate, IR).

1.1.6 System



If a non-matching pump-cannula-combination (see section 12.1.10: Pump-cannula combinations, page 114) was chosen, use only the connector sets provided with the system in order to minimize the risk of clots at the junctions. Be aware of increased risk of thrombosis and hemolysis.

The cannula diameter may be adapted only once (either by using a staged cannula or a connector set.) Multiple staging could result in limited pump performance and compromised hemodynamics.

If the *Ikus* is operating in emergency pulse mode, immediately visually check whether the blood pump(s) is (are) filling and ejecting completely. If one pump is not filling and/or ejecting completely, the patient must be supported immediately using the manual pump (see section 10.2: Driving blood pump(s) with the manual pump, page 92). Otherwise there is a risk that the patient will not be supported sufficiently.

1.1.7 Procedures to minimize risk of thrombosis



Ensure complete filling/ejection of the pump.

When using staged cannulae or a connecting set, the pumping rate may not be greater than the respective value found in Tab. 12-9, as the pump will not eject its full volume at higher rates.

At least every 4 hours, visually check of blood pump(s) for deposit formation.

1.1.8 Cleaning the components



Cleaning the pump and the drive line: Do not use any acetone or petroleum based products near the pump or drivelines. We recommend using only water or alcohol to clean the pump and the drive line.
IMPORTANT: Do not use any corrosive or colored solutions or organic solvents to clean the blood pump or drivelines as they may alter the surface of the product.

Cleaning the cannulae and transcutaneous exit site: Do not use any acetone or petroleum based products near the the cannulae and the transcutaneous exit site. We recommend using chlorhexidine to clean the cannulae and transcutaneous exit site.
IMPORTANT: Do not use any corrosive or colored solutions or organic solvents to clean the cannulae and transcutaneous exit site as they may alter the surface of the product.

1.1.9 Errors and corrective measures



Any time an error message has occurred, visually check that the blood pump(s) is (are) filling and ejecting completely over a period of several pump cycles, then address the error message with the appropriate corrective action.

If the emergency pulse mode is activated while the backup system is already active, the *Ikus* is no longer able to drive both pumps. In this case, the patient must immediately be supported using the manual pump (see section 10.2: Driving blood pump(s) with the manual pump, page 92). Otherwise there might not be sufficient pump output.

In order for a driving tube to be replaced, the pump must be stopped for a short time. If the left driving tube is being replaced in a driving unit providing biventricular support, the right pump must also be stopped while the driving tube is being replaced in order to avoid overloading of the pulmonary circulation (danger of pulmonary edema).

If the left pump is being replaced in a VAD providing biventricular support, the right pump must also be stopped while the pump is being replaced in order to avoid overloading the pulmonary circulation (danger of pulmonary edema).

If the *Ikus* is operating in emergency pulse mode, the user must immediately visually check the blood pump(s) to determine whether the pump(s) are filling and ejecting completely. If one pump is not filling and/or ejecting completely the patient must be supported immediately with the replacement *Ikus*. Use the manual pump while securing the replacement *Ikus* (see IFU and section 10.2: Driving blood pump(s) with the manual pump, page 92 resp.). Otherwise there is the risk that the patient will not be supported sufficiently.

If the emergency pulse mode is activated while the backup system is already active, the *Ikus* is no longer able to drive both pumps. In this case the patient must be supported immediately with the replacement *Ikus*. Use the manual pump while securing the replacement *Ikus* (see IFU and section 10.2: Driving blood pump(s) with the manual pump, page 92 resp.). Otherwise there is the risk that the patient will not be supported sufficiently.

1.1.10 Replacing the blood pump(s)



When replacing a blood pump, follow the instruction given here. Otherwise the duration of the pump stop will be prolonged and the patient might suffer from inadequate support.

The blood pump may only be replaced under sterile conditions!

When connecting the blood pump(s), pay attention to the direction of the arrows on the inflow and outflow stubs! These show the direction of the blood flow.

The cable tie covering the *EXCOR* cannula on the stub of the blood pump should be removed carefully. Important: never use a sharp instrument, for example, a scalpel or scissors, to remove the cable tie. This may cause damage to the cannula.

If the left pump is being replaced in a VAD providing biventricular support, the right pump must also be stopped while the pump is being replaced in order to avoid overloading the pulmonary circulation (danger of pulmonary edema).

1.1.11 Driving blood pump(s) with the manual pump



The use of the manual pump is only permitted for medical personnel trained in the use of it.

Pay attention to the colored markings on the driving tubes and on the connectors of the manual pump. Otherwise, there is a risk of lung edema.

Always keep manual pump attached to the *Ikus*. Otherwise in an emergency situation the adequate support of the patient is not guaranteed.

Call one or more persons to assist. Otherwise in an emergency situation the adequate support of the patient is not guaranteed.

The driving tubes and cannulae should be arranged in a bend-free position. Otherwise in an emergency situation the adequate support of the patient is not guaranteed.

When operating the manual pump with 1 hand, do not block the valves with your feet (see valve "c" in Fig. 10-3, page 94 and in Fig. 10-4, page 94).

1.1.12 Ambient conditions



Protect the *Ikus* from exposure to moisture and wetness. Never store or operate the *Ikus* in a damp environment (e.g. bathroom, etc.). Otherwise there is a risk of functional limitation and/or *Ikus* malfunction.

In terms of electromagnetic compatibility (EMC) the *Ikus* is subject to special precautions! Avoid exposure to strong electromagnetic radiation (as generated by mobile/cell phones and cordless phones when switched on, electromagnetic security systems etc.), see IFU. Otherwise there is a risk of electromagnetic disturbances and fault-free functioning of the *Ikus* cannot be guaranteed.

When using a cell phone in the immediate environment of an *Ikus* in operation please make sure to keep a distance of at least 0.77 m. For further information please refer to IFU.

When using an RFID device in the immediate environment of an *Ikus* in operation please make sure to keep a distance of at least 1 m. For further information please refer to IFU.

If an ambient temperature of +30°C is continuously exceeded during operation, the lifetime of the batteries is reduced. Therefore, a person trained to use the manual pump should always be present in this case. This should ensure patient care in case of emergency.

Use the *Ikus* as far away as possible from environments containing flammable gases and use extreme caution. Otherwise there is a risk of explosion or gas ignition. The *Ikus* would be severely limited in function or malfunction altogether as a result of this damage.

Also see IFU.

1.1.13 Interaction with other procedures and therapies



The following procedure is not possible:

- Magnetic resonance imaging
-

EXCOR patients with prosthetic aortic valves may have increased risk of thromboembolism.

If EXCOR is used in interaction with other procedures and therapies, observe the movement of the membrane to determine whether the blood pump is filling and ejecting completely. If a pump is not filling and/ or ejecting completely, stop the interacting procedure or therapy and institute the appropriate corrective action.

In terms of electromagnetic compatibility (EMC) the *Ikus* is subject to special precautions! When exposing *Ikus* to the procedures and therapies listed below please observe EMC regulations given in the IFU.

For the following procedures and therapies, the manufacturer does not expect any harmful interaction with the *Ikus* due to the general electromagnetic shielding of the device (see IFU). However, these procedures and therapies must only be applied after consultation with the treating physician.

- Radiotherapy
 - Nuclear diagnostics / nuclear therapy
 - Electro-stimulation therapy
 - Therapeutic ultrasonic treatment (e.g. lithotripsy)
 - External defibrillation
-

The following procedures and therapies have been tested in regard to their interaction with the *Ikus* and no harmful effects were found, however, these procedures and therapies must only be applied after consultation with the treating physician. Additionally the manufacturer does not guarantee that equivalent devices will not interfere.

- Diathermy
 - X-rays
 - Computed tomography
-

1.2 Precautions

1.2.1 VAD placement technique



Implantation - anesthesia: There should be an adequate supply of pre-matched stored blood, fresh frozen plasma and platelet concentrates available for immediate transfusion if required.

Implantation - anesthesia: Keep blood product transfusions to a minimum. Blood transfusions may lead to the development of antibodies, which are known to promote coagulation and inflammatory response.

The titanium connectors of the blood pumps have sharp edges designed to minimise the risk of clot formation at the junction. Be careful to avoid cutting yourself while connecting the pump and the cannulae.

1.2.2 Ambient conditions



The *Ikus* is intended solely for use in a hospital setting.

1.2.3 Caution while using the *Ikus*



At least daily, the *EXCOR* cannulae should be inspected for signs of wear or damage. **ADVICE:** To avoid needless kinking of the cannulae use a mirror for inspection of the bottom side of the blood pump.

At least every 4 hours, check visually that the blood pump(s) is (are) filling and ejecting completely over a period of several pump cycles. If a pump is not filling and/ or ejecting completely, then take the appropriate corrective action.

Under certain circumstances, the message **left/right pump is not filling adequately** in some circumstances is not generated with the 10 ml EXCOR blood pump due to the low volume of air which is moved in the pump. Therefore in pumps of this size, pay special attention to the movement of the membrane and ensure that each pump fills and empties completely.

After changing over to biventricular operation the device is operating in separate mode. All parameters are reset to the default parameters (see IFU). The patient-customized parameters have to be adjusted again.

Replacing the blood pump due to growth of the patient: In children, plan to replace the pump(s) with a larger pump(s) in good time, to prevent the possibility of inadequate support due to an insufficient discharge rate.

1.3 Obligations of the operator



Only qualified medical personnel trained specifically in the use of the system are permitted to work with EXCOR. Training courses can be arranged with *Berlin Heart, Inc.*



The operator (i.e. the hospital using the system) is responsible for instruction and care of the patient. The patient must be instructed on safety risks and cautionary measures (moisture, temperature, electromagnetic fields, etc.).

A replacement *Ikus* and replacement equipment must always be available in the hospital.

2 General Information

2.1 Device description

EXCOR is an extracorporeal, pneumatically driven ventricular assist device. It is designed to support the right and/or left ventricle when the native heart is unable to maintain normal blood flows and pressures even with help of drug therapy and intraaortic balloon counterpulsation. The device is designed for mid to long term mechanical support.

The EXCOR consists of 1 or 2 extracorporeal, pneumatically driven blood pumps and cannulae which connect the blood pump(s) to the atrium or ventricle and to the great arteries. The *Ikus* provides alternating air pressure to the blood pumps through driving tubes.

The blood pump is divided into an air chamber and a blood chamber by a multi-layer flexible polyurethane membrane. The alternating air pressure provided by the *Ikus* moves the membrane, thus filling and emptying the blood pump. Both the blood chamber and the polyurethane connectors are transparent to allow for visual detection of deposits and for monitoring the filling and emptying of the blood pump.

Valves (three-leaflet polyurethane valves) are located at the inlet and outlet positions of the blood pump connector stubs, thus ensuring the unidirectional blood flow.

Pulse rate, systolic drive pressure, diastolic suction pressure and the relative systolic duration can all be monitored and adjusted on the driving unit.

2.2 Indications for use

The EXCOR is intended to provide mechanical support as a bridge to cardiac transplantation for pediatric patients. Pediatric patients with severe isolated left ventricular or biventricular dysfunction who are candidates for cardiac transplant and require circulatory support may be treated using the EXCOR.

2.3 IDE Clinical Study Summary

See chapter 3: IDE Clinical Study Summary, page 23.

2.4 Intended operation environment

Ikus is intended for use in a clinical setting. It can be used in any kind of hospital unit, e.g. OR, ICU, intermediate care unit or general care unit. It may be moved between clinical units using the built-in wheels, however in this case the patient must always be accompanied by a person trained in the use of the manual pump and emergency procedures. Thus, the patient shall be guaranteed care in case of an emergency.

Transporting the device during operation by any vehicles (e.g. ambulance, aircraft, etc.) is not allowed.

During movement of the device in operation within the clinic all electromagnetic compatibility precautions (EMC precautions) must be observed. See IFU. Otherwise there is a risk of electromagnetic disturbances and the fault-free operation of *Ikus* could not be guaranteed.

2.5 Contraindications

Patients unable to tolerate systemic anticoagulation therapy should not be implanted.

Magnetic Resonance Imaging (MRI) is contraindicated in patients after being implanted with the EXCOR.

2.6 Storage and durability

⚠ WARNING

The expiration date of each *EXCOR* product is found on the product labels located on both the outer and inner packaging. The pumps, cannulae and accessories must not be used after the expiration date and even not be re-sterilized. Otherwise there is a risk of patient infection.

An *EXCOR* blood pump may not be used on a patient for more than 1 year. After this it shall be replaced with new products.

IMPORTANT: EXCOR must be stored at room temperature and be protected against extreme temperature fluctuations and moisture. Otherwise there is a risk of functional limitation and/or damage to the Iku.

3 IDE Clinical Study Summary

3.1 Indications for use

EXCOR® Pediatric Ventricular Assist Device (referred to as EXCOR) is intended to provide mechanical circulatory support as a bridge to cardiac transplantation for pediatric patients. Pediatric candidates with severe isolated left ventricular or biventricular dysfunction who are candidates for cardiac transplant and require circulatory support may be treated using the EXCOR.

3.2 Contraindications

Patients unable to tolerate systemic anticoagulation therapy should not be implanted.

Magnetic Resonance Imaging (MRI) is contraindicated in patients after being implanted with the EXCOR.

3.3 Alternative Practices or Procedures

FDA approved therapies include the DeBakey Child device for left ventricular support for body surface area $> 0.7 \text{ m}^2$ and $< 1.5 \text{ m}^2$. EXCOR is the only ventricular assist device approved for univentricular and biventricular support in children from 3-60 kg.

3.4 Marketing History

EXCOR was approved to apply the CE Mark in 1996. Since that authorization, EXCOR has been distributed to the following countries: Germany, Austria, Belgium, Bulgaria, Estonia, Switzerland, Denmark, Spain, Finland, France, Great Britain, Greece, Hungary, Italy, Lithuania, Netherlands, Poland, Portugal, Romania, Sweden, Slovakia, Turkey, Argentina, Australia, Azerbaijan, Brazil, Canada, Chile, Taiwan, China, Hong Kong, Israel, Iran, New Zealand, Serbia, Russia, Saudi Arabia, and South Africa. The EXCOR has not been removed from the market in any country.

3.5 Potential Adverse Effects

Serious adverse events (SAEs) for all primary cohort patients were reported in the primary study analysis for events per patient-day. The total time on device for Cohort 1 (BSA $< 0.7 \text{ m}^2$) subjects of 1411 days yielded a rate of 0.068 SAEs per patient-day. The total time on device for Cohort 2 (BSA > 0.7 to $< 1.5 \text{ m}^2$) subjects was 1376 days yielded a rate of 0.079 SAEs per patient-day.

The following table details each SAE with the number of events experienced and the number and percent of subjects experiencing each SAE. Some of the SAEs have subcategories (see indented descriptions) which provide additional detail regarding the type of SAE.

Rates for subjects enrolled in the Cohorts 1 CAP (Continued Access Protocol which allowed continued access to the device following the conclusion of enrollment in the primary cohorts) and Compassionate / Emergency Use Cohorts 3A and 3B are included to support the assessment of reasonable assurance of safety as specified in the IDE Investigational Plan.

Serious Adverse Event Summary per Cohort

EVENT	COHORT										
	1 Total	Per Subject (% of 24)	1 CAP Total	Per Subject (% of 20)	3A Total	Per Subject (% of 35)	2 Total	Per Subject (% of 24)	3B Total	Per Subject (% of 6)	
Major Bleeding	15	10 (41.7%)	12	7 (35.0%)	25	18 (51.4%)	22	12 (50.0%)	3	3 (50.0%)	
Cardiac Arrhythmia	1	1 (4.2%)	2	2 (10.0%)	3	3 (8.6%)	6	4 (16.7%)	2	1 (16.7%)	
Sustained VT	1	1 (4.2%)	0	0 (0.0%)	2	2 (5.7%)	2	2 (8.3%)	2	1 (16.7%)	
Sustained SVT	0	0 (0.0%)	2	2 (10.0%)	1	1 (2.9%)	4	3 (12.5%)	0	0 (0.0%)	
Pericardial Fluid Collection	3	3 (12.5%)	5	5 (25.0%)	4	4 (11.4%)	4	3 (12.5%)	1	1 (16.7%)	
With Tamponade	1	1 (4.2%)	3	3 (15.0%)	2	2 (5.7%)	2	2 (8.3%)	0	0 (0.0%)	
Without Tamponade	2	2 (8.3%)	2	2 (10.0%)	2	2 (5.7%)	2	2 (8.3%)	1	1 (16.7%)	
Hemolysis	1	1 (4.2%)	1	1 (5.0%)	1	1 (2.9%)	1	1 (4.2%)	1	1 (16.7%)	
Hemolysis-Early	0	0 (0.0%)	0	0 (0.0%)	0	0 (0.0%)	0	0 (0.0%)	1	1 (16.7%)	
Hemolysis-Late	1	1 (4.2%)	1	1 (5.0%)	1	1 (2.9%)	1	1 (4.2%)	0	0 (0.0%)	
Hepatic Dysfunction	1	1 (4.2%)	0	0 (0.0%)	6	5 (14.3%)	1	1 (4.2%)	3	2 (33.3%)	
Hypertension	12	12 (50.0%)	15	13 (65.0%)	9	9 (25.7%)	8	8 (33.3%)	1	1 (16.7%)	
Major Infection	35	15 (62.5%)	15	7 (35.0%)	39	16 (45.7%)	24	12 (50.0%)	8	4 (66.7%)	
Infection-Localized Non-Device	25	12 (50.0%)	10	6 (30.0%)	20	11 (31.4%)	18	10 (41.7%)	7	3 (50.0%)	
Infection-Percutaneous Site or Pocket	4	4 (16.7%)	1	1 (5.0%)	0	0 (0.0%)	0	0 (0.0%)	0	0 (0.0%)	
Infection-Septis	6	5 (20.8%)	4	2 (10.0%)	19	9 (25.7%)	6	6 (25.0%)	1	1 (16.7%)	
Psychiatric Episode	0	0 (0.0%)	0	0 (0.0%)	0	0 (0.0%)	1	1 (4.2%)	0	0 (0.0%)	

Tab. 3-1 Serious adverse event summary per cohort

Serious Adverse Event Summary per Cohort, continued

EVENT	COHORT										
	1 Total	Per Subject (% of 24)	1 CAP Total	Per Subject (% of 20)	3A Total	Per Subject (% of 35)	2 Total	Per Subject (% of 24)	3B Total	Per Subject (% of 6)	
Neurological Dysfunction	8	7 (29.2%)	6	5 (25.0%)	6	6 (17.1%)	9	7 (29.2%)	4	3 (50.0%)	
TIA	0	0 (0.0%)	1	1 (5.0%)	0	0 (0.0%)	0	0 (0.0%)	1	1 (16.7%)	
Ischemic CVA	8	7 (29.2%)	5	5 (25.0%)	4	4 (11.4%)	7	7 (29.2%)	3	3 (50.0%)	
Hemorrhagic CVA	0	0 (0.0%)	0	0 (0.0%)	2	2 (5.7%)	2	2 (8.3%)	0	0 (0.0%)	
Renal Dysfunction	3	2 (8.3%)	0	0 (0.0%)	7	7 (20.0%)	4	3 (12.5%)	2	1 (16.7%)	
Acute	3	2 (8.3%)	0	0 (0.0%)	7	7 (20.0%)	2	2 (8.3%)	2	1 (16.7%)	
Chronic	0	0 (0.0%)	0	0 (0.0%)	0	0 (0.0%)	2	2 (8.3%)	0	0 (0.0%)	
Respiratory Failure	3	3 (12.5%)	8	8 (40.0%)	6	5 (14.3%)	9	6 (25.0%)	6	5 (83.3%)	
Right Heart Failure	2	2 (8.3%)	2	2 (10.0%)	8	7 (20.0%)	3	3 (12.5%)	1	1 (16.7%)	
Arterial Non-CNS Thromboembolism	1	1 (4.2%)	1	1 (5.0%)	2	2 (5.7%)	0	0 (0.0%)	0	0 (0.0%)	
Venous Thromboembolism Event	1	1 (4.2%)	1	1 (5.0%)	0	0 (0.0%)	0	0 (0.0%)	0	0 (0.0%)	
Wound Dehiscence	0	0 (0.0%)	0	0 (0.0%)	1	1 (2.9%)	0	0 (0.0%)	0	0 (0.0%)	
Other	10	6 (25.0%)	6	5 (25.0%)	17	12 (34.3%)	15	6 (25.0%)	7	4 (66.7%)	
Other Ischemic w/o symptoms	0	0 (0.0%)	0	0 (0.0%)	1	1 (2.9%)	0	0 (0.0%)	0	0 (0.0%)	
Other Covert Stroke	0	0 (0.0%)	0	0 (0.0%)	0	0 (0.0%)	0	0 (0.0%)	1	1 (16.7%)	

Tab. 3-2 Serious adverse event summary per cohort (table continued)

The rates of SAEs per patient-day were calculated separated by whether the subjects were supported with ECMO pre-implant and are summarized in the following table.

In Cohort 1, those supported with ECMO pre-implant had twice as many events per patient-day of support. For Cohort 2, those supported with ECMO pre-implant had 1.5 times as many events per patient-day of support.

Serious Adverse Events per Patient-day by pre-implant ECMO

Group	ECMO Pre-Implant	# Events	Total Time on Support (Days)	Rates Success Criterion <0.25	
				Events per Patient-Day	Upper bound of CI
Cohort 1	Yes	38	345	0.110	0.151
	No	58	1066	0.054	0.070
Cohort 2	Yes	43	450	0.096	0.129
	No	64	926	0.069	0.088

Tab. 3-3 Serious adverse events per patient-day pre-implant ECMO

3.6 Summary of Clinical Studies

3.6.1 IDE Clinical Study Summary

Berlin Heart Inc. conducted a prospective, multi-center, single arm study to assess the safety and probable benefit of the EXCOR.

The purpose of the study was to determine whether use of the EXCOR for bridge-to-transplantation is associated with reasonable assurance of safety and probable benefit such that the EXCOR merits approval by the Food and Drug Administration (FDA) under a Humanitarian Device Exemption (HDE).

3.6.2 Study Cohorts

The primary study population of 48 subjects aged 0-16 years consisted of 24 subjects with a body surface area (BSA) < 0.7 m² (Cohort 1) and 24 subjects with a body surface area (BSA) ≥ 0.7 m² to < 1.5 m² (Cohort 2).

A third cohort of subjects was enrolled under Compassionate / Emergency Use and is classified as Cohort 3. These subjects followed the study protocol unless otherwise noted within the approval documentation for the subject. This cohort is further divided into groups based on the subject’s BSA similar to Cohorts 1 and 2 and is labeled Cohort 3A if the subject’s BSA is < 0.7 m² and Cohort 3B if the BSA is ≥ 0.7 m² and <1.5 m².

For the primary effectiveness endpoint, the protocol prescribed an ECMO historical control group. The historical ECMO control group was compiled from the Extracorporeal Life Support Organization (ELSO) registry, the most extensive registry of patients treated with ECMO in North America. The database was filtered to best match the EXCOR IDE study population. Patients included for comparison to the EXCOR cohorts included patients from both genders, age 0-16 years, with weight greater than 3 kg, cardiac only ECMO support, support initiation from 2000 onward who met critical eligibility criteria. The dataset for the ELSO registry included baseline and outcomes data comparable to the EXCOR dataset. The control group was then created by matching the EXCOR subjects to the patients in the subset using a propensity score analysis (PSA).

3.6.3 Inclusion/Exclusion Criteria

Subjects of both genders who satisfy all inclusion and exclusion criteria were eligible for entrance into the primary cohorts of the clinical study.

Inclusion Criteria

1. Severe NYHA Functional Class IV (or Ross Functional Class IV for subjects ≤ 6 years) heart failure refractory to optimal medical therapy, and has met at least one of the following criteria:
 - A INTERMACS™ profile status 1 or 1A, i.e. critical cardiogenic shock (low BP unresponsive to support, compromised end organ perfusion, < 24 hour survival expected without mechanical support; may be due to VT/VF (1A)
 - B INTERMACS profile status 2 or 2A (i.e. progressive decline): not in imminent danger, but worsening despite optimal inotropic therapy; may be due to VT/VF (2A) AND at least one of the following criteria
 - a Decline in renal function as defined by a 50 % reduction in estimated GFR despite optimization of subject volume status
 - b Decline in nutritional status as defined by a sustained (≥ 7 days) inability to tolerate an enteral nutritional intake sufficient to provide at least 75 % of the prescribed caloric needs for the subject, or signs of nutritional compromise (cachexia, nutritional weight loss) despite appropriate intervention
 - c Decline in mobility/ambulation as defined by sustained bed confinement (≥ 7 days without prospect for improvement) attributable to heart failure symptoms or its treatment (e.g. intubation for pulmonary edema)
 - C Support with extra-corporeal membrane oxygenation (ECMO) or other mechanical circulatory support device OR
 - D Unable to separate from cardiopulmonary bypass (must be listed for heart transplantation at time of transfer to the operating room)
2. Listed (UNOS status 1A or equivalent) for cardiac transplantation
3. Two-ventricle circulation, including cardiomyopathy, repaired structural heart disease (e.g. ALCAPA, aortic stenosis) or acquired heart disease (e.g. myocarditis, Kawasaki disease)
4. Age 0 to 16 years; corrected gestational (CGA) at least 37 weeks
5. Weight ≥ 3 kg and ≤ 60 kg
6. Legal guardian (and subject if age-appropriate) understands the nature of the procedure, are willing to comply with associated follow-up evaluations, and provide written informed consent and assent prior to the procedure

Exclusion Criteria

1. Support on ECMO for ≥ 10 days
2. Cardiopulmonary resuscitation (CPR) duration ≥ 30 minutes within 48 hours prior to device implantation
3. Body weight < 3.0 kg or BSA > 1.5 m²
4. Presence of mechanical aortic valve
5. Unfavorable or technically-challenging cardiac anatomy including single ventricle lesions, complex heterotaxy, and restrictive cardiomyopathy
6. Evidence of intrinsic hepatic disease as defined by a total bilirubin level or AST/ALT greater than five times the upper limit of normal for age, except in association with acute heart failure as determined by the principal investigator
7. Evidence of intrinsic renal disease as defined by a serum creatinine greater than 3 times the upper limit of normal for age, except in association with acute heart failure as determined by the principal investigator
8. Hemodialysis or peritoneal dialysis (not including dialysis or Continuous Veno-Venous Hemofiltration (CVVH) for volume removal
9. Evidence of intrinsic pulmonary disease (e.g. chronic lung disease, RDS) as defined by need for chronic mechanical ventilation, except in association with acute heart failure as

- determined by the principal investigator
10. Moderate or severe aortic and/or pulmonic valve insufficiency considered technically challenging to repair at the time of the device implantation as determined by the principal investigator
 11. Apical VSD or other hemodynamically-significant lesion considered technically challenging to repair at the time of device implantation as determined by the principal investigator
 12. Documented heparin induced thrombocytopenia (HIT) or idiopathic thrombocytopenia purpura (ITP) or other contraindication to anticoagulant/antiplatelet therapy
 13. Documented coagulopathy (e.g. Factor VIII deficiency, disseminated intravascular coagulation) or thrombophilic disorder (e.g. Factor V Leiden mutation)
 14. Hematologic disorder causing fragility of blood cells or hemolysis (e.g. sickle cell disease)
 15. Active infection within 48 hours of implant demonstrated by:
 - A Positive blood culture OR
 - B Temperature >38 degrees C and WBC >15, 000/ ml
 16. Documented human immunodeficiency virus (HIV) infection or acquired immunodeficiency syndrome (AIDS)
 17. Evidence of recent or life-limiting malignant disease
 18. Stroke within past 30 days prior to enrollment, or congenital CNS malformation syndrome associated with increased risk of bleeding (e.g. arteriovenous malformation, moya moyo)
 19. Psychiatric or behavioral disease (e.g. antisocial disorder) with a high likelihood for non-compliance
 20. Currently participating in another investigational device or drug trial and has not completed the required follow-up period for that study
 21. Subject is pregnant or nursing

3.6.4 Study Enrollment

The following table summarizes the complete enrollment (including the subjects enrolled at non IDE sites) by subject's body size. As of the data cutoff for the final HDE report (February 2011 report with January 17, 2011 data cutoff), there were 151 smaller sized subjects (BSA < 0.7m²) enrolled and 53 larger sized subjects (BSA ≥ 0.7 to <1.5 m²) enrolled.

Subject Enrollment

Cohort	IDE Site Implants	Non-IDE Site Implants	Total
BSA < 0.7 m²			
Cohort 1	24	<i>n/a</i>	24
Cohort 1 CAP	20	<i>n/a</i>	20
Cohort 3A	35	72	107
<i>Subtotal</i>	79	72	151
BSA ≥ 0.7 m² to < 1.5 m²			
Cohort 2	24	<i>n/a</i>	24
Cohort 3B	6	23	29
<i>Subtotal</i>	30	23	53
TOTAL	109	95	204

Tab. 3-4 Subject enrollment

Note: Enrollment in Cohorts 1 CAP, 3A, 3B (IDE and non-IDE) are supportive data and are included only in the safety summary tables.

Study Enrollment and Outcome

<p>Total Enrollment June 21, 2007 -- December 1, 2010 n=204</p>						
<p>BSA < 0.7m² n=151</p> <p>Transplant n=88 Weaned n=10 Death n=45 On device n= 8</p>				<p>BSA ≥ 0.7 m²- < 1.5 m² n=53</p> <p>Transplant n=42 Weaned n= 2 Death n= 6 On device n= 3</p>		
<p>Cohort 1 n=24</p>	<p>Cohort 1 CAP n=20</p>	<p>Cohort 3A IDE Sites n=35</p>	<p>Cohort 3A Non-IDE Sites n=72</p>	<p>Cohort 2 n=24</p>	<p>Cohort 3B IDE Sites n=6</p>	<p>Cohort 3B Non-IDE Sites n=23</p>
<p>TX n=21 Weaned n=1 Death n=2 On Device n=0</p>	<p>TX n=16 Weaned n=0 Death n=1 On Device n=3</p>	<p>TX n=20 Weaned n=3 Death n=10 On Device n=2</p>	<p>TX n=31 Weaned n=6 Death n=32 On Device n=3</p>	<p>TX n=21 Weaned n=1 Death n=2 On Device n=0</p>	<p>TX n= 4 Weaned n=1 Death n=1 On Device n=0</p>	<p>TX n=17 Weaned n=0 Death n=3 On Device n=3</p>

Fig. 3-1 Study enrollment and outcome

Enrollment in Cohorts 1 CAP, 3A, 3B (IDE and non-IDE) are supportive data and are only included in the safety summary tables.

3.6.5 Subject Demographics

The following table summarizes the demographic data for Cohorts 1 and 2. Males comprised the majority of the subjects in Cohort 2 (54%) and half (50%) of Cohort 1. The smaller group of subjects ranged in age from 2.6 to 45.6 months while the larger group ranged in age from 51 to 192 months (or 4.2 to 16 years). The weight range for Cohort 1 was 3.6 to 13.6 kilograms with a BSA range of 0.23 to 0.62 m² and the weight range for Cohort 2 was 16.0 to 58.1 kilograms with a BSA range of 0.71 to 1.66 m².

The most predominant cardiac diagnosis for Cohort 1 was dilated cardiomyopathy (79.2%) and the majority of this group, 54.2%, presented with progressive decline. The most predominant cardiac diagnosis for Cohort 2 was also dilated cardiomyopathy (70.8%) and most (54.2%) were listed as in critical cardiogenic shock.

Demographic Data Summary			
Variable	Category	Cohort 1 n=24	Cohort 2 n=24
Gender	Female	12 (50.0%)	11 (45.8%)
	Male	12 (50.0%)	13 (54.2%)
Age (months)	Mean \pm Std (N)	15.4 \pm 12.4 (24)	113.2 \pm 37.6 (24)
	Median	11.7	111.2
	Min – Max	2.6 - 45.6	50.8 - 191.8
BSA (m ²)	Mean \pm Std (N)	0.43 \pm 0.10 (24)	1.09 \pm 0.29 (24)
	Median	0.44	1.08
	Min – Max	0.23 - 0.62	0.71 - 1.66
Weight (kg)	Mean \pm Std (N)	9.1 \pm 2.7 (24)	32.2 \pm 12.5 (24)
	Median	9.2	30.7
	Min – Max	3.6 - 13.6	16.0 – 58.1
Race	African-American	7 (29.2%)	6 (25.0%)
	American Indian/Alaska Native	1 (4.2%)	0 (0.0%)
	Asian	0 (0.0%)	1 (4.2%)
	Hawaiian/other Pacific Islander	0 (0.0%)	1 (4.2%)
	White	13 (54.2%)	15 (62.5%)
	Other/none of the above	3 (12.5%)	1 (4.2%)
Ethnicity: Hispanic or Latino	Yes	7 (29.2%)	1 (4.2%)

Tab. 3-5 Demographic data summary (a)

Demographic Data Summary, continued

Variable	Category	Cohort 1 n=24	Cohort 2 n=24
Patient Profile/Status	1 Critical Cardiogenic Shock	11 (45.8%)	13 (54.2%)
	2 Progressive decline	13 (54.2%)	11 (45.8%)
	3 Stable but Inotrope dependent	0 (0.0%)	0 (0.0%)
Modifier A Arrhythmia (# Yes)		4 (16.7%)	4 (16.7%)
Primary Cardiac Diagnosis	Congenital Heart Disease	3 (12.5%)	6 (25.0%)
	Dilated Myopathy	19 (79.2%)	17 (70.8%)
	Hypertrophic cardiomyopathy	1 (4.2%)	0 (0.0%)
	Restrictive Myopathy	1 (4.2%)	1 (4.2%)
Secondary Cardiac Diagnosis (multiple Choices)	Congenital Heart Disease	2 (8.3%)	3 (12.5%)
	Coronary Artery Disease	0 (0.0%)	2 (8.3%)
	Dilated Myopathy: Familial	1 (4.2%)	0 (0.0%)
	Dilated Myopathy: Idiopathic	0 (0.0%)	2 (8.3%)
	Dilated Myopathy: Ischemic	0 (0.0%)	1 (4.2%)
	Dilated Myopathy: Myocarditis	0 (0.0%)	2 (8.3%)
	Dilated Myopathy: Viral	1 (4.2%)	0 (0.0%)
	Dilated Myopathy: Other	1 (4.2%)	2 (8.3%)
	Restrict Myopathy: Secondary to Radiation/Chemo	0 (0.0%)	1 (4.2%)
	Valvular Heart Disease	0 (0.0%)	1 (4.2%)
	CHD/Dilated Myopathy Familial	1 (4.2%)	0 (0.0%)
	None	18 (75.0%)	10 (41.7%)
	Heart Rate	Mean ± Std (N)	126.3 ± 25.5 (24)
Min – Max		91.0 - 175.0	85.0 - 168.0
Systolic Blood Pressure	Mean ± Std (N)	85.3 ± 16.0 (24)	95.2 ± 13.5 (24)
	Min – Max	45.0 - 110.0	60.0 - 112.0
Diastolic Blood Pressure	Mean ± Std (N)	56.0 ± 14.1 (24)	65.9 ± 14.8 (24)
	Min – Max	38.0 - 89.0	46.0 - 100.0
Previous Cardiac operations (# Yes)		5 (20.8%)	8 (33.3%)

Tab. 3-6 Demographic data summary (b)

Pre-implant support for the subjects is detailed in the following table. ECMO support was used pre-implant for 25% of Cohort 1 subjects and 33.3% of Cohort 2 subjects.

Pre-Implant Support

Variable	Category	Cohort 1	Cohort 2
		n=24	n=24
Prior support within 48 hours	No support	0 (0.0%)	0 (0.0%)
	Ventilator	20 (83.3%)	12 (50.0%)
	ECMO	6 (25.0%)	8 (33.3%)
	Ultrafiltration	3 (12.5%)	1 (4.2%)
	VAD	2 (8.3%)	0 (0.0%)
	Dialysis	0 (0.0%)	0 (0.0%)
	Feeding Tube	10 (41.7%)	7 (29.2%)
	IABP	0 (0.0%)	0 (0.0%)
	Inotropes	22 (91.7%)	21 (87.5%)

Tab. 3-7 Pre-implant support

3.6.6 Results**3.6.6.1 Probable Benefit**

Efficacy for the IDE trial was assessed by comparing survival (defined by the interval of time from initiation of mechanical support as a bridge to transplant or recovery) to the historical ECMO control. Subjects who were transplanted were censored at the time of explant. Subjects who were explanted due to recovery of their ventricular function and survived to 30 days or discharged with acceptable neurologic status were censored at the time of explant. Subjects who were explanted due to recovery of their ventricular function and died within 30 days or discharge (whichever was longer) were counted as a failure with time to failure being the explant date.

For the 2 primary cohorts, the rate of successfully bridging the subjects to transplant was 87.5% for Cohort 1 (21/24) and 91.7% for Cohort 2 (22/24) or 89.6% overall (43/48). The following table summarizes the survival to transplant/successful recovery for each primary Cohort ITT and PP as well as their matched ECMO control groups.

Three (3) of the Cohort 1 subjects (12.5%) failed (2 deaths and 1 weaned subject with unacceptable neurological outcome at 30 days post-explantation) compared to 12 of the 48 (25%) patients in the matched ECMO control group. The 3 subjects from Cohort 1 who died or were considered failures were all supported with ECMO at the time of implant. The failures occurred at day 0 (death), day 38 (death) and day 146 (weaned-failure).

The control group for Cohort 1 was on ECMO for a median of 4.9 days and a maximum of 20.5 days compared to the primary cohort subjects who were supported a median of 27.5 days and maximum of 174 days. Seventeen (17) of the 24 (71%) Cohort 1 subjects were supported longer than the entire ECMO control group (i.e. longer than 20.5 days).

Two of the Cohort 2 subjects (8.3%) failed compared to 16 of the 48 (33.3%) patients in the matched ECMO control group. One of the subjects who died in Cohort 2 was supported with ECMO at the time of implant. The deaths occurred at day 19 and day 144.

The control group for Cohort 2 was on ECMO for a median of 4.7 days and a maximum of 27.5 days compared to the primary cohort subjects who were supported a median of 42.5 days and a maximum of 192 days. Seventeen (17) of the 24 (71%) subjects in Cohort 2 were supported longer than the entire ECMO control group (i.e. longer than 27.5 days).


Primary Efficacy Study and Control Groups							
Group	Total	Max Time on Device (days)	# Successes	# Failures	Survival Time		
					30 days	60 days	90 days
Cohort 1 ITT	24	174	21 (87.5%)	3 (12.5%)	95.8%	87.1%	87.1%
Cohort 1 Per-Protocol	22	174	19 (86.4%)	3 (13.6%)	95.5%	86.8%	86.8%
ECMO Control Group	48	20.5	36 (75.0%)	12 (25.0%)	NA	NA	NA
Cohort 2 ITT	24	192	22 (91.7%)	2 (8.3%)	94.7%	94.7%	94.7%
Cohort 2 Per-Protocol	22	144	20 (90.9%)	2 (9.1%)	94.1%	94.1%	94.1%
ECMO Control Group	48	27.5	32 (66.7%)	16 (33.3%)	NA	NA	NA

Tab. 3-8 Primary Efficacy Study and Control Groups

Comparison of the ITT groups to their respective matched ECMO control group survival rates were both statistically significant (log-rank p value <0.0001). Therefore, there is a significantly higher survival rate of Cohort 1 and 2 subjects as compared to their respective ECMO control group.

The following figures display the Kaplan-Meier curves for the endpoint of death/weaned with unacceptable outcome for both Cohort 1 ITT and Cohort 2 ITT and their respective ECMO control groups.

**Survival to Death/Weaned with Unacceptable Neurological Outcome:
Cohort 1 versus ECMO**




COHORT 1	Interval Ending (Days Post Implant)									
	0	1	7	14	30	45	60	90	120	150
# Left	24	21	21	20	12	10	9	6	5	1
Total # Failed	0	1	1	1	1	2	2	2	2	3
Survival	100%	95.8%	95.8%	95.8%	95.8%	87.1%	87.1%	87.1%	87.1%	58.1%
Std Error	0%	4.1%	4.1%	4.1%	4.1%	9.1%	9.1%	9.1%	9.1%	24.5%

ECMO CONTROL	Interval Ending (Days Post Implant)				
	0	1	7	14	30
# Left	48	46	16	4	0
Total # Failed	0	2	7	10	12
Survival	100%	95.8%	80.1%	52.0%	17.3%
Std Error	0%	2.9%	7.1%	14.2%	14.9%

Fig. 3-2 Cohort 1 Survival

**Survival to Death/Weaned with Unacceptable Neurological Outcome:
Cohort 2 versus ECMO**



COHORT 2	Interval Ending (Days Post Implant)									
	0	1	7	14	30	45	60	90	120	150
# Left	24	23	21	20	17	11	9	6	3	1
Total # Failed	0	0	0	0	1	1	1	1	1	2
Survival	100%	100%	100%	100%	94.7%	94.7%	94.7%	94.7%	94.7%	47.4%
Std Error	0%	0%	0%	0%	5.1%	5.1%	5.1%	5.1%	5.1%	33.6%

ECMO CONTROL	Interval Ending (Days Post Implant)				
	0	1	7	14	30
# Left	48	41	12	3	0
Total # Failed	0	5	10	15	16
Survival	100%	89.4%	69.6%	42.2%	0%
Std Error	0%	4.5%	8.9%	12.2%	.

Fig. 3-3 Cohort 2 Survival

Because the Kaplan-Meier analysis censors subjects at time of transplant, “Competing Outcomes” curves were constructed to show a more complete picture of the endpoints.

The following figure shows the “Competing Outcomes” for Cohort 1. The curves represent each of the outcomes and at any time point the sum of the proportions of outcomes equals 100%.

Of the 24 Cohort 1 subjects, 21 were transplanted between 1 to 174 days of support. The 2 deaths in this Cohort occurred at 0 and 38 days post implant. One subject was weaned after 146 days due to poor prognosis.

Competing Outcomes – Cohort 1




Fig. 3-4 Cohort 1 Competing outcomes

The next figure shows the “Competing Outcomes” for the ECMO control group for Cohort 1. The longest support time was 20.5 days at which time 75% were weaned from ECMO for recovery or transplant.

Competing Outcomes – ECMO Control group for Cohort 1



Fig. 3-5 Cohort 1 control group competing outcomes

The following figure shows the “Competing Outcomes” for Cohort 2. Of the 24 Cohort 2 subjects, 21 were transplanted between 3 to 192 days of support. The 2 deaths in this Cohort occurred at 19 and 144 days post implant. One subject was successfully weaned to recovery after 9 days.

Competing Outcomes – Cohort 2




Fig. 3-6 Cohort 2 competing outcomes

The next figure shows the “Competing Outcomes” for the ECMO control group for Cohort 2. The longest support time was 27.5 days at which time 67% were weaned from ECMO for recovery or transplant.

Competing Outcomes – ECMO Control group for Cohort 2




Fig. 3-7 Cohort 2 Control Group Competing Outcomes

a) Secondary Efficacy Results

There were two secondary efficacy objectives of the study. The first was to summarize the days of transplant eligible support.

Only one subject was removed from the transplantation listing at any point during their support. The subject (in Cohort 2) was first listed on day 3 of support (10/03/09) and then was delisted from 01/15/10 to 02/22/10 due to a neurological event. The subject was successfully transplanted on 04/10/10. The summary statistics of time of eligible support are detailed in the following table.

Days of Transplant Eligible Support

Cohort	N	Median	Mean ± Std	Range
Cohort 1	24	27.5	58.8 ± 56.1	0 – 174
Cohort 2	24	42.5	55.6 ± 44.3	3 – 151

Tab. 3-9 Days of transplant eligible support

The second objective was to show the ability to de-intensify concomitant hemodynamic support. At each visit, the subject’s status was recorded with the following choices: sedated, intubated, on ECMO, awake, ambulating or eating. The following table summarizes those choices pre-implant, and at 2 weeks and 1 month post-implant. A subject could have more than one status subcategory checked.

Prior to implant, 22 of the 24 Cohort 1 subjects (92%) and 16 of 24 Cohort 2 subjects (67%)

were sedated and/or intubated and over 30% were supported by ECMO immediately prior to device implant.

In Cohort 1 there were 7 subjects (7/20=35%) who were sedated and intubated at 2 weeks with 1 sedated and awake (1/20=5%). The other 12 (12/20=60%) were awake with some of those also ambulating and eating.

In Cohort 2, 6 subjects (6/20=30%) were still sedated and intubated at 2 weeks with 1 awake and intubated (1/20=5%) and the remaining 13 awake (13/20=65%). At 1 month post, those numbers drop to only 3 of the Cohort 1 and 4 of the Cohort 2 subjects remaining sedated and intubated.

Support Status at each Follow-up Visit

Time Point	Status (more than 1 could be checked)	Cohort 1 n=24	Cohort 2 n=24
Pre-implant N=24 In each cohort	Sedated	21 (87.5%)	16 (66.7%)
	Intubated	21 (87.5%)	14 (58.3%)
	On ECMO/other	8 (33.3%)	9 (37.5%)
	Awake	3 (12.5%)	12 (50.0%)
	Ambulating	0 (0.0%)	5 (20.8%)
	Eating	0 (0.0%)	8 (33.3%)
2 Weeks N=20 In each cohort	Sedated	8 (40.0%)	6 (30.0%)
	Intubated	7 (35.0%)	6 (30.0%)
	Awake	13 (65.0%)	14 (70.0%)
	Ambulating	3 (15.0%)	4 (20.0%)
	Eating	6 (30.0%)	12 (60.0%)
1 Month N=12 Cohort 1 N=17 Cohort 2	Sedated	4 (33.3%)	5 (29.4%)
	Intubated	3 (25.0%)	5 (29.4%)
	Awake	9 (75.0%)	13 (76.5%)
	Ambulating	3 (25.0%)	8 (47.1%)
	Eating	4 (33.3%)	9 (52.9%)

Tab. 3-10 Support status at each follow-up visit

3.6.6.2 Primary Safety

The total time on device of the Cohort 1 subjects was 1411 days. There were 96 serious adverse events (SAEs) for this cohort yielding a rate of **0.068 events per patient-day**. The 95% Poisson confidence interval was calculated as: [0.055, 0.083]. The total time on device for Cohort 2 was 1376 days. There were 109 SAEs for this cohort yielding a rate of **0.079 events per patient-day** with the confidence interval as [0.065, 0.096]. A summary of SAEs rates for each cohort is included in the first table of this clinical study section.

a) Infection Serious Adverse Events

Major Infection events were reported according to the Investigational Plan definition (which is the same as the INTERMACS definition). Any time an additional medication was added for treating a different organism a new SAE was reported (or adjudicated as an event). The study

design was intentionally broad with regard to setting a low threshold for calling an event an infection. Fever was defined at 38 degrees, WBC > 15,000, positive cultures from any source, or decision to start antibiotics with or without positive cultures were listed as an SAE and subsequently adjudicated. Each infection was counted as a separate event even when occurring concurrently in one patient, ensuring that the infection rate would not be under-reported.

In Cohort 1, 15 subjects had 35 total infectious events reported. In Cohort 1, a majority of subjects had pre-existing risks for infection including ventilation (83%), pre-implant ECMO support (33%), and previous cardiac surgery (21%).

In the larger subjects (Cohorts 2) there were fewer events (12 subjects with 24 events) which is as expected based on age and body size.

Outcomes of any of the subjects did not appear to be affected by infections as the deaths that occurred were not solely related to infection, even when one was present. These cases tended to have multi-factorial contributors such as stroke, end-organ failure, arrhythmias, or thromboembolism. All other subjects with a noted infectious SAE were transplanted or weaned. Infection had little impact on the transplant wait time since 99.3% of the total time the subjects were on support was considered transplant eligible time.

b) Major Bleeding Serious Adverse Events

Major Bleeding was the third most frequently reported SAE in Cohort 1 (10 subjects with at least one event). All bleeding events for Cohort 1 occurred in subjects less than 2 years old. Five of the 10 subjects in Cohort 1 with bleeding events were younger than 9 months old. Young infants have some degree of ineffective erythropoiesis. Hemoglobin subsequently falls to a nadir at around 2–3 months of age due to decreased RBC production. Anemia in acute or critical illness may be exacerbated by numerous factors including blood loss (due to hemorrhage or sampling), reduced RBC production (due to nutritional deficits, inflammatory processes or low erythropoietin levels) and increased RBC turnover due to hemolysis.

Cohort 1 subjects had a pre-implant history of transfusion in 92% (22/24), history of ECMO or previous VAD in 33% (8/24), and 21% (5/24) of subjects had previous cardiac surgeries. These factors along with the strict Major Bleeding definition could have contributed to the percentage of events reported.

Major Bleeding was one of most prevalent events in Cohort 2 with 12 of 24 (50%) subjects experiencing a bleeding event.

c) Hypertension Serious Adverse Events

Hypertension was reported per the protocol definition (consistent with the INTERMACS definition). An event was logged each time a subject's blood pressure reached the 95th percentile for age and was treated with an IV agent. Several hypertension events were reported in the early post-op periods. However, 75% (15/20) of the hypertension events were in Cohort 1 and 2 subjects who only received LVAD support. This is not surprising as it is common for patients supported only with left sided devices to require pharmacological support in order to optimize right ventricular function with agents that can cause hypertension, resulting in the concomitant need for agents to lower the blood pressure in the early post-operative period. Additionally, hypertension is one of the leading post operative cardiac surgical events for children, especially the younger children, possibly due to their reactive vasculature. In order to follow the event definition, hypertension events were reported when the values met the definition even if the subject was also on a pressor or in a period where the site was trying to optimize the overall hemodynamic status of the subject in the early post-op period. There did not appear to be a correlation between Hypertension and Major Bleeding.

d) Neurological Dysfunction Serious Adverse Events

Four of the 48 (8.3%) Cohort 1 and 2 subjects experienced a neurological dysfunction with long term severe results (PSOM scores ≥ 2) and another 2 (4.2%) were withdrawn from support due to the neurological injury.

In Cohort 1, 7 of the 24 subjects experienced a neurological event. One subject experienced 2

ischemic events. Of the 7 subjects, 1 was withdrawn from support as a result of the neurological injury. Of the remaining 6 subjects, PSOM exams were performed post explant and 1 had no deficit (assessed 17 days post explant); 2 had mild deficits (23 and 221 days post explant), 1 had moderate deficit (82 days post) and 2 had severe deficits (PSOM score of 3 at 34 days post and score 4 at 54 days post).

In Cohort 2, 7 of the 24 subjects experienced a neurological event. Two of those subjects experienced both an ischemic and hemorrhagic event. Of the 7 subjects, 1 was withdrawn from support as a result of the neurological injury. Of the remaining 6 subjects, PSOM exams were performed post explant and 1 had no deficit (50 days post explant); 2 had mild deficits (27 and 49 days post explant), 1 had moderate deficit (357 days post) and 2 had severe deficits (PSOM scores of 10 at 29 and 38 days post).

This table summarizes the status information.

Summary of Neurological Event Status

Long term Result	Cohort 1 N=24	Cohort 2 N=24	Total N=48
No Deficit (PSOM 0.0)	1 (4.2%)	1 (4.2%)	2 (4.2%)
Mild (PSOM 0.5-1.0)	2 (8.3%)	2 (8.3%)	4 (8.3%)
Moderate (PSOM 1.5-2.0)	1 (4.2%)	1 (4.2%)	2 (4.2%)
Severe (PSOM ≥ 2.5)	2 (8.3%)	2 (8.3%)	4 (8.3%)
Support withdrawn	1 (4.2%)	1 (4.2%)	2 (4.2%)
TOTAL	7 (29.2%)	7 (29.2%)	14 (29.2%)

Tab. 3-11 Summary of neurological event status

Pump Replacement Due to Thrombus

During the course of the support, a clinician may have identified that a pump required replacement due to visualized thrombus within the blood pump. These replacements were not considered adverse events. However, these were nonetheless regarded as sentinel events due to their frequency and association with thromboemboli.

In the primary cohorts, 24 (50%) of the subjects had at least one pump replacement due to suspected thrombus (11 Cohort 1, 13 Cohort 2). The number of pump replacements ranged from 0 to 4 per subject. The average number of replacements per subject was 0.9 ± 1.2 . However, subjects were supported on the device for varying lengths of time therefore it may be more informative to consider the replacements per length of time on device. The average replacements-per-day on device was 0.02 ± 0.03 per day.

At the IDE sites, 57 (52.3%) of the 109 subjects had at least one pump replacement due to thrombus (11 Cohort 1, 14 Cohort 1 CAP, 13 Cohort 2, and 19 Cohort 3). The number of pump replacements ranged from 0 to 6 per subject. The average number of replacements per subject was 1.1 ± 1.4 and the average replacements-per-day on device was 0.02 ± 0.03 per day.

Additionally, 95 subjects were enrolled at non-IDE sites. Of the 204 subjects, 93 (45.6%) subjects had at least one pump replacement due to thrombus (11 Cohort 1, 14 Cohort 1 CAP, 13 Cohort 2, and 19 Cohort 3, 36 Cohort 3 Non-IDE). The number of pump replacements ranged from 0 to 6 per subject. The average number of replacements per subject was 1.1 ± 1.4 and the average replacements-per-day on device was 0.02 ± 0.03 per day.

Cohort	N	# Subjects with at least 1 replacement	Total number of replacements	Replacements per Subject	Total Days on Device	Replacements per Days on Support	Time to first replacement (days)
primary Cohorts *	48	25 (50.0%)	43	0.9 ± 1.2 0 - 4	2787	0.02 ± 0.03 0.00 - 0.13	24.1 ± 19.7 4 - 105
IDE Cohorts	109	57 (52.3%)	114	1.1 ± 1.4 0 - 6	6350	0.02 ± 0.03 0.00 - 0.18	19.1 ± 16.9 2 - 105
Non-IDE Cohorts	95	36 (37.9%)	58	0.6 ± 1.0 0 - 4	7240	0.01 ± 0.03 0.00 - 0.27	41.9 ± 44.6 2 - 198
Total	204	93 (45.6%)	172	0.8 ± 1.2 0 - 6	13590	0.02 ± 0.03 0.00 - 0.27	27.8 ± 32.3 2 - 198

* Note: the 48 subjects in the “Primary Cohorts” group are a subset of the “IDE Cohorts” group (n=109)

Tab. 3-12 Pump replacement

3.6.6.3 Death information

Two subjects in each of the primary cohorts died after support was withdrawn. The 4 subjects were supported a median time of 28.5 days ranging from 0 to 144 days (mean ± std: 50.3 ± 64.4 days). Of the 4 subjects who died, 75% (3/4) were supported with ECMO at the time of EXCOR implant.

The CEC reviewed all deaths at the IDE sites and assigned primary and secondary causes of death. These causes are summarized by subject in the following table.

Patient	Days on Device	Primary Cause	Secondary Cause(s)
COHORT 1 (2 deaths/ 24 subjects)			
#1	0	Pulmonary Respiratory Failure	Cardiovascular: Left A-V valve regurgitation
#2	38	CNS: Multiple ischemic strokes	None
COHORT 2 (2 deaths/ 24 subjects)			
#3	144	Other: Arterial CNS and non-CNS Thromboembolism	Infection
#4	19	CNS: Large ischemic strokes with hemorrhagic conversion	Other: Tonsillar herniation

Tab. 3-13 Primary and secondary cause of death

3.6.7 Conclusion

Despite the reported SAEs, 42 of the 48 subjects supported by the EXCOR were adequately supported to transplant and 1 subject was able to be weaned successfully from the device after 9 days of support yielding an 89.6% success rate (43/48). The device supported children safely

to cardiac transplantation for a median transplant eligible time of 27.5 and 42.5 days for cohort 1 and 2 respectively. Only one subject was temporarily removed from transplant eligibility during their support and was eventually relisted and transplanted.

Data that strongly supports the consideration for probable benefit is summarized for both Cohort 1 and 2 subjects as shown in the following tables.

Probable Benefit

Cohort	N	Outcome				Success (Transplant or Weaned- Recovered)
		Transplant	Weaned- Recovered	Weaned- Failure	Died	
Cohort 1	24	21	0	1	2	21/24 (87.5%)
Cohort 2	24	21	1	0	2	22/24 (91.7%)
Total	48	42	1	1	4	43/48 (89.6%)

Tab. 3-14 Probable Benefit

Post-Explant/Transplant Follow-up

Cohort	N	Outcome	30 days post-explant		1 year post-explant	
		# Explanted	# (%) alive 30 days	Lost to Follow-up	# (%) alive 1 Year	Lost to Follow-up
Cohort 1	24	22	22/22 (100%)	n/a	17/22 (77%)	0
Cohort 2	24	22	21/22 (95%)	1*	16/17 (94%)**	1
Total	48	44	43/44 (97.7%)	1	33/39 (85%)	1

* 1 subject was weaned and returned to home

** 5 subjects have regular contact with the site for post transplant care but are not 1 year post-explant as of this report: 3 subjects are due in June (last report alive at 313, 257 and 250 days), 1 subject is due in July (last report alive at 170 days) – verbal report; denominator includes 1 LTF

Tab. 3-15 Post-explant/transplant status follow up

Beyond the primary endpoint of survival to transplant, the majority of subjects remain alive at 1 year post-explant/transplant as noted in the previous table.

HDE regulations require the device under study to show **reasonable safety and probable benefit**. In the EXCOR® Pediatric IDE trial the device demonstrated probable benefit as a bridge to transplantation in patients who are transplant eligible with severe left ventricular or biventricular dysfunction. The majority of patients implanted with the EXCOR were transplant eligible during device support with adequate end organ function and decreasing need for hemodynamic support such as intubation, sedation or ECMO support. While the concomitant support decreased, the subjects were able to spend more time awake, eating and ambulating.

The benefits offered to subjects implanted with the EXCOR® Pediatric include additional time to await transplant and improved hemodynamics allowing removal of pre-implant hemodynamic support allowing for increase time awake, ambulating and eating contributing to post implant transplant eligible wait times. These far-reaching benefits outweigh the risks associated with the adverse events that occurred.

This page was left blank intentionally

4 Description: blood pump, cannulae and accessories

EXCOR is an extracorporeal electro-pneumatically driven ventricular assist device. It can be used for either univentricular or biventricular support. EXCOR is comprised of the following permanently active components:

- extracorporeal blood pump(s)
- inflow and outflow cannula(e)
- 1 driving tube for each blood pump
- *Ikus*




Fig. 4-1 EXCOR shown in situ as a biventricular assist device in pediatric application

Overview

The blood flows from the atrium or the ventricle through the inflow cannula into the blood chamber of the pump and then from this blood chamber through the outflow cannula into the aorta or into the pulmonary artery. A driving tube is used to connect the air chamber of the pump to the electro-pneumatic *Stationary Driving Unit Ikus*. *Ikus* generates the suction and driving pressures required to move the triple-layer membrane separating the blood chamber from the air chamber.

4.1 EXCOR blood pumps




Fig. 4-2 60 ml blood pump

EXCOR blood pumps have a transparent polyurethane (PU) housing which is divided into an air chamber and a blood chamber by a triple-layer membrane.

The blood chamber has an inflow and an outflow stub to which the inflow and outflow cannula, respectively, are connected. The pump stubs themselves are made of polyurethane, the end of each stub is fitted with a titanium connector to which the cannula will be connected. The valves located in the pump stubs keep the blood flowing in one direction. EXCOR blood pumps are available with three-leaflet valves made of polyurethane (10 - 60 ml stroke volume).

All surfaces of the pump coming into contact with the blood are coated with a *Carmeda® BioActive Surface (CBAS®)* coating. The transparent casing of the blood pump allows easy visual monitoring of the filling and emptying of the blood chamber.

The blood pump is equipped with a de-airing nipple which is used for de-airing the blood chamber when the pump is being commissioned.

The air chamber of the pump is equipped with a driving tube connector. This connector is used to connect the blood pump to the driving tube through which air is pumped from the *Ikus*. *Ikus* generates the suction and driving pressures required to move the blood pump's triple-layer membrane. A graphite powder layer is located between the membrane layers in order to minimize friction.

4.2 EXCOR cannulae

3 different types of cannulae are available for EXCOR in various sizes for each type:

- atrial cannulae (as inflow cannulae)
- LV apex cannulae (as inflow cannulae)
- arterial cannulae (as outflow cannulae)

The cannulae are made of tissue-friendly silicone. Polyester-velour suture rings enable convenient and safe anastomosis of the cannulae. The mid section of all cannulae is covered with polyester-velour in order to promote good ingrowth of the cannulae where they pass through the skin.

Some vascular cannulae have a shaping wire which allows the cannulae to be adapted to each individual patient's anatomic conditions.



Fig. 4-3 Cannula heads: 1) atrial cannula, 2) LV apex cannula, 3) arterial cannula

4.3 EXCOR accessories

The following *EXCOR* accessories are required in order to commission and operate *EXCOR*:

- 1 driving tube (PVC) for each blood pump
- 2 tank units
- 1 accessory (T00L-002) set which includes:
 - membrane set
 - de-airing set (2 x trocar, 2 x de-airing tube)
 - de-airing hammer
 - tube connecting set (cable ties, cable-tie gun)

There is enough material in 1 accessory set (T00L-002) to commission 2 *EXCOR* blood pumps.

This page was left blank intentionally

5 Implantation: Preparations in the operating room

NOTE: This chapter omits safety instructions, information and procedures that refer to the *Ikus* exclusively. Please refer also to the IFU.

5.1 Preparing the components and materials required

NOTICE

Selection of blood pump(s): see section 12.1: Overview: Product range and possible combinations, page 111.

ADVICE

It is advantageous to provide a sterile table on which to place the prepared sterile components.

General (all sterile)

- 500 ml sterile injectable saline
- 2 small sterile basins
- 50 ml disposable syringe with luer lock connector
- suture (to secure the trocar to the de-airing nipple and the de-airing tube to the trocar)
- heavy scissors
- towel clamp, tube clamp
- other instruments and equipment as required for open-heart surgery

EXCOR components and accessories

- blood pump(s), each with a pump seal
- 1 driving tube for each blood pump
 - univentricular: driving tube, red
 - biventricular: 1 red driving tube and 1 blue driving tube
- inflow cannula(e) (atrial or LV apex cannula)
- outflow cannula(e)
- accessory set (T00L-002) for blood pumps with PU valves
 - membrane set
 - de-airing set (2 x trocar, 2 x de-airing tube)
 - de-airing hammer
 - tube connecting set (cable ties, cable-tie gun)

5.2 Checking and adjusting the settings of the cable tie gun


Before using the cable tie gun contained in the EXCOR *Tube connecting set* the accuracy of settings has to be checked and if necessary to be corrected.



Fig. 5-1 Cable tie gun

► INSTRUCTION

1. Check if the following values are set:
 - coarse adjustment on STD (2)
 - fine adjustment on 5 (1)



2. In the case of deviations loosen the screw (4) and disassemble the locking cap (3).



3. Adjust the above-mentioned values with the adjusting wheels (6 and 5). Begin with adjusting wheel 6.



4. Assemble the locking cap (3) and secure it with the screw (4).



5.3 Unpacking the sterile components

WARNING

Only use sterile components which have been delivered in undamaged sterile condition (sterile packaging intact, expiration date not expired).

Only use blood pumps which have an undamaged aluminum-coated outer packaging.

INSTRUCTION

1. Pump: a non-sterile person opens the aluminum-coated package and removes the pump in its double sterile packaging.
2. The non-sterile person opens the outer sterile package.
3. A sterile person takes out the inner sterile package, opens it and places the components on the prepared sterile field.

5.4 Moving the membrane to the end-of-diastole position



- a de-airing nipple (blood chamber)
- b driving tube connector (air chamber)

Fig. 5-2 De-airing nipple and driving tube connector


INSTRUCTION

1. Pick up adapter tube, disposable syringe (membrane set) and the pump.
2. Connect the adapter tube to the disposable syringe.
3. Connect the free end of the adapter tube to the driving tube connector of the blood pump.
4. Remove all air from the air chamber of the pump. The blood pump membrane is now in the end-of-diastole position.
5. Seal the adapter tube with a tube clamp in order to keep the membrane in this position.

5.5 De-airing the blood pump

Prepare and place the following ready for use:

- blood pump(s) with pump seal(s)
- 1 de-airing set (trocar and a de-airing tube) for each blood pump
- 50 ml disposable syringe for each blood pump



- a trocar (de-airing needle with obturator in place)
- b de-airing nipple
- c blood chamber
- d membrane in end-of-diastole position
- e air chamber

Fig. 5-3 Pump with trocar in place (de-airing needle with inserted obturator)

5.5.1 Inserting the de-airing needle

WARNING

The membrane must be kept in the end-of-diastole position. Keep the clamped membrane set connected to the blood pump.

INSTRUCTION

1. Take hold of the trocar (de-airing needle with obturator) and remove the protective silicone cap.

De-airing the blood pump

2. Push the trocar as pictured above as far as it will go through the center of the blood pump's de-airing nipple. Never turn the trocar when inserting it, this increases the risk of removing a large piece of the silicone material in the de-airing nipple.
3. Remove the obturator.
4. Withdraw the de-airing needle by approx. 2 mm. Important: The tip of the cannula should still be visible in the blood chamber.
5. Use the suture to fix the de-airing needle to the de-airing nipple.
6. Remove the adapter tube from the pump.

5.5.2 Rinsing and filling the blood pump

ADVICE

Before commencing surgery, mark the points for the exit sites of the cannulae. The aim is to achieve a stable final position of the cannulae without exerting any tension on the skin. Caution: with biventricular support, 2 of the 4 cannulae will cross each other. This crossing point should be outside of the thorax as far as possible.

INSTRUCTION

1. Fill and empty the pump once or twice with sterile injectable saline.
2. Push the free end of the de-airing tube onto the trocar as far as it will go. Secure the de-airing tube to the trocar with a suture tie.
3. Fill the syringe with sterile injectable saline.
4. Connect the syringe to the stopcock end of the de-airing tube.
5. Slowly fill the pump with sterile injectable saline. Rock the pump back and forth to move any bubbles to the outflow stub.
6. Close the stopcock on the de-airing tube.
7. Tap the blood pump body gently in order to free all remaining bubbles. Remove all air from the pump through the outflow connector.
8. Use the seal caps to close the titanium cannula connectors.
9. Place the pump ready for connection with the connectors pointing up.

6 Implantation - surgical procedure

This chapter describes the product-specific measures to be observed when implanting an EXCOR blood pump.

NOTE: This chapter omits safety instructions, information and procedures that refer to the *Ikus* exclusively. Please refer also to the IFU.

Unless any specific instructions to the contrary are given, the same protocol as for any other major cardiothoracic surgical procedure should be followed. Implantation is accomplished using a CPB with bicaval cannulation. Implantation can be achieved with induced ventricular fibrillation or on a beating heart, hypothermia is usually not required.

WARNING

After implantation each cannulae and all connections must be inspected for it's solidity, safeness and tightness.

Do not start pump operation until the blood pump is completely free of air!

Do not touch or manipulate the blood pump with pointed or sharp-edged objects (e. g. surgical instruments)!

If a cannula is bent with flexible metal reinforcement to adjust it to the anatomical conditions: determine by visual inspection that the blood flow in the cannula is not restricted.

NOTICE

For the suture use an appropriate suture material. It should be a nonabsorbable monofilament, not traumatizing material.

ADVICE

For BVAD, carry out anastomosis of the cannulae in the following order:

- apical cannulation
 1. LV apex
 2. right atrium
 3. pulmonary artery
 4. aorta
 - atrial cannulation
 1. left atrium
 2. right atrium
 3. pulmonary artery
 4. aorta
-

6.1 Cannula exit sites




Fig. 6-1 Cannula position following implantation




Fig. 6-2 Suggested cannulae exit sites (Example: BVAD with LV apex cannulation)

6.2 Use of the cannula tunneling tip

The cannula tunneling tip is a sterile disposable product and is supplied with each cannula. Sizes available: see figure Fig. 6-3: Available sizes of cannula tunneling tips, page 55. Staged cannulae are supplied with 2 different tunneling tips.

► INSTRUCTION

1. Push the cannula tunneling tip firmly into the distal end of the cannula.
2. Advance the forceps through the subcostal incision and the cannula tunnel into the

- mediastinum, so that the cannula tunneling tip can be gripped.
3. Use the forceps to firmly grip the flat end piece, pull it through the cannula tunnel and the skin incision and position it.
 4. Carefully remove the tunneling tip from the cannula by bending it back and forth.

Refer to the respective cannula type as described in sections 6.3 to 6.6 of the instruction for use to determine the sequence of cannulae anastomosis and tunneling.




Fig. 6-3 Available sizes of cannula tunneling tips




Fig. 6-4 Use of cannula tunneling tip

6.3 Cannulae and connector set

To avoid damages of cannulae careful attention should be paid to the following safety precautions.

⚠ WARNING

During implantation the *Cannula Tunneling Tip* (provided with each cannula) should be used during implantation of the *EXCOR* system.

If it is necessary to apply a clamp directly to the cannula in order to pull the cannula through the skin, the following procedures should be observed:

- Position the clamp at the distal end of the cannula
 - After the cannula has been pulled through the skin, cut off and discard the part of the cannula where the clamp was applied.
-

If it is necessary to clamp any other part of the cannula that is not covered with velour, cover the part of the cannula that will be clamped with a gauze sponge.

If the connector set is being used: Secure each connection between blood pump and cannula with at least one cable tie.

If an *EXCOR connecting set* is required for implantation and the length of the tube part needs to be reduced, the tube part should be cut but only to achieve the following minimum lengths:

Part Number	Diameter Reduction	Minimum Length
A12-016	16 to 12 mm	90 mm
A09-012	12 to 9 mm	75 mm
A06-009	9 to 6 mm	60 mm

Tab. 6-1 Connector set: minimum length of connector tube

⚠ WARNING

If replacement of an *EXCOR* blood pump is required, the following procedures should be observed:

- The cable tie covering the *EXCOR* cannula on the stub of the blood pump should be removed carefully. Important: never use a sharp instrument, for example, a scalpel or scissors, to remove the cable tie. This may cause damage to the cannula.
- If a connecting set needs to be cut for a pump replacement, ensure that there will be sufficient length of the tube part remaining to meet the minimum length recommendations.

Do not kink the drivelines. Otherwise there might not be sufficient pump output

Do not kink the cannulae needlessly. Otherwise there might not be sufficient pump output. Moreover, cannulae might be damaged.

At least daily, the *EXCOR* cannulae should be inspected for signs of wear or damage. ADVICE: To avoid needless kinking of the cannulae use a mirror for inspection of the bottom side of the blood pump.

In no case should the cannulae either be kinked directly at the connector to the blood pump or at the transition area between velour and silicone.

6.4 Access

➤ INSTRUCTION

1. Median sternotomy. Make sure that there is absolutely no bleeding.
2. Insert standard cardiopulmonary bypass cannulae (bicaval cannulation).
3. Initiate extracorporeal circulation.
4. Place a vent in the left atrium, if necessary.

6.5 LV apex cannula

Refer to section 6.2: Use of the cannula tunneling tip, page 54.

6.5.1 Anastomosis of inflow cannula with LV apex

⚠ WARNING

During anastomosis of the LV apex cannula, make sure that the cannula head is facing in the right direction: the long side of the head should be parallel to the lateral wall. This prevents the ventricular lateral wall from being sucked into the tip of the cannula. After the cannula head has been placed, its position can be checked by means of the flow direction arrow on the cannula body (except LV apex cannulae C10A-030, C14A-040, C18A-020). The arrow is aligned with the long side of the cannula head (see figure Fig. 6-6: Ideal position of the LV apex cannula, page 57).

➤ INSTRUCTION

1. If indicated, initiate ventricular fibrillation as needed.
2. Apical excision of the LV: The ideal implant position of the LV cannula is slightly off-center of the LV apex toward the lateral wall. The distance from LAD/ septum to the center of the excised muscle core is about 2 cm for children.
3. We recommend to excise a circular apical core with a diameter slightly smaller than the size of the cannula head.
4. Start with muscle core incision on the side away from the septum/ LAD (see b in figure Fig. 6-6, page 57) to avoid septal injury.
5. Check left ventricle for thrombi and excise the excess trabeculae.




Fig. 6-5 Anastomosis of LV apex cannula




Fig. 6-6 Ideal position of the LV apex cannula



a Long side of LV apex cannula head

Fig. 6-7 Head of LV apex cannula

6.5.2 Creating a transcutaneous tunnel for the LV apex cannula

⚠ WARNING

Always use the cannula tunneling tip provided (see section 6.2: Use of the cannula tunneling tip, page 54) to advance the cannula through the prepared transcutaneous tunnel. Never use a sharp surgical instrument directly on the cannula.

Make sure that the blood pump and cannulae come to rest in a stable position without tension.

Do not touch or manipulate the silicone cannulae with pointed or sharp-edged objects (e. g. surgical instruments).

If it is necessary to apply a clamp directly to the cannula in order to pull the cannula through the skin, the following procedures should be observed:

- Position the clamp at the distal end of the cannula
- After the cannula has been pulled through the skin, cut off and discard the part of the cannula where the clamp was applied.
- If it is necessary to clamp any other part of the cannula that is not covered with velour, cover the part of the cannula that will be clamped with a gauze sponge.

The skin incision must be slightly smaller than the cannula diameter (to ensure good ingrowth) but large enough to prevent necrosis.

Plan the cannula exit sites appropriately. Leave an adequate bridge of skin and subcutaneous tissue between the cannula exit incisions to prevent breakdown and necrosis of the skin and tissue. If possible, the cannula exit sites should be on different planes (see fig. Fig. 6-2, page 54).

➤ INSTRUCTION

1. Prepare the transcutaneous tunnel. Ensure that the incision is large enough.
2. Incise the pericardium widely in a lateral direction. Prepare the cannula tunnel by blunt dissection. Important: Do not tunnel transperitoneally.
3. Tunnel the LV apex cannula through the transcutaneous passage by using a pair of forceps to firmly grip the flat end piece of the tunneling tip and pull it through the cannula tunnel and the skin incision.
Important: Do not rotate the cannula while pulling it through the tunnel. At the end of this procedure, the apex of the heart should be in its native position without torsion.
4. Terminate ventricular fibrillation if necessary.

6.6 Atrial cannula(e)

Refer to section 6.2: Use of the cannula tunneling tip, page 54.

ADVICE

For atrial cannulae supplied with a forming wire, the transcutaneous tunnel should be created and the cannula advanced through the tunnel and skin incision prior to the anastomosis.
For all other atrial cannulae, the sequence is arbitrary.

6.6.1 Creating a transcutaneous tunnel for atrial cannula(e)

WARNING

If possible, always use the cannula tunneling tip provided (see section 6.2: Use of the cannula tunneling tip, page 54) to advance the cannula through the prepared transcutaneous tunnel.

If it is necessary to apply a clamp directly to the cannula in order to pull the cannula through the skin, the following procedures should be observed:

- Position the clamp at the distal end of the cannula
- After the cannula has been pulled through the skin, cut off and discard the part of the cannula where the clamp was applied.
- If it is necessary to clamp any other part of the cannula that is not covered with velour, cover the part of the cannula that will be clamped with a gauze sponge.

Care must be taken to ensure that the cannulae come to rest in a stable position free of tension.

Do not touch or manipulate the silicone cannulae with pointed or sharp-edged objects (e. g. surgical instruments).

Using a pair of forceps, firmly grip the flat end piece of the tunneling tip and pull it through the cannula tunnel and the skin incision.
Important: Do not rotate the cannula while pulling it through the tunnel.

The incision must be slightly smaller than the cannula diameter (to ensure good ingrowth) but large enough to prevent necrosis.

Plan the cannula exit sites appropriately. Leave an adequate bridge of skin and subcutaneous tissue between the cannula exit incisions to prevent breakdown and necrosis of the skin and tissue. If possible the cannula exit incisions should be on different planes.

INSTRUCTION

1. Prepare the transcutaneous tunnel. Ensure that the incision is large enough.
2. Prepare the cannula tunnel by blunt dissection. Important: Do not tunnel transperitoneally.
3. Using a pair of dressing forceps, tunnel the cannula through the transcutaneous tunnel. Important: Do not rotate the cannula while pulling it through the tunnel.

6.6.2 Anastomosis of atrial cannulae

Right atrium



Create the anastomosis laterally, directly above the tricuspid valve.

a) closed technique




Fig. 6-8 Cannulation of right atrium




Fig. 6-9 Suture technique, right atrium

INSTRUCTION

1. Make a running (purse-string) suture with monofilament, secured with pledgets at 4 positions.
2. Place 4 single U-sutures secured with pledgets on each side of the purse string suture.
3. Make a sufficiently long incision inside of the suture circle and extend it as required.
4. Push the cannula down on the sutures, at the same time slightly reduce the venous inflow to the CPB while inflating the lung in order to prevent negative pressure in the left atrium.
5. Remove all air from the cannula and use a tube clamp to clamp the cannula below the anastomosis.

b) open technique with bicaval cannulation

With bicaval cannulation, the right atrial cannula can be inserted in an open technique.

Left atrium

The procedure for anastomosis of the left atrium corresponds to the procedure applied to the right atrium.



Place anastomosis at the junction of the right upper pulmonary vein and the left atrium. The atrial wall is the recommended implantation location. The pulmonary vein should be left intact.

6.7 Arterial cannula(e)

Refer to section 6.2: Use of the cannula tunneling tip, page 54.



For cannulae supplied with a forming wire, the transcutaneous tunnel should be created and the cannula advanced through the tunnel and skin incision prior to the anastomosis.

6.7.1 Creating a transcutaneous tunnel for arterial cannula

⚠ WARNING

Care must be taken to ensure that the blood pump and cannulae come to rest in a stable position.

Do not touch or manipulate the silicone cannulae with pointed or sharp-edged objects (e. g. surgical instruments).

Using a pair of forceps, firmly grip the flat end piece of the tunneling tip and pull it through the cannula tunnel and the skin incision. Important: Do not rotate the cannula while pulling it through the tunnel.

The incision must be smaller than the cannula diameter (to ensure good ingrowth) but large enough to prevent skin necrosis.

Plan the cannula exit sites appropriately. Leave an adequate bridge of skin and subcutaneous tissue between the cannula exit incisions to prevent breakdown and necrosis of the skin and tissue. If possible the cannula exit incisions should be on different planes (see fig. Fig. 6-2, page 54).

➤ INSTRUCTION

1. Prepare the transcutaneous tunnel. Ensure that the incision is large enough.
2. Prepare cannula tunnel by blunt dissection. Important: Do not tunnel transperitoneally.
3. Using a pair of forceps, firmly grip the flat end piece of the tunneling tip and pull it through the cannula tunnel and the skin incision. Important: Do not rotate the cannula while pulling it through the tunnel.

6.7.2 Anastomosis of the arterial cannula

Aorta

➤ INSTRUCTION

1. Tangentially clamp the ascending aorta and make a longitudinal opening of a length which is suitable for the cannula diameter. If necessary, offset the incision laterally to the right by up to 45°.
2. Anastomose the cannula using ten teflon-backed double-reinforced individual monofilament (e. g. 4-0 EB) U-sutures. (If simpler conditions are encountered, a running suture can be made instead.)
3. Remove all air from the cannula and use a tube clamp to clamp the cannula below the anastomotic site. If it is necessary to clamp any other part of the cannula that is not covered with velour, cover the part of the cannula that will be clamped with a gauze sponge.




Fig. 6-10 Anastomosis of the aortic cannula

Pulmonary artery

INSTRUCTION

1. Make a longitudinal incision of a size suitable for the cannula diameter in the pulmonary artery.
2. Anastomose the cannula using 10 teflon-backed, double-reinforced individual monofilament (e. g. 4-0 EB) U-sutures. (If simpler conditions are encountered, a running suture can be made instead.)
3. Remove all air from the cannula and use a tube clamp to close it below the anastomosis. If it is necessary to clamp any other part of the cannula that is not covered with velour, cover the part of the cannula that will be clamped with a gauze sponge




Fig. 6-11 Cannulation of the pulmonary artery

6.8 Shortening the cannulae if necessary

WARNING

If an *EXCOR connecting set* is required for implantation and the length of the tube part needs to be reduced, the tube should be cut but only to achieve the following minimum lengths:

Part Number	Diameter Reduction	Minimum Length
A12-016	16 to 12 mm	90 mm
A09-012	12 to 9 mm	75 mm
A06-009	9 to 6 mm	60 mm

Tab. 6-2 Connector set: minimum length of connector tube

INSTRUCTION

1. Cut the cannulae to the required length. Make the cut perpendicular to the cannula axis and ensure that the cut is straight.
2. Make sure that the lengths of the 2 cannulae leading to the same pump match. It must be possible to connect the cannulae to the pump without having to exert any tension.

6.9 Connecting the blood pumps to the cannulae

WARNING

Ensure that cannulae, blood pump(s) and driving tubes are not subject to external forces and are free of kinks or sharp bends.

When connecting the blood pump(s), pay attention to the direction of the arrows on the inflow and outflow stubs. These show the direction of the blood flow.

Type of support	Anastomosis of inflow cannula to	Points upwards...
Univentricular		
LVAD	apex	blood chamber
LVAD	atrium	air chamber
Biventricular		
LVAD	apex	blood chamber
LVAD	atrium	air chamber
RVAD	atrium	air chamber

Tab. 6-1 Anastomosis and direction of the blood chambers

NOTICE

Finally, the driving tube is connected to the *Ikus*. The *Ikus* is started and the parameters are gradually adjusted (see section 6.10.4: Checking the parameters when the pump is started and adjusting them, page 66).

INSTRUCTION

1. Bring the patient into the Trendelenburg position.
2. Release the tube clamps, flush the cannulae and then use tube clamps to clamp the cannulae below the exit sites. If it is necessary to clamp any other part of the cannula that is not covered with velour, cover the part of the cannula that will be clamped with a gauze sponge.
3. First connect the inflow cannula to the pump, then connect the outflow cannula. When doing so, add sterile injectable saline with a bulb syringe in order to connect the pump air free. Be careful to avoid damaging the gloves and the inner cannula (lumen) and pump surfaces.
4. Release the tube clamps, de-air the pump(s) and the cannulae.
5. Connect the driving tube to the blood pump. Biventricular: use the red driving tube for the left blood pump and the blue driving tube for the right blood pump. Univentricular: always use the red driving tube.




Fig. 6-12 Final position of the blood pumps, for example: BVAD with LV apex cannulation

6.10 Intraoperative drive management

NOTE: This section omits safety instructions, information and procedures that refer to the *Ikus* exclusively. Please refer also to the IFU.

6.10.1 Connecting the blood pump(s) to the *Ikus*

WARNING

Do not kink either the driving tubes or the cannulae.

NOTICE

State of the blood pumps when they are initially connected: filled with sterile injectable saline, de-airing needle in place. To allow easier handling, the driving tubes are not connected until the inflow and outflow cannulae have been connected to the pump (see section 6.9: Connecting the blood pumps to the cannulae, page 63).

INSTRUCTION

1. Open the driving tube connector marked in red (univentricular) or both connectors (biventricular). To do so, pull the seal plugs out of the connector(s).
2. Connect the driving tube to the *Ikus*. To do so, push the plug of the driving tube into the connector. The sound of the plug snapping into place is clearly audible. Check that the plug is securely connected. To do so, grip the plug body above the release sleeve and pull on it. Do not pull from the release sleeve, and never from the tube!
3. In biventricular mode: observe the color of the markings.
4. In biventricular mode: repeat the procedure for the second pump.

Operating mode	<i>Ikus</i> connector
biventricular	LVAD: connector marked red RVAD: connector marked blue
univentricular	connector marked red

Tab. 6-2 Assignment: operating mode, blood pump, connector


6.10.2 De-airing the blood pumps in single-step mode

NOTICE

Each de-airing step (**Step left** / **Step right**) carries out half a pump cycle (systole or diastole), the 1th step being a diastole. Normally, several de-airing steps are required for each pump. In single-step mode, the pumps will operate using the pressures shown in the parameter table. It will not be possible to switch to the standard view unless at least 1 de-airing step has been completed for each connected pump.

INSTRUCTION

1. Bring the patient into the Trendelenburg position.
2. Move the cursor to the field marked **Step left**.
3. Lift the pump. The de-airing nipple is the highest point.
4. To trigger a single step, press the <Enter> key. If necessary, use the de-airing needle to vent the air from the pump (see section 5.5: De-airing the blood pump, page 51). After consulting the surgeon: If necessary, press <Enter> repeatedly to trigger further single steps until all air has been removed from the pump(s). If the the blood pump is not filling sufficiently, ensure there is sufficient preload and if necessary, increase the diastolic pressure.
5. In biventricular mode: Move the cursor to the field **Step right**. Repeat the procedure for the 2nd pump.



Transition to continuous pumping mode is not yet possible because **OK** field is still inactive.

Fig. 6-13 Single-step mode

6.10.3 Starting the blood pump (changing to standard view)

WARNING

Do not start the pump(s) until all air has been removed.

Once the de-airing needle has been removed it cannot be re-inserted.

Only remove the de-airing needle after all air has been removed from the blood pump, the blood pump is running and the parameters have been adjusted (see section 6.10.4: Checking the parameters when the pump is started and adjusting them, page 66 and section 6.11: Removing the de-airing needle, page 69).

INSTRUCTION

1. Move cursor to the **OK** field and press <Enter> to confirm. The system now starts with the parameter values visible in the parameter table.

6.10.4 Checking the parameters when the pump is started and adjusting them

⚠ WARNING

In order to avoid air being sucked into the blood pump through the cannula anastomosis, adjust the parameters gradually. If air does enter the system, disconnect the driving tubes from the *Ikus* and de-air the system using the de-airing needle.

Continuously monitor all settings.

Once the de-airing needle has been removed it cannot be re-inserted.

NOTICE

If the pump is not filling adequately at this stage, increase the pre-load by adding volume from the CPB circuit. After adding volume, adjust the parameters on the laptop of the *Ikus* as described in the following table.

➤ INSTRUCTION

1. Observe the left blood pump. Is the pump ejecting completely? If not: increase the left driving pressure if necessary.
2. Observe the right blood pump. Is the pump ejecting completely? If not: increase the right driving pressure if necessary.

Observe	Action / measure
Right pump Is the pump filling properly? (see below)	If not: check the filling pressure (central venous pressure; CVP) CVP too low: substitute volume CVP too high: increase suction pressure If no improvement occurs: check the position of the cannulae via echographic monitoring!
Left pump Is the pump ejecting properly?	If not: check mean arterial pressure (Guideline value: 70mmHg)
Compare left and right pump. Is left pump filling considerably worse than right pump?	If yes: increase suction pressure on left side If no improvement occurs: check the position of the cannulae via echographic monitoring!

Tab. 6-3 Pump filling criteria

Keep the following points in mind with regard to filling of the right pump:

The aim is to reduce the right ventricle’s load to a large extent but not completely. Signs that the RV load has been reduced completely are:

- filling of the pump depends largely on the respiratory cycle
- ventricle is empty/limp
- membrane stops abruptly during filling

Important: If the three above-mentioned phenomena are observed, do one of the following:

- reduce the diastolic pressure
- substitute volume

Adjusting parameters

➤ INSTRUCTION

1. Use the <←>/<→> keys to move the cursor to the desired field in the parameter table. The selected field is given a colored background.

- Use the <↓>,<↑> or <Bild-↓>, <Bild-↑> keys to adjust the value, then press <Enter> to confirm the input.

Parameter	Range possible	<↓>/<↑> changes value by	<Bild-↓>/<Bild-↑> changes value by
Systolic pressure [mmHg]; driving pressure	60 to 350	2.5	25
Diastolic pressure [mmHg]; suction pressure	0 to -100	2.5	25
Rate [bpm]	30 to 150	1	10
Relative systolic duration [%]	20 to 70	1	10

Tab. 6-4 Parameter's possible adjustments

In biventricular operation: adjusting the operating mode

To run the pumps in the asynchronous mode or separate mode instead of the synchronous mode the appropriate mode must be selected.

- asynchronous mode is recommended for patients who have a small thorax volume in comparison to the pump volume. In asynchronous mode, the intrathoracic blood volume remains unchanged.
- separate mode is useful, under some circumstances, for patients with intracardiac shunts.

INSTRUCTION

- Use the <<->/<->> keys to move the cursor to the field showing the current operating mode. A pop-up menu showing the available operating modes is opened (see table Tab. 6-2: Assignment: operating mode, blood pump, connector, page 64).
- Select the desired operating mode with <↓>,<↑> and confirm with <Enter>. The system will now work in the selected mode.

Guideline values

The most important criteria when selecting drive parameters is that they ensure a good filling and emptying of the pump; the parameters must be set to achieve this goal.

NOTICE

The systolic driving pressure must be higher than the patient's physical systolic pressure. Important: If the systolic duration (% systole) is reduced or if very small cannulae are used, it may be necessary in some cases to select a higher value than recommended here.

The actual driving pressures achieved are influenced by the diameter of the cannulae used.

The following values are merely guideline values; they may not be appropriate in each individual case

Systolic pressure [mmHg], left/ right	Diastolic pressure [mmHg], left/ right	Rate [bpm]	Rel. systolic duration [%], left/ right
220/150	-40/-40	80	40/40

Tab. 6-5 Recommended guideline values for normal operation

ADVICE

Remove the de-airing needle after all air has been removed from the blood pump, the blood pump is running and the parameters have been adjusted (see section 6.10.4: Checking the parameters when the pump is started and adjusting them, page 66 and section 6.11: Removing the de-airing needle, page 69).
 Important: Once the de-airing needle has been removed it cannot be re-inserted.

6.10.5 Switching from CPB support to VAD support

The aim here is to reduce the CPB flow and in doing so to shift the volume from the CPB to the patient (i.e. to the VAD).

WARNING

Secure the driving tubes and cannulae to the blood pump(s) as soon as the proper function of the EXCOR is established (see section 6.12: Securing the connections, page 70).

NOTICE

When using staged cannulae or a connecting set, do not set a pumping rate > 100 bpm, as the pump will not eject its full volume at higher rates. With these cannulae rates > 100 bpm are to be avoided.

INSTRUCTION

1. When the blood pump(s) starts to fill, reduce the CPB flow and gradually increase the EXCOR rate from an initial 30 bpm until CPB has been terminated and the required flow is achieved. Important: In doing so, make sure that the pump fills adequately, and if necessary regulate the driving pressure.
2. If necessary, adjust the systolic pressure, diastolic pressure and the systolic percent.

6.10.6 Possible complications

Decreased filling after stable filling conditions

If a good filling behavior was achieved at first (filling pressures LA/CVP < 10 mmHg and diastolic pulmonary artery pressure < 15 mmHg) with good drainage and nominal rate (normally 80 bpm), but the filling has deteriorated over time, it usually will not help to increase the diastolic pressure.

Deterioration in the filling behavior despite stable inflow conditions may indicate hypovolemia or obstruction of the inflow cannula. The cause of deterioration in filling behavior must be identified and addressed.

NOTICE

Manipulations during implantation can severely influence the inflow temporarily – wait for the situation to stabilize before adjusting the values.

INSTRUCTION

1. Evaluate volume status and transfuse if necessary. Evaluate and if necessary correct the cannula position.

Pump filling deteriorates when thorax is closed

If atrial cannulation is used, a slight decrease in the filling may be observed in some cases when the thorax is closed. This may be caused by compression of the atria or a slight shift in the position of the cannulae.

INSTRUCTION

1. Evaluate volume status and transfuse if necessary. Important: Observe the effect volume replacement on the pump filling!
2. Increase suction pressure.

Distinct decrease in filling or generally poor inflow conditions on right side

INSTRUCTION

1. Make sure that there is no inflow obstruction.
2. If a suction pressure of less than -50 mmHg is necessary, increase the relative diastolic duration as an additional measure. At the same time, reduce the relative systolic duration. Important: Increase the driving pressure accordingly!

Incomplete ejection right/left

INSTRUCTION

1. Observe the arterial blood pressure, and at the same time observe the ejection movement of the pump membrane.
2. If complete emptying of the pump is no longer achieved, adjust the driving pressure accordingly. Important: Do not respond to extreme – temporary – increases in the arterial blood pressure (due to manipulation, catecholamine, etc.).

6.11 Removing the de-airing needle

WARNING

When removing the de-airing needle, never pull on the de-airing tube, but on the de-airing needle itself.

Before removing the de-airing needle, be sure that the de-airing tube is secured to the de-airing needle. Important: Once the de-airing needle has been removed it cannot be re-inserted.

NOTICE

Do not remove the de-airing needle until all air is removed, the blood pump is running, all parameters have been adjusted and the chest has been closed. (see section 6.10.4: Checking the parameters when the pump is started and adjusting them, page 66).

INSTRUCTION

1. Cut the suture material between the de-airing needle and the de-airing nipple (see image 1 in figure Fig. 6-14, page 70). Important: Leave the ligature around the de-airing nipple (see image 2 in figure Fig. 6-14, page 70).
2. Pull the de-airing needle out of the de-airing nipple.




Fig. 6-14 Removing the de-airing needle

After the patient has been weaned from the CPB and the proper function of the EXCOR is established, the connections of the driving tubes and cannulae to the blood pump(s) have to be secured.

6.12 Securing the connections

WARNING



All connections have to be secured by at least 1 cable tie. 2 cable ties may be used. Exception: connection between drive line and drive line connector of the blood pump: 1 cable tie only!

INSTRUCTION

1. Pick up the Tube connecting set.
2. Secure the following connections:
 - inflow cannula on the connector
 - outflow cannula on the connector
 - drive line on the drive line connector (1 cable tie only!)

3. The 1. cable tie must be positioned exactly on the groove profile of the connector (1).
Important: the heads of the cable ties have to be directed away from the patient's body.
4. Fasten the cable ties by the cable tie gun. Important: pay attention to 5.2: Checking and adjusting the settings of the cable tie gun, page 49.

5. A 2nd cable tie can be used optionally. If a 2nd cable tie shall be used (2) it has to be positioned above the 1st cable tie. IMPORTANT: the heads of the cable tie straps should both be staggered and directed away from the patient's body.



6. If an EXCOR Connecting set is required for implantation after that secure also those connections with cable ties. Proceed thereby as described in the instruction steps 3 to 5.

6.13 Postoperative drive management

NOTICE

The patient should receive the same treatment as is usual after any other major cardiac surgical procedure.

6.13.1 After transfer to the ward

If a good filling and stable ejection of the blood pump(s) is observed in the immediate post-operative period, it is normally not necessary to adjust the driving and suction pressures.

- Good filling means that the suction pressure is adequate.
- Stable ejection (at normal arterial blood pressure) means that the driving pressure is adequate.

WARNING

At least every 4 hours, visually check that the pump(s) is (are) filling and ejecting completely over a period of several pump cycles. If a pump is not filling and/ or ejecting completely, appropriate measures are to be taken.

NOTICE

For further details on regular monitoring of pump(s) and cannulae, see section 8.5: Regular checks of blood pump(s) and cannulae, page 79.

6.13.2 Follow-up treatment

Guideline values and criteria for adjusting the parameter settings: see table Tab. 6-5: Recommended guideline values for normal operation, page 68.

It is only necessary to adjust the left driving pressure when

- the arterial blood pressure increases (e. g. after lifting sedation, when the patient wakes up)
- when the patient is mobilized (moving to an upright position, sitting, standing – in order to compensate for the additional hydrostatic pressure component).

This page was left blank intentionally

7 Implantation - anesthesia

The following risk factors should be closely monitored for anesthetic and hemodynamic management:

- right heart function during LVAD implantation
- coagulopathy
- renal insufficiency
- abnormal reactions to inotrope administration
- pulmonary hypertension

CAUTION

There should be an adequate supply of pre-matched stored blood, fresh frozen plasma and platelet concentrates available for immediate transfusion if required.

Keep blood product transfusions to a minimum. Blood transfusions may lead to the development of antibodies, which are known to promote coagulation and inflammatory response.

ADVICE

Medication for right ventricular afterload reduction should be available for use in the operating room (nitric oxide NO, phosphodiesterase inhibitor, prostaglandin, etc)

Auto-transfusion equipment (e. g. Cellsaver) should be available for use in the operating room.

For patients with an LVAD, start ventilation with nitric oxide or administer the appropriate medication to treat pulmonary hypertension and reduce afterload for right ventricle 15 minutes before weaning from the CPB. This can help to prevent or lower the risk of right ventricular failure.

Monitoring procedure

Intraoperative monitoring should include the same monitoring procedures applied during major cardiothoracic surgery:

- central venous line
- Swan-Ganz catheter (if appropriate)
- arterial line
- ECG
- pulseoximetry
- central temperature monitor
- urine catheter

Additional recommended monitoring procedures

- cardiac output calculation (if appropriate)
- intraoperative transesophageal echocardiogram (inflow cannula position, heart valve function, intracardial shunts, volume status)
- right heart function in case of LVAD

Any other monitoring processes can be used (e. g. neurological monitoring) at the anesthesiologist's discretion.

This page was left blank intentionally

8 Wound care and treatment

Cannula exit sites should be treated like open wounds. The patient's wounds should always be attended to by a small group of nurses in the inpatient area.

The only way to ensure there is a minimum risk of infection is to provide good wound care.

WARNING

Before cleaning the wound (see 8.3: Cleaning of the wound, page 77), put on sterile disposable gloves, cap and mask.

Cleaning the pump and the drive line: Do not use any acetone or petroleum based products near the pump or drivelines. We recommend using only water or alcohol to clean the pump and the drive line.

IMPORTANT: Do not use any corrosive or colored solutions or organic solvents to clean the blood pump or the drive line as they may alter the surface of the product.

Cleaning the cannulae and transcutaneous exit site: Do not use any acetone or petroleum based products near the cannulae and the transcutaneous exit site.

We recommend using chlorhexidine to clean the cannulae and transcutaneous exit site.

IMPORTANT: Do not use any corrosive or colored solutions or organic solvents to clean the cannulae and the transcutaneous exit site as they may alter the surface of the product.

NOTICE

Do not stick bandages to the cannulae. Over time, remnants of adhesive contaminate the cannulae and increase the risk of infection.

Do not use any adhesive on the velour coating of the cannula as it is difficult to remove and may adversely manipulate the cannula.

Do not use organic solvents near the EXCOR Pediatric such as petroleum ether or turpentine oil, as they could damage the cannulae and the pumps. The plastic parts must not get in contact with chlorinated hydrocarbon (e.g. chloroform), thinners (e.g. acetone, naphtha, toluol, xylene, heptane) or similar compounds.

Do not mark or write on the plastic parts.

Material required (with biventricular access):

- Sterile dressing tray
- Disinfectant i.e. 2% chlorhexidine solution
- Clean gloves
- Mask
- Sterile gloves and towel
- *Metalline*[®] drain compress
- 2X2 gauze, 4X4 gauze
- Adhesive dressing (i.e. *Mepore*[®])
- Adhesive remover
- Non sting barrier film sticks
- Abdominal pads
- Tape
- Tubular bandage (i.e. *Burnnet*)



Fig. 8-1 Materials for dressing change

How often to change the dressing

If the wound is dry and not infected:

- POD 1- once a day
- POD 11-28 every second day, if the wound is dry and not infected
- POD >28 twice a week, if the wound is dry and not infected

If the wound shows signs of infection: clean wound and change dressing twice a day

8.1 Removing the old dressings

➤ INSTRUCTION

1. Unpack all the material required to dress the wound and place this within reach on a sterile sheet.
2. Put on disposable gloves, remove old dressings.
3. Take off the disposable gloves, put on the sterile gloves.
4. Remove old dressing using no-touch technique.
5. Examine the places where the cannulae pass through the skin and if changes are apparent take appropriate measures if necessary.
6. Use adhesive remover to remove any adhesive dressing.
Important: adhesive remover (depending on contents) might damage cannula and the pump, use only on skin.

8.2 Cleaning the blood pump



Fig. 8-2 Cleaning the blood pump

INSTRUCTION

1. Cleanse the exposed cannula and the pump head with disinfectant (i.e. 2% chlorhexidine solution) then place on sterile towel.
2. Observe cannulae and cannulae exit sites.
3. Remove gloves.

8.3 Cleaning of the wound

INSTRUCTION

1. Hand hygiene, prepare sterile dressing tray, put on sterile gloves. If assistance is necessary notify Berlin Heart.
2. 4X4 gauze soaked in 2% chlorhexidine cleanse each cannula exit site in a circular motion outward to a radius of approximately 10 cm.
3. Using a new soaked 4X4 repeat 2 more times beginning at the exit site and clean in larger circles each time.



Fig. 8-3 Cleanse each cannula exit site

4. Wrap 4X4 gauze soaked in 2% chlorhexidine around cannula and gently cleanse with back/forth motion.
5. Repeat with each cannula exit site.
6. Cleanse entire cannula (upper and bottom side).
7. 4X4 gauze soaked in 2% chlorhexidine solution.
8. Starting at the exit site moving down cannula approximately 10 cm from exit site.
9. Repeat for each cannula exit site.
10. Allow chlorhexidine to dry completely.



Fig. 8-4 Cleanse with back/forth motion



Fig. 8-5 Cleanse entire cannula

8.4 The new dressing

8.4.1 Preparing a new dressing

INSTRUCTION

1. Apply non sting barrier film to skin around cannulae. Non sting barrier prevents skin maceration around cannula exit sites.



Fig. 8-6 Non sting barrier film

8.4.2 Applying a new dressing

INSTRUCTION

1. Wrap a Metalline drain compress around each cannula (from right to left, slit always facing upwards (see figure Fig. 8-7).
2. Attach the Metalline drain compresses above the cannulae using sterile bandages. First secure the outer compresses, then the inner compresses (see figure Fig. 8-8).



Fig. 8-7 Metalline drain compress



Fig. 8-8 Secure with a sterile bandage

3. Pass a gauze compress folded lengthwise beneath the 2 left cannulae. The open end of the folded compress should point in the direction of the wound. Pull the cannulae into place by tugging the compress slightly (see figure Fig. 8-9).
4. Fold the left end of the compress upwards, diagonally to the right and secure with a sterile bandage (see figure Fig. 8-9).
5. Fold the right end of the compress upwards, diagonally to the left and secure with a sterile bandage (see figure Fig. 8-11).



Fig. 8-9 Gauze compress under the cannulae



Fig. 8-10 Fold the left end of compress and secure



Fig. 8-11 Fold the right end of compress and secure

6. Repeat this procedure for the 2 right cannulae. In this way, the 4 cannulae are padded so that they do not press on the skin or wound (see figure Fig. 8-12).

7. Cover the entire wound broadly with gauze compresses (see fig. Fig. 8-13).



Fig. 8-12 Cannulae are padded



Fig. 8-13 Cover with sterile gauze compresses

8. Secure the upper part of the dressing with a sterile bandage (see figure Fig. 8-14).
9. Finally, seal the dressing at the left and right side, below the cannulae and between the individual cannulae with strips of adhesive bandage (e. g. Leukoplast), see figure Fig. 8-15.



Fig. 8-14 Secure with a sterile bandage



Fig. 8-15 Seal with strips of adhesive bandage

10. Place tubular bandage (i.e. *Burnnet*) around patient (see figure Fig. 8-16).
11. Tie in front to secure dressing.



Fig. 8-16 Tubular bandage

8.5 Regular checks of blood pump(s) and cannulae

Frequency of inspection: every 4 hours



Everyone involved in caring for an EXCOR patient must be trained to carry out a visual check, to evaluate the filling behavior of the blood pump(s) and to detect deposits.



At least daily, the EXCOR cannulae should be inspected for signs of wear or damage. ADVICE: To avoid needless kinking of the cannulae use a mirror for inspection of the bottom side of the blood pump.

At least every 4 hours, check visually that the blood pump(s) is (are) filling and ejecting completely over a period of several pump cycles. If a pump is not filling and/ or ejecting completely, then take the appropriate corrective action.

Under certain circumstances, the message **left/right pump is not filling adequately** in some circumstances is not generated with the 10 ml EXCOR blood pump due to the low volume of air which is moved in the pump. Therefore in pumps of this size, pay special attention to the movement of the membrane and ensure that each pump fills and empties completely.

8.5.1 Visual inspection: pump filling and ejection

The filling and ejection behavior of a blood pump is optimal when the membrane surface is completely smooth at the end-of-systole and end-of-diastole positions. Check visually that the pump(s) is (are) filling and ejecting completely over a period of several pump cycles. If a pump is not filling and/ or ejecting completely, take the appropriate corrective action.

Cautionary measures

For all blood pumps: check the position and condition of the driving tube and the cannulae (inflow deterioration due to kinks in cannulae/driving tubes is rather rare).

For all blood pumps: check the membrane movement.

Medical examination of patient

Check CVP, mean arterial pressure and adjust therapy if necessary.

Check the volume status:

- amount of bleeding
- increased urine output (use of diuretics?)
- tamponade
- *Important:* Increasing the suction pressure will not bring about any distinct improvement if there is not sufficient volume available.

LVAD: observe the functions of the right ventricle.

Adjusting the parameter values

Only adjust the parameters if the measures listed above have no effect or in case of:

- *Mobilization of patient:* adjust the systolic pressure, both left and right. When pressures have increased, do not reduce these again, even when the patient is lying down.
- *Signs of low cardiac output:* the membrane is moving properly while at the same time a decrease in urine output, lactate increase and dyspnea (shortage of breath) can be observed. In this case, increase the rate and adjust other settings as required.

INSTRUCTION

1. Use the <<->/<->> keys to move the cursor to the desired field in the parameter table. The selected field is given a colored background.
2. Use the <↓>,<↑> and <Bild↓>, <Bild↑> keys to adjust the value, then press <Enter> to confirm the input. The system will now operate using the new settings.

Cautionary measure

Confirm each changed parameter value by pressing <Enter>. The system does not take over the new, changed value until it has been confirmed with <Enter>.


ADVICE

Enter all the changes to the parameter values into the parameter log. (see section 12.4: Sample copy: EXCOR parameter log, page 120).

8.5.2 Visual inspection: deposits

Check the blood pump(s) for visible deposits (fibrin, clots) every 4 hours. If deposits develop, check the pump(s) every hour.

Checking the pump areas which come in contact with blood



- 1 transition inflow cannula – inflow connector
- 2 inflow stub in front of inflow valve (only for pumps with PU valves)
- 3 inflow valve
- 4 inflow stub behind inflow valve
- 5 area between inflow and outflow stubs
- 6 remaining area of blood chamber
- 7 transition blood chamber - membrane (directly above the reinforcement ring)
- 8 outflow stub in front of outflow valve
- 9 outflow valve
- 10 (only for pumps with PU valves) outflow stub behind outflow valve
- 11 transition outflow connector – outflow cannula

Fig. 8-17 Diagram of EXCOR blood pump (top view of blood chamber)

ADVICE

During the visual check, first clean the blood pump then illuminate the blood chamber with a flashlight. This makes it easier to detect deposits. Enter all of the findings into the blood pump log. (see section 12.3: Sample copy: EXCOR pump log, page 118).

Cautionary measures

Initial signs of deposits: check anticoagulation therapy and adjust therapy if necessary.
 Floating deposits inside the pump: replace the pump!

8.5.3 Checks using the monitor program

Record all drive parameters and adjust if necessary.

Objective: the blood pump(s) must fill and eject completely in each pumping cycle, the diastolic pressure should be as low as possible.

ADVICE

Record the parameter values once a day.

To record the parameters use the sample copy in section 12.4: Sample copy: EXCOR parameter log, page 120.

8.5.4 Replacing the blood pump due to growth of the patient



In children, plan to replace the pump(s) with a larger pump(s) in good time, to prevent the possibility of inadequate support due to an insufficient discharge rate.

The pump selected at the time of transplantation may not be adequate for the entire period of cardiac support. Growth and/or weight gain can result in the patient not receiving adequate support. Use the chart in section 12.1.2: Overview: Relationship: body weight – pump size, page 111, to plan, in good time, which pump(s) the patient may need to change over to. This chart is for guideline purposes only and is not binding for each individual case. This decision must be taken by the surgeon in consultation with Berlin Heart GmbH.



Notify Berlin Heart! 866.249.0128

The blood pump(s) must be replaced as described in section 10.1: Replacing the blood pump(s), page 89.

9 Anticoagulation therapy

9.1 Before Implantation of the EXCOR

9.1.1 General considerations

Patients with an EXCOR system must be maintained on anticoagulation therapy.

Anti-Xa levels should be specific to the drug being used, either unfractionated heparin or enoxaparin.

The TEG[®] may be useful in managing unfractionated heparin and antiplatelet therapy. Please contact *Berlin Heart*, Clinical Affairs for further information.

9.1.2 Pre implantation

The following laboratory tests should be considered prior to implantation.

- Platelet Function Studies, INR, PTT, fibrinogen, antithrombin III, and platelet count to establish a baseline. Assessment for thrombophilia by measuring Protein C, S, Factor V Leiden, Prothrombin 20210 defect, as well as Heparin Induced Thrombocytopenia (HIT) is recommended.

9.2 During Implantation - Cardiopulmonary Bypass

9.2.1 Cardiopulmonary Bypass (CPB)

Use unfractionated heparin as per institutional protocol for cardiopulmonary bypass.

9.2.2 Post CPB

Completely reverse heparin with protamine sulphate as per institutional protocol.

The goal post-CPB is to achieve normal (institution specific) coagulation parameters (INR, PTT, fibrinogen, platelet count).

In the early post-operative period, the possibility of surgical bleeding, GI bleeding, internal bleeding in the retro-peritoneum or other bleeding diathesis is possible and must be monitored.

If the patient is bleeding despite normal coagulation parameters consider:

- Von Willebrand's
- Surgical bleeding

9.3 Postoperative anticoagulation therapy

9.3.1 General Considerations

Primary tests used to evaluate anticoagulation in the patient include antifactor Xa levels and/or PTT.

9.3.2 Starting anticoagulation therapy

During the first 24 hours following implantation, no anticoagulants should be administered.

Approximately 24 - 48 hours after implantation, commence unfractionated heparin therapy (i.v.) if the following criteria are met:

- Platelet count >20,000/ μ l
- Normal Platelet Function Studies
- Minimal bleeding in infants and young children.

9.3.3 Unfractionated heparin therapy (i.v.) Patient < 12 months

- Initial dose 15 IU/kg/hour.
- Do not use a bolus
- After 6 hours if the patient does not have increased bleeding, increase the heparin infusion to 28 IU/kg/hour (therapeutic dose).

6 hours after increasing the heparin to the therapeutic dose, obtain a PTT and an antifactor Xa level.

If the anti factor Xa level is desired range (0.35-0.5 U/ml) and the PTT is in the therapeutic range (institution dependent), then either the PTT or anti factor Xa level may be used to follow the heparin therapy.

If the anti factor Xa level is <0.35 U/ml or >0.5 U/ml, increase or decrease the heparin infusion, respectively until the anti factor Xa level is the therapeutic range (see Tab. 9-1, page 87).

Anti factor Xa levels should be obtained daily. Important: hyperbilirubinemia may result in falsely low anti factor Xa levels. If anti Xa levels do not correlate with the PTT in this setting, consider using the PTT to monitor heparin therapy.

Antithrombin should be >70%. If the antithrombin is <70%, treat according to institutional protocol.

9.3.4 Unfractionated heparin therapy (i.v.) Patient ≥ 12 months

Initial dose 10 IU/kg/hour.

Do not use a bolus.

After 6 hours if the patient does not have increased bleeding, increase the heparin infusion to 20 IU/kg/hour (therapeutic dose).

6 hours after increasing the heparin to the therapeutic dose, obtain a PTT and an anti factor Xa level.

If the anti factor Xa level is desired range (0.35-0.5 U/ml) and the PTT is in the therapeutic range (institution dependent), then either the PTT or anti factor Xa level may be used to follow the heparin therapy.

If the anti factor Xa level is < 0.35 U/ml or > 0.5 U/ml, increase or decrease the heparin infusion, respectively until the anti factor Xa level is the therapeutic range (see Tab. 9-1, page 87).

Anti factor Xa levels should be obtained daily. Important: hyperbilirubinemia may result in falsely low anti factor Xa levels. If anti Xa levels do not correlate with PTT in this setting, consider using the PTT to monitor heparin therapy.

Antithrombin should be >70%. If the antithrombin is <70%, treat according to institutional protocol.

NOTICE	<p>If during standard unfractionated heparin therapy:</p> <ol style="list-style-type: none">1. Platelet count is < 40,000/μl revert to the Stage I heparin dose for continuous infusion (see Tab. 9-1, page 87)2. Platelets <20,000/ul discontinue heparin and consider evaluation for heparin induced thrombocytopenia (HIT). <hr/> <p>If the anti factor Xa or PTT is too low or too high during heparin therapy, never use a bolus of heparin or protamine. Instead, increase or decrease the heparin dose, IU/hour, as required (see Tab. 9-1, page 87).</p>
---------------	---

9.3.5 Thrombelastography (TEG®)

TEG® analysis may be useful in managing the anticoagulation and anti-platelet therapy. Please contact *Berlin Heart Inc.*, Clinical Affairs for further information.

9.4 Low Molecular Weight Heparin

At 48 hours following surgery if all bleeding has stopped, the creatinine is within normal limits, and the patient is hemodynamically stable, switching from unfractionated heparin to low molecular weight heparin (LMWH) is recommended.

- Patient < 3 months start administration of Enoxaparin at 1.5 mg/kg subcutaneously every 12 hours.
- Patient > 3 months start administration of Enoxaparin at 1 mg/kg subcutaneously every 12 hours.
- Stop heparin infusion and administer LMWH (subcutaneously) simultaneously.
- Obtain the first anti factor Xa level at 4 hours after the 2nd LMWH dose is administered. See Tab. 9-2, page 88 for monitoring and dosing.
- Anti factor Xa therapeutic range: 0.6 to 1.0 U/ml.
- Anti factor Xa should be monitored along with platelet count, and creatinine
- When using LMWH, monitor Anti factor Xa daily. Once the Anti Factor Xa level is in the therapeutic range at a stable dose, monitor twice a week for 2 weeks, and then weekly.

9.5 Oral Anticoagulation Therapy (only for patients \geq 12 months of age who are taking a full oral diet)

ADVICE

This section only applies to patients \geq 12 months. Oral anticoagulation in children < 12 months of age is not recommended due to difficulties with monitoring the warfarin effect.

When the patient's condition has been fully stabilized (e.g. hemodynamically stable, no evidence of bleeding, etc), switch to oral anti-coagulation therapy with a vitamin K antagonist (target INR: 2.7 to 3.5), with an initial loading dose of 0.2 mg/kg/day. Do not exceed maximum loading dose of 5mg/day. The INR must be checked daily in the first 4 weeks, twice a week for the next 4 weeks (if INR is stable), and once a week there after (see Tab. 9-3, page 88 and Tab. 9-4, page 88).

Until the target INR is achieved, simultaneous administration of warfarin and heparin is necessary (approximately 4 days). Once the target INR is achieved, heparin therapy can be discontinued. If the INR decreases to < 2.7, administer LMW heparin immediately and then q12h until an INR of \geq 2.7 is achieved. (Table 5, Appendix 2) If INR is 2.0- 2.7 use an enoxaparin dose of 0.5 mg/kg targeting an anti factor Xa level of 0.3-0.5, if INR is <2.0 use an enoxaparin dose of 1 mg/kg targeting an anti factor Xa level of 0.5 - 1.0.

When unable to achieve a stable INR with warfarin, LMWH should be used instead. Discontinue the warfarin and administer LMWH as per previously discussed age related dosing (see Tab. 9-2, page 88).

9.6 Monitoring of Blood Count and Anticoagulation Status

Monitoring the anticoagulation status as well as infection risk, and renal and hepatic function is important and should be monitored with the following frequency:

- Daily while on UFH, twice a week while on enoxaparin/coumadin for 4 weeks then once week : Fibrinogen, D-dimer, aPTT, PT/INR, Platelet Count, TEG[®], Antithrombin, WBC, HgB, HCT, BUN/SCr, AST/ALT, bilirubin T/D, prealbumin, CRP.
- While on UFH obtain anti factor Xa level daily.
- While on enoxaparin obtain anti factor Xa daily until in therapeutic range and on a stable dose, then twice a week for two weeks and then weekly.

If infection is suspected, appropriate measures must be taken immediately (antibiotic therapy, adjustment of the anticoagulation and platelet inhibition therapy) and increased monitoring of the coagulation system. In addition, in the setting of hemodynamic instability, organ dysfunction, and inadequate anticoagulation daily monitoring should be performed until any of

these issues are resolved

9.7 Postoperative platelet inhibition therapy

As individual patient responses vary to the anti-platelet agents, the optimum dosage for each patient will be that which minimizes both the risk of thromboembolic complications when the dose is too low and the risk of hemorrhagic complications when the dose is too high. Acetylsalicylic acid (ASA) and dipyridamole are the anti-platelet agents recommended.

9.7.1 Start of therapy

Dipyridamole

At 48 hours after surgery, start dipyridamole, 4mg/kg/day p.o. divided into 4 doses (1 mg/kg Q6) (maximum dose 15mg/kg/day). If the following are present:

- All bleeding has stopped, AND
- The patient is hemodynamically stable AND,
- Platelet studies do not show significantly decreased function,
- Platelet count is > 40,000/ μ l,

Acetylsalicylic Acid

At 4 days post implantation, following the removal of all drainage tubes, start acetylsalicylic acid (ASA) 1mg/kg/day p.o., divided into 2 doses (0.5 mg/kg Q 12), if the following are present:

- Platelet studies show platelet inhibition in the presence of AA < 70%

The ASA dose should split and be administered two times daily (0.5 mg/kg Q 12) due to the short half life and the high turnover of the platelets (approximately 10 % new platelets per day).

9.8 Adjunctive Medication

The inflammation parameters (Tissue factor pathway inhibitor, prothrombin fragment 1-2, fibrinogen, Factor VIII) for patients on ventricular assist device support are often elevated above normal. Accordingly, the physician may choose to administer the following medications at his/her discretion to facilitate the overall anticoagulation/anti-platelet management of the patient:

- Omega-3 fatty acids (e.g. DHA/EPA), have been shown to have an anti-inflammatory effect and also decrease premature activation of platelet membrane. Omega-3-fatty acids are composed of long chain polyunsaturated long chain carbons. Only alpha-linolenic acid (ALA) of the omega-3 family is truly essential.

Antioxidants (Vitamin C and E) also have been shown to have an anti-inflammatory effect, and may be considered.

9.9 Anticoagulation Therapy

9.9.1 Therapeutic Heparin administration and adjustment

NOTICE

This table assumes the site therapeutic PTT is 60 to 85 seconds (Monagle, P, et al.). Each site should use their hospital calculated therapeutic range.

Stage	Description	Anti factor Xa [u/ml]/PTT	Infusion	Hold heparin	Rate Change [%]	Repeat PTT
I	Initial Dose (first 6 hours)					
	Infant < 12 mo		15 IU/kg/h			
	Child ≥12mo		10 IU/kg/h			
II	Therapeutic Dose					
	Infant < 12 mo		28 IU/kg/h			after 6h
	Child ≥12mo		20 IU/kg/h			after 6h
III	Adjustment					
		<0.1/<50	0	0	+15%	4h
		0.1-0.34/50-60	0	0	+10%	6h
		0.35-0.50/60-85	0	0	0	next day
		0.51-0.70/86-95			-10%	6h
		0.71-0.89/96-120		30 min.	-10%	4h
		= 0.90/ >120		60 min.	-15%	4h

Tab. 9-1 Unfractionated Heparin adjusted to maintain an anti factor Xa level of 0.35 to 0.50 U/ml.

Anti Factor Xa level U/ml?	Hold Next Dose?	Dose Change?	Repeat Anti Factor Xa?
< 0.35	no	increase dose by 25%	4h after next dose
0.36 - 0.45	no	increase dose by 15%	4h after next dose
0.46 - 0.59	no	increase dose by 10%	4h after next dose
0.6 - 1.0	no	no	4h after next dose
1.1 - 1.25	no	decrease dose by 20%	4h after next dose
1.26 - 1.5	no	decrease dose by 30%	4h after next dose
1.6 - 2.0	yes for 3h	decrease dose by 40%	Before next dose then 4h after next dose

Anti Factor Xa level U/ml?	Hold Next Dose?	Dose Change?	Repeat Anti Factor Xa?
> 2.0	yes, until anti factor Xa level is <0.5 U/ml	decrease dose by 50%	Before next dose is administered, if >0.5U/ml (therapeutic level), do not give next enoxaparin dose & repeat anti Xa level in 12h. When level <0.5 U/ml, administer 50% original dose.

Tab. 9-2 Enoxaparin, low molecular weight heparin dosing (⁶Monagle, P, et al.)

9.9.2 Oral Anticoagulation Therapy

Stage	INR	Action
Day 1	1.0 - 1.8	0.2 mg/kg orally
Day 2-4	1.1 - 1.3	repeat day 1 loading dose
	1.4 - 1.9	50% of day 1 loading dose
	2.0 - 3.0	50% of day 1 loading dose
	3.1 - 3.5	25 % of day 1 loading dose
	> 3.5	hold dosing until INR is < 3.5

Tab. 9-3 Warfarin loading dose to maintain an INR of 2.7 - 3.5 (⁶Monagle, P, et al.)

Stage	INR	Action
Maintenance : = Day 5 and long term	1.1 - 1.9	increase dose by 40 -50%
	2.0 - 2.4	increase dose by 10%
	2.7 - 3.5	no change
	3.6 - 4.0	administer next dose at 50% then restart at 20% less maintenance dose
	4.1- 5.0	hold one dose then 20% less maintenance dose

Tab. 9-4 Warfarin Maintenance Dosing for Day 5 and longer to maintain INR 2.7-3.5

INR 2.7 to 3.5	use only warfarin p.o.
INR < 2.7	use warfarin plus enoxaparin as outlined in section 5 until INR ≥ 2.7

Tab. 9-5 Drugs and Dose for specific INR range

10 Troubleshooting and correcting faults

NOTE: This chapter omits safety instructions, information and procedures that refer to the *Ikus* exclusively. Please refer also to the IFU.

For information on *Ikus* error messages and corrective measures please refer to the IFU.

 **HOTLINE**

Notify Berlin Heart! 866.249.0128

Problem	Cause of problem / action to be taken
Visible blood pump faults	Replace the pump, see section 10.1: Replacing the blood pump(s), page 89.
Deposits in the pump	Initial deposits: check anticoagulation status and adjust therapy if necessary. If floating deposits are detected (may cause thromboembolic complication): replace the pump, see section 10.1: Replacing the blood pump(s), page 89.
Pump is filling or ejecting blood incorrectly	Assess the condition of the patient and the hemodynamic status. If necessary, adjust the system parameters.

Tab. 10-1 Possible problems

10.1 Replacing the blood pump(s)

 **WARNING**

When replacing a blood pump, follow the instruction given here. Otherwise the duration of the pump stop will be prolonged and the patient might suffer from inadequate support.

The blood pump may only be replaced under sterile conditions!

When connecting the blood pump(s), pay attention to the direction of the arrows on the inflow and outflow stubs! These show the direction of the blood flow.

The cable tie covering the *EXCOR* cannula on the stub of the blood pump should be removed carefully. Important: never use a sharp instrument, for example, a scalpel or scissors, to remove the cable tie. This may cause damage to the cannula.

BVAD: If the left pump is being replaced, the right pump must also be stopped while the pump is being replaced. Otherwise there is the risk of pulmonary edema.

NOTICE

If the replacement pump has a larger volume than the one being replaced,

- the use of a connector set must be considered
- the corresponding parameter in the view *Pump size and single-step mode* must be updated

IMPORTANT: When 2 blood pumps need to be replaced, replace the right blood pump in the first place, subsequently replace the left blood pump.

IMPORTANT: Sedate the patient if necessary and administer a bolus of Heparin according to the the anticoagulation protocol.

10.1.1 Preparing a replacement blood pump

Material

- 1 replacement blood pump of appropriate type and size
- 1 driving tube, red or blue
- 1 accessory set (for blood pumps with PU valves) with tube connecting set;
IMPORTANT: Only the cable ties and cable tie guns provided should be used.

► INSTRUCTION

1. Bring membrane to the end-of-diastole position, position de-airing needle, rinse and fill pump with sterile injectable saline (see section 10.2: Driving blood pump(s) with the manual pump, page 92).
2. Connect the driving tube to the respective driving tube connector of the pump.
3. Place the pump, ready for connection, with the titanium connectors pointing upwards.

10.1.2 Replacing the right blood pump (RVAD/ BVAD)

Material

- 1 prepared replacement blood pump (see section 10.1.1: Preparing a replacement blood pump, page 90)
- 1 tube connecting set (cable tie, cable-tie gun), included in the accessory set. Only the cable ties and cable tie guns provided should be used.

Stopping the right blood pump and detaching the blood pump from Ikus

► INSTRUCTION

1. Bring the patient into the Trendelenburg position.
2. The cable tie covering the EXCOR cannula on the stub of the blood pump should be removed carefully. Important: never use a sharp instrument, for example, a scalpel or scissors, to remove the cable tie. This may cause damage to the cannula. Check cannulae immediately to make sure they are not damaged.
3. If necessary log into the monitor program by entering user ID and password, confirming the password with <Enter>.
4. BVAD: Reduce rate of left blood pump to 30 bpm. Use <<->/ <->> to navigate cursor to the respective field of the parameter table, then use <↓> to adapt value. Confirm with <Enter>.
5. In the monitor program, select the option **Pause left** respectively **Pause right** and press <Enter> to confirm. Respond to the prompt in the dialog window by pressing the <X> key or the <1> key. The right blood pump will stop.
RVAD: **Pause left**
BVAD: **Pause right**
The view *Pump size and single-step* mode is displayed.
6. As soon as the right pump has stopped, clamp off the cannulae beneath the right pump to be replaced and slide the cannulae off the pump. If it is necessary to clamp any other part of the cannula that is not covered with velour, cover the part of the cannula that will be clamped with a gauze sponge.
7. Check cannulae for visible deposits. If necessary, remove these deposits carefully.
8. Remove the driving tube of the pump to be replaced from the connector. To do so, take hold of the release sleeve and pull this out of the connector.

Connect new right blood pump to the IkuS**INSTRUCTION**

1. Fill the free ends of the cannulae with sterile saline solution. Make sure that all air has been removed. Connect the prepared replacement pump to the cannulae.
2. Plug the new driving tube into the freed connector. The plug snaps into place clearly audible.
3. Check that the plug is securely connected. To do so, grip the plug body above the release sleeve and pull on it. Do not pull from the release sleeve, and never from the tube!
4. Release the tube clamps from the cannulae.

Starting the IkuS**INSTRUCTION**

1. Move the cursor to the field **step left** (RVAD) respectively **step right** (BVAD).
2. RVAD: Confirm **Step left** with <Enter> to trigger a single step.
BVAD: Confirm **Step right** with <Enter> to trigger a single step.
3. If any air bubbles are visible remove them via the de-airing needle. When all air has been completely removed from the left pump: remove the de-airing needle.
4. Move cursor to the **OK** field and press <Enter> to confirm. The driving unit starts up again using the defined parameters.
5. Check whether the pump is filling correctly and, if necessary, adjust the parameters.
6. Secure all connections with cable ties. See 6.12: Securing the connections, page 70.

10.1.3 Replacing the left blood pump (LVAD/ BVAD)**WARNING**

BVAD: If the left pump is being replaced, the right pump must also be stopped while the pump is being replaced. Otherwise there is the risk of pulmonary edema.

Material

- 1 prepared replacement blood pump (see section 10.1.1: Preparing a replacement blood pump, page 90)
- 1 tube connecting set (cable tie, cable-tie gun), included in the accessory set. Only the cable ties and cable tie guns provided should be used.

Stopping the left blood pump and detaching the blood pump from IkuS**INSTRUCTION**

1. Bring the patient into the Trendelenburg position.
2. The cable tie covering the *EXCOR* cannula on the stub of the blood pump should be removed carefully. Important: never use a sharp instrument, for example, a scalpel or scissors, to remove the cable tie. This may cause damage to the cannula. Check cannulae immediately to make sure they are not damaged.
3. If necessary log into the monitor program by entering user ID and password, confirming the password with <Enter>.
4. In the monitor program, select the option **Pause left** respectively **Drive pause** and press <Enter> to confirm. Respond to the prompt in the dialog window by pressing the <X> key or the <1> key. The left blood pump respectively both blood pumps will stop.
LVAD: **Pause left**. The view *Pump size and single-step* mode is displayed.
BVAD: **Drive pause**. The view *Select operating mode* is displayed.
5. As soon as the right pump(s) has/ have stopped, clamp off the cannulae beneath the left pump to be replaced and slide the cannulae off the pump. If it is necessary to clamp any other part of the cannula that is not covered with velour, cover the part of the cannula that will be clamped with a gauze sponge.
6. Check cannulae for visible deposits. If necessary, remove these deposits carefully.
7. Remove the driving tube of the left pump to be replaced from the connector. To do so,

Driving blood pump(s) with the manual pump

take hold of the release sleeve and pull this out of the connector.

Connect new left blood pump to the Ikus

► INSTRUCTION

1. Fill the free ends of the cannulae with sterile saline solution. Make sure that all air has been removed. Connect the prepared replacement pump to the cannulae.
2. Plug the new driving tube into the freed connector. The plug snaps into place clearly audible.
3. Check that the plug is securely connected. To do so, grip the plug body above the release sleeve and pull on it. Do not pull from the release sleeve, and never from the tube!
4. Release the tube clamps from the cannulae.

Starting the Ikus

► INSTRUCTION

1. Move the cursor to the field **step left**.
2. Confirm **Step left** with <Enter> to trigger a single step.
3. If any air bubbles are visible remove them via the de-airing needle. When all air has been completely removed from the left pump: remove the de-airing needle.
4. Move cursor to the **OK** field and press <Enter> to confirm. The driving unit starts up again using the defined parameters.
5. Check whether the pump is filling correctly and, if necessary, adjust the parameters.
6. Secure all connections with cable ties.
7. .

10.2 Driving blood pump(s) with the manual pump




Fig. 10-1 Patient on *Ikus*



Fig. 10-2 Patient on manual pump

This is necessary if ...

- the power supply to the *Ikus* cannot be ensured
- the *Ikus* has to be restarted (e.g. emergency operating mode) and there is no replacement *Ikus* available

WARNING

The use of the manual pump is only permitted for medical personnel trained in the use of it.

Pay attention to the colored markings on the driving tubes and on the connectors of the manual pump. Otherwise, there is a risk of lung edema.

Always keep manual pump attached to the *Ikus*. Otherwise in an emergency situation the adequate support of the patient is not guaranteed.

Call one or more persons to assist. Otherwise in an emergency situation the adequate support of the patient is not guaranteed.

The driving tubes and cannulae should be arranged in a bend-free position. Otherwise in an emergency situation the adequate support of the patient is not guaranteed.

When operating the manual pump with 1 hand, do not block the valves with your feet (see valve "c" in Fig. 10-3, page 94 and in Fig. 10-4, page 94).

NOTICE


Seal the connector(s) on the *Ikus* immediately after removing the driving tube(s) in order to avoid contaminates from entering the system.

IMPORTANT: In biventricular mode: the blood pumps are driven asynchronously by the manual pump.

INSTRUCTION


1. The patient is lying down.
2. Disconnect the driving tube(s) from the *Ikus*. To do so, take hold of the release sleeve and pull this out of the connector.
3. Connect the driving tube(s) to the manual pump.
Important: Observe the colored markings.
4. Check that the plug is securely connected. To do so, grip the plug body above the release sleeve and pull on it. (see f in figure Fig. 10-4, page 94) Do not pull from the release sleeve, and never from the tube!
5. Pump steadily and rhythmically at roughly 60 to 80 strokes per minute. *Important:* Move the piston so far that the membrane reaches its final position. The piston need not necessarily be moved to its end position.
6. Perform a visual check of the blood pump to verify that the membrane is moving and that blood is being pumped.

Driving blood pump(s) with the manual pump



- a Connector for driving tube with blue marking
- b Connector for driving tube with blue marking
- c Valve

Fig. 10-3 Manual pump



- d Base plate
- e Release sleeve
- f Area above release sleeve

Fig. 10-4 Plug on the driving tube




Fig. 10-5 Examples to operate the manual pump

The manual pump can be operated with both hands or with one hand (placing the pump between the feet). Alternating between two-handed or one-handed pumping, as well as using the left or right hand, is allowed. When doing so, care of the patient must remain ensured.

11 Weaning and Explantation for BTR and BTT

11.1 Weaning Procedure

11.1.1 Introduction

This document summarizes the clinical guideline for weaning and explantation of the EXCOR. The decision to wean the EXCOR should be made cautiously after careful review of all available clinical and laboratory data. This document should be considered a guideline only. As always treatment must be individualized to each patient based on his/her unique clinical circumstances.

It is important to recognize that prolonged pump stoppage and operation of the device at lower beat rates is not recommended due to the risks of blood stagnation and thrombus formation. This risk increases with the smaller blood pumps (e.g. 10, 25 and 30 ml devices) where the luminal sizes and flow rates are the lowest. Therefore, a size-based guideline has been developed to test the adequacy of the native circulation without a prolonged pump stoppage using a combination of gradual weaning, brief pump stoppages, careful anticoagulation monitoring, invasive hemodynamic testing, and a brief afterload challenge. It is not recommended that weaning proceed unless all parameters especially those pertaining to anticoagulation have been fully optimized. This protocol reflects the most recent understanding of the safest possible weaning strategy based on the collective US and European experience to date. Consultation with *Berlin Heart, Inc.* prior to weaning and explantation is strongly recommended.

11.1.2 Indication

Weaning may be considered in children supported with the EXCOR judged to have sufficient evidence of myocardial recovery to provide adequate systemic perfusion independent of VAD support.

11.1.3 Eligibility Criteria



WARNING

Continuous reassessment of eligibility criteria is critical to reducing the risks associated with weaning of VAD support. At all times each of the weaning criteria must be satisfied in order to proceed with the weaning protocol.

Special attention must be taken to ensure the patient's anticoagulation status remains within the targeted range.

Weaning of the EXCOR may be considered in subjects who meet the following eligibility criteria:

- LVEDD within normal limits (<98th percentile, or Z-score of +2)
- EF = 45% (i.e. no less than mild dysfunction)
- Lactate <3 mmol/L
- No clinical evidence of thromboembolism or bleeding
- Anticoagulation markers within target parameters

11.1.4 Weaning Protocol



WARNING

Rates < 60 bpm are intended to be used only for implantation and explantation. Never use the *Ikus* with a rate < 60 bpm without constant supervision.

If the patient does not meet the eligibility criteria at any time during the weaning process: Resume pumping at rate prior to any weaning (initial rate, IR).

The weaning protocol can be divided into 5 steps and generally takes one week to complete.

- Day 0 (and throughout the weaning process). Confirmation of eligibility criteria for weaning.
- Day 0. Acute weaning challenge
- Day 1-4. Graduated weaning challenge with non-invasive assessment (echo).
- Day 5. Pump stoppage with invasive hemodynamic assessment with afterload challenge.
- Day 6. Pump stoppage with invasive hemodynamic assessment in OR (full anticoagulation).

This size-based weaning protocol accounts for physiologic differences in heart rate and stroke volume observed in children of varying ages.

11.1.5 10 ml pump

The individual weaning progress is based upon the following parameters:

Parameter	Explanation	Abbr.	Value
initial rate	rate prior to any weaning	IR	Please enter: IR = _____ bpm
weaning rate	lowest rate achieved during weaning process, depends on pump size	WR	50 bpm
total weaning interval	Difference between initial rate and explanation rate: TWI = IR – WR	TWI	Please enter: IR ____ bpm – WR 50 bpm = TWI ____ bpm
reduced rate	rate resumed at the end of day 1 to 3	RR ₁ to RR ₃	Please refer to table Tab. 11-2.

Tab. 11-1 Important parameters for weaning progress

Reduced rate (RR _x)	Calculation
RR ₁	Please enter: $RR_1 = WR\ 50\ bpm + 0.75 \times TWI\ (\ ___ \text{ bpm}) = ____ \text{ bpm}$
RR ₂	Please enter: $RR_2 = WR\ 50\ bpm + 0.50 \times TWI\ (\ ___ \text{ bpm}) = ____ \text{ bpm}$
RR ₃	Please enter: $RR_3 = WR\ 50\ bpm + 0.25 \times TWI\ (\ ___ \text{ bpm}) = ____ \text{ bpm}$

Tab. 11-2 Reduced rate day 1 to day 3

10 ml pump Weaning Sequence

10 ml pump Weaning Sequence		
Day 0	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> Administer unfractionated heparin (UFH) 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. After 5 minutes, reduce the pump rate step-wise from IR (_____ bpm) to 30 bpm in increments of 5 bpm q5 min. After 5 minutes at 30 bpm, reassess LV size and function. After an additional 5 minutes (i.e. total time = 10 min at 30 bpm), stop the pump for 3 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. After 3-minute pump stop, reconnect pump to Ikus and resume pump speed at IR(____ bpm). 	<p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p>	<p><input type="checkbox"/> NO - STOP</p> <p><input type="checkbox"/> YES - Proceed</p> <p>MD: _____</p>
Day 1	<p>After confirmation of eligibility criteria, the following steps should be performed sequentially under echo guidance ¹:</p> <ol style="list-style-type: none"> Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. After 5 minutes, reduce the pump rate step-wise by from the IR (_____ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. After a total time of 10 min at 30 bpm, stop the pump for 3 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. After 3-minute pump stop, reconnect pump to Ikus and resume pumping at rate RR₁ (____ bpm). 	<p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p>	<p><input type="checkbox"/> NO - STOP</p> <p><input type="checkbox"/> YES - Proceed</p> <p>MD: _____</p>
Day 2	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. After 5 minutes, reduce the pump rate step-wise from RR₁ (_____ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. After a total time of 20 min at 30 bpm, stop the pump for 3 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. After 3-minute pump stop, reconnect pump to Ikus and resume pumping at RR₂ (____ bpm). 	<p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p>	<p><input type="checkbox"/> NO - STOP</p> <p><input type="checkbox"/> YES - Proceed</p> <p>MD: _____</p>

10 ml pump Weaning Sequence		
Day 3	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 2. After 5 minutes, reduce the pump rate step-wise from RR₂ (____ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. 3. Initiate exercise with gentle age-appropriate play tasks (e.g. rattle, clapping) as clinically appropriate, where possible 4. After a total time of 30 min at 30 bpm, stop the pump for 3 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. After 3-minute pump stop, reconnect pump to Ikus and resume pumping at RR₃ (____ bpm). 	<p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p>	<p><input type="checkbox"/> NO - STOP</p> <p><input type="checkbox"/> YES - Proceed</p> <p>MD: _____</p>
Day 4	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 2. After 5 minutes, reduce pump rate step-wise from RR₃ (____ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. 3. Initiate exercise with gentle age-appropriate play tasks (e.g. rattle, clapping) as clinically appropriate, where possible. 4. After a total time of 30 min at 30 bpm, stop the pump for 3 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. 5. After a 3-minute pump stop: If the patient meets all eligibility criteria, reconnect pump to Ikus and resume pumping at WR (50 bpm). If the patient does not meet all criteria, reconnect Ikus and resume pumping at IR. 	<p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p>	<p><input type="checkbox"/> NO - STOP</p> <p><input type="checkbox"/> YES - Proceed</p> <p>MD: _____</p>

10 ml pump Weaning Sequence	
Day 5	<p>After confirmation of eligibility criteria, the following steps should be performed in the cath lab under echo guidance¹:</p> <ol style="list-style-type: none"> 1. Obtain standard access for RHC (if possible with out sedation). 2. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 3. After 5 minutes, reduce the pump rate step-wise from WR (50 bpm) to 30 bpm in increments of 5 bpm q5 min. After 5 minutes at 30 bpm, reassess LV size and function. 4. After a total time of 10 min at 30 bpm, stop the pump for 3 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. 5. After 3 minutes, initiate norepinephrine infusion at 0.01 mcg/kg/min IV gtt titrated to MAP 20% above baseline x 5 min. While doing so, proceed pumping manually twice q30 seconds. 6. If LV size and function acceptable, proceed pumping manually twice q30 seconds for 3 min. While doing so, reassess LV size & function, and record RAP, PAP, PCWP and MVS. 7. After 6-minute pump stop: If the patient meets all eligibility criteria, reconnect pump to Ikus and resume pumping at 50 bpm until the actual surgical procedure of explantation takes place. If the patient does not meet all criteria, reconnect Ikus and resume pumping at IR.
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p> <p><input type="checkbox"/> NO - STOP <input type="checkbox"/> YES - Proceed MD:_____</p>

¹ TTE unless echo windows insufficient. The last weaning increment may be less than 5 bpm if the wean interval is not a multiple of 5.

11.1.6 25/ 30 ml pump

The individual weaning progress is based upon the following parameters:

Parameter	Explanation	Abbr.	Value
initial rate	rate prior to any weaning	IR	Please enter: IR = _____ bpm
weaning rate	lowest rate achieved during weaning process, depends on pump size	WR	40 bpm
total weaning interval	Difference between initial rate and explanation rate: TWI = IR – WR	TWI	Please enter: IR ____ bpm – WR 40 bpm = TWI ____ bpm
reduced rate	rate resumed at the end of day 1 to 3	RR ₁ to RR ₃	Please refer to table Tab. 11-4.

Tab. 11-3 Important parameters for weaning progress

Reduced rate (RR _x)	Calculation
RR ₁	Please enter: $RR_1 = WR\ 40\ bpm + 0.75 \times TWI\ (\ ___ \text{ bpm}) = ____ \text{ bpm}$
RR ₂	Please enter: $RR_2 = WR\ 40\ bpm + 0.50 \times TWI\ (\ ___ \text{ bpm}) = ____ \text{ bpm}$
RR ₃	Please enter: $RR_3 = WR\ 40\ bpm + 0.25 \times TWI\ (\ ___ \text{ bpm}) = ____ \text{ bpm}$

Tab. 11-4 Reduced rate day 1 to day 3

25/ 30 ml pump Weaning Sequence

25/ 30 ml pump Weaning Sequence		
Day 0	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> Administer unfractionated heparin (UFH) 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. After 5 minutes, reduce the pump rate step-wise from IR (_____ bpm) to 30 bpm in increments of 5 bpm q5 min. After 5 minutes at 30 bpm, reassess LV size and function. After an additional 5 minutes (i.e. total time = 10 min at 30 bpm), stop the pump for 5 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. After 5-minute pump stop, reconnect pump to Ikus and resume pump speed at IR(____ bpm). 	<p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p>	<p><input type="checkbox"/> NO - STOP</p> <p><input type="checkbox"/> YES - Proceed</p> <p>MD: _____</p>
Day 1	<p>After confirmation of eligibility criteria, the following steps should be performed sequentially under echo guidance ¹:</p> <ol style="list-style-type: none"> Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. After 5 minutes, reduce the pump rate step-wise by from the IR (_____ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. After a total time of 10 min at 30 bpm, stop the pump for 5 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. After 5-minute pump stop, reconnect pump to Ikus and resume pumping at rate RR₁ (____ bpm). 	<p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p>	<p><input type="checkbox"/> NO - STOP</p> <p><input type="checkbox"/> YES - Proceed</p> <p>MD: _____</p>
Day 2	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. After 5 minutes, reduce the pump rate step-wise from RR₁ (_____ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. After a total time of 20 min at 30 bpm, stop the pump for 10 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. After 10-minute pump stop, reconnect pump to Ikus and resume pumping at RR₂ (____ bpm). 	<p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p>	<p><input type="checkbox"/> NO - STOP</p> <p><input type="checkbox"/> YES - Proceed</p> <p>MD: _____</p>

25/ 30 ml pump Weaning Sequence		
Day 3	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 2. After 5 minutes, reduce the pump rate step-wise from RR₂ (____ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. 3. Initiate exercise with gentle age-appropriate play tasks (e.g. patty cake) as clinically appropriate, where possible 4. After a total time of 30 min at 30 bpm, stop the pump for 10 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. 5. After 10-minute pump stop, reconnect pump to Ikus and resume pumping at RR₃ (____ bpm). 	<p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p>	<p><input type="checkbox"/> NO - STOP</p> <p><input type="checkbox"/> YES - Proceed</p> <p>MD: _____</p>
Day 4	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 2. After 5 minutes, reduce pump rate step-wise from RR₃ (____ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. 3. Initiate exercise with gentle age-appropriate play tasks (e.g. patty cake) as clinically appropriate, where possible. 4. After a total time of 30 min at 30 bpm, stop the pump for 15 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. 5. After a 15-minute pump stop: If the patient meets all eligibility criteria, reconnect pump to Ikus and resume pumping at WR (40 bpm). If the patient does not meet all criteria, reconnect Ikus and resume pumping at IR. 	<p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p>	<p><input type="checkbox"/> NO - STOP</p> <p><input type="checkbox"/> YES - Proceed</p> <p>MD: _____</p>

25/ 30 ml pump Weaning Sequence	
Day 5	<p>After confirmation of eligibility criteria, the following steps should be performed in the cath lab under echo guidance¹:</p> <ol style="list-style-type: none"> 1. Obtain standard access for RHC (if possible with out sedation). 2. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 3. After 5 minutes, reduce the pump rate step-wise from WR (50 bpm) to 30 bpm in increments of 5 bpm q5 min. After 5 minutes at 30 bpm, reassess LV size and function. 4. After a total time of 30 min at 30 bpm, stop the pump for 15 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. 5. After 15 minutes, initiate norepinephrine infusion at 0.01 mcg/kg/min IV gtt titrated to MAP 20% above baseline x 5 min. While doing so, proceed pumping manually twice q30 seconds. 6. If LV size and function acceptable, proceed pumping manually twice q30 seconds for 5 min. While doing so, reassess LV size & function, and record RAP, PAP, PCWP and MVS. 7. After 20-minute pump stop: If the patient meets all eligibility criteria, reconnect pump to Ikus and resume pumping at 50 bpm until the actual surgical procedure of explantation takes place. If the patient does not meet all criteria, reconnect Ikus and resume pumping at IR.
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p> <p><input type="checkbox"/> NO - STOP <input type="checkbox"/> YES - Proceed MD: _____</p>

¹ TTE unless echo windows insufficient. The last weaning increment may be less than 5 bpm if the wean interval is not a multiple of 5.

11.1.7 50/ 60 ml pump

The individual weaning progress is based upon the following parameters:

Parameter	Explanation	Abbr.	Value
initial rate	rate prior to any weaning	IR	Please enter: IR = _____ bpm
weaning rate	lowest rate achieved during weaning process, depends on pump size	WR	30 bpm
total weaning interval	Difference between initial rate and explanation rate: TWI = IR – WR	TWI	Please enter: IR ____ bpm – WR 30 bpm = TWI ____ bpm
reduced rate	rate resumed at the end of day 1 to 3	RR ₁ to RR ₃	Please refer to table Tab. 11-6.

Tab. 11-5 Important parameters for weaning progress

Reduced rate (RR _x)	Calculation
RR ₁	Please enter: $RR_1 = WR\ 30\ bpm + 0.75 \times TWI\ (\ ___ \text{ bpm}) = ____ \text{ bpm}$
RR ₂	Please enter: $RR_2 = WR\ 30\ bpm + 0.50 \times TWI\ (\ ___ \text{ bpm}) = ____ \text{ bpm}$
RR ₃	Please enter: $RR_3 = WR\ 30\ bpm + 0.25 \times TWI\ (\ ___ \text{ bpm}) = ____ \text{ bpm}$

Tab. 11-6 Reduced rate day 1 to day 3

50/ 60 ml pump Weaning Sequence

50/ 60 ml pump Weaning Sequence		
Day 0	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> Administer unfractionated heparin (UFH) 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. After 5 minutes, reduce the pump rate step-wise from IR (_____ bpm) to 30 bpm in increments of 5 bpm q5 min. After 5 minutes at 30 bpm, reassess LV size and function. After an additional 5 minutes (i.e. total time = 10 min at 30 bpm), stop the pump for 5 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. After 5-minute pump stop, reconnect pump to Ikus and resume pump speed at IR(____ bpm). 	<p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p>	<p><input type="checkbox"/> NO - STOP</p> <p><input type="checkbox"/> YES - Proceed</p> <p>MD: _____</p>
Day 1	<p>After confirmation of eligibility criteria, the following steps should be performed sequentially under echo guidance ¹:</p> <ol style="list-style-type: none"> Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. After 5 minutes, reduce the pump rate step-wise by from the IR (_____ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. After a total time of 10 min at 30 bpm, stop the pump for 10 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. After 10-minute pump stop, reconnect pump to Ikus and resume pumping at rate RR₁ (____ bpm). 	<p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p>	<p><input type="checkbox"/> NO - STOP</p> <p><input type="checkbox"/> YES - Proceed</p> <p>MD: _____</p>
Day 2	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. After 5 minutes, reduce the pump rate step-wise from RR₁ (_____ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. After a total time of 15 min at 30 bpm, stop the pump for 15 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. After 15-minute pump stop, reconnect pump to Ikus and resume pumping at RR₂ (____ bpm). 	<p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p>	<p><input type="checkbox"/> NO - STOP</p> <p><input type="checkbox"/> YES - Proceed</p> <p>MD: _____</p>

50/ 60 ml pump Weaning Sequence		
Day 3	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 2. After 5 minutes, reduce the pump rate step-wise from RR₂ (____ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. 3. Initiate exercise with gentle age-appropriate play tasks (e.g. ambulate) as clinically appropriate, where possible 4. After a total time of 30 min at 30 bpm, stop the pump for 20 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. 5. After 20-minute pump stop, reconnect pump to Ikus and resume pumping at RR₃ (____ bpm). 	<p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p>	<p><input type="checkbox"/> NO - STOP</p> <p><input type="checkbox"/> YES - Proceed</p> <p>MD: _____</p>
Day 4	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 2. After 5 minutes, reduce pump rate step-wise from RR₃ (____ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. 3. Initiate exercise with gentle age-appropriate play tasks (e.g. ambulate) as clinically appropriate, where possible. 4. After a total time of 30 min at 30 bpm, stop the pump for 30 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. 5. After a 30-minute pump stop: If the patient meets all eligibility criteria, reconnect pump to Ikus and resume pumping at WR (30 bpm). If the patient does not meet all criteria, reconnect Ikus and resume pumping at IR. 	<p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p>	<p><input type="checkbox"/> NO - STOP</p> <p><input type="checkbox"/> YES - Proceed</p> <p>MD: _____</p>

50/ 60 ml pump Weaning Sequence		
Day 5	<p>After confirmation of eligibility criteria, the following steps should be performed in the cath lab under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Obtain standard access for RHC (if possible with out sedation). 2. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 3. Assess LV size and function to obtain data for comparison. 4. Stop the pump for 15 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. 5. After 15 minutes, initiate norepinephrine infusion at 0.01 mcg/kg/min IV gtt titrated to MAP 20% above baseline x 5 min. While doing so, proceed pumping manually twice q30 seconds. 6. If LV size and function acceptable, proceed pumping manually twice q30 seconds for 15 min. While doing so, reassess LV size & function, and record RAP, PAP, PCWP and MVS. 7. After 30-minute pump stop: If the patient meets all eligibility criteria, reconnect pump to Ikus and resume pumping at 50 bpm until the actual surgical procedure of explantation takes place. If the patient does not meet all criteria, reconnect Ikus and resume pumping at IR. 	<hr/> <hr/> <hr/> <hr/> <hr/> <hr/> <hr/> <hr/>
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p>	<input type="checkbox"/> NO - STOP <input type="checkbox"/> YES - Proceed MD: _____

¹ TTE unless echo windows insufficient. The last weaning increment may be less than 5 bpm if the wean interval is not a multiple of 5.

11.1.8 Explantation Criteria

NOTICE	<p>ASA and dipyridamole should be discontinued 24-hours prior to device explantation; coumadin/Enoxaparin should be transitioned back to unfractionated heparin (titrated to therapeutic levels).</p> <hr/> <p>Milrinone 0.75 µg/kg/min should be started 12 hours prior explantation. ACE inhibitor, β-Blocker and Spirinolactone should be not stopped.</p>
---------------	---

In the operating room, explantation should be considered if the following criteria are met with the pump stopped for 20 minutes (after anticoagulation has been established in the target range for cardiopulmonary bypass):

- LVEDD less than 98th percentile (Z-score less than +2)
- EF ≥45% (i.e. no more than mild ventricular dysfunction)
- Normotensive on only Milrinone (no other inotropes)
- Lactate <3 mmol/L
- LVEDP < 12 mm Hg
- Resting CI of > 2.8 L/min/m²

Surgery should be performed without Cardiopulmonary Bypass. Control all bleeding immediately during and post implantation.

11.2 Explantation for BTR

11.2.1 Explantation with univentricular support

The procedure is analogous to that used after BTT (see 11.3: Explantation for BTT, page 109). Sew over all anastomosis areas where cannulae were placed.

11.2.2 Explantation after biventricular support

Stopping the right pump

INSTRUCTION

1. Select **Pause right** (see Fig. 11-1), then press <Enter> to confirm. Respond to the prompt in the dialog window by pressing the <X> key or the <1> key. The right pump will stop. The view *Pump size and single-step mode* is shown (see Instructions for Use Rev. 6, Fig. 6-7 on page 71). The cursor is located on the **OK** field.




Fig. 11-1 Pause right

2. Unplug the driving tube of the right pump from the connector on the *Ikus*. Use the seal plug to seal the connector.
3. To confirm the **OK** selection, press <Enter>. The *Ikus* continues running. The screen shows the standard view.

Switching the *Ikus* off

WARNING

The *Ikus* power switch (toggle switch) should always be in the [I] position, even if the main switch (key switch) is in the [O] position!. Otherwise there is a risk that the drive may fail in future due to the *Ikus* batteries being totally discharged.

Always follow the above sequence of operations. Always use the key switch to switch off the *Ikus*.

Do not switch the *Ikus* off unless the batteries are fully charged. Leave the *Ikus* switched on until all yellow LEDs light up, then switch off the *Ikus* with main switch (key switch).

Keep all driving tube connectors covered at all times when not in use.

INSTRUCTION

1. Put the patient on cardiopulmonary bypass (CPB).
2. Disconnect the driving tubes and connect both tank units to the *Ikus*.
3. Leave the *Ikus* running with the tank units until the patient is stable on CPB and the blood pumps have been explanted.
4. Next in the monitor program, select the option **Drive OFF** (see figure Fig. 11-2, page 109) and press <Enter> to confirm.
5. Respond to the prompt in the dialog window by pressing the <X> key or the <1> key. The system stops operation immediately and writes an operating log.
6. Disconnect the driving tube(s) from the connector(s). To do so, take hold of the plug's release sleeve and pull the plug out of the connector.

7. Use the seal plugs to seal the driving tube connector sockets.
8. Wait until the log has been completed. When the message **Switch drive off with main switch!** appears, press <F10> to shut down the monitor program. Confirm by pressing the <X> key or the <1> key.
9. Select **3. End** (<3>, see figure Fig. 11-3, page 110) in the start menu and switch off the laptop.
10. Switch the *Ikus* off, provided that the batteries are fully charged. To do so, turn the key switch to [0] position.

11.3 Explanation for BTT

NOTICE

When planning and timing the transplantation, be aware that massive adhesions may exist in the transplant recipient.

Preparing the donor organ

ADVICE

Leave adequate lengths of the aorta and the pulmonary artery attached to the donor organ in order to be able to continue using those parts of the original vessels used for anastomosis of the VAD cannulae.

Leave the *Ikus* running with the tank units until the patient is stable on CPB and the blood pumps have been explanted.

Switching the *Ikus* off

WARNING

The *Ikus* power switch (toggle switch) should always be in the [I] position, even if the main switch (key switch) is in the [0] position! Otherwise there is a risk that the drive may fail due to the *Ikus* batteries being totally discharged.

CAUTION

Always follow the above sequence of operations. Always use the key switch to switch off the *Ikus*.

Do not switch the *Ikus* off unless the batteries are fully charged. To do this leave the *Ikus* switched on until all yellow LEDs light up, then switch the *Ikus* off using the key switch.

INSTRUCTION

1. Put the patient on cardiopulmonary bypass.
2. Disconnect the driving tubes and connect both tank units to the *Ikus*.
3. Leave the *Ikus* running with the tank units until the patient is stable on CPB and the blood pumps have been explanted.
4. Next in the monitor program, select the **Drive OFF** option and press <Enter> to confirm (see figure Fig. 11-2).




Fig. 11-2 Drive OFF

5. Respond to the prompt in the dialog window by pressing the <X> key or the <1> key. The system stops operation immediately and writes an operating log.

6. Disconnect the driving tube(s) from the connector(s). To do so, take hold of the release sleeve and pull this out of the connector.
7. Use the seal plugs to seal the driving tube connectors.
8. Wait until the log has been completed. When the message **Switch drive off with main switch!** appears, press <F10> to shut down the monitor program. Confirm by pressing the <X> key or the <1> key.
9. Select **3. End** (<3>, see figure Fig. 11-3) in the start menu and switch off the laptop.
10. Switch the *Ikus* off, provided that the batteries are fully charged. To do so, turn the key switch to [0] position.

```
          Berlin Heart <R>
          Ikus2000 <R> Rev. 2.1
          Build 2009.06
          Copyright <C> 1997-2009

1. Start Program
2. Entry codes
3. End
4. Save data
5. Change date or time
6. Change language

Input:
```

Fig. 11-3 Start menu

Removing the VAD cannulae

➤ INSTRUCTION

1. Clamp off the cannulae.
2. Disconnect the pump from the cannulae.
3. Remove the cannulae. Sew over the anastomosis areas of the atrium.

The remaining procedure is the same as for any primary orthotopic heart transplantation.

12 Appendix

12.1 Overview: Product range and possible combinations

12.1.1 Blood pumps

Designation	Article number	Inflow / outflow [mm]
Blood pump PU valves 10 ml in/out ø 6 mm	P10P-001	6/6
Blood pump PU valves 25 ml in/out ø 9 mm	P25P-001x01	9/9
Blood pump PU valves 30 ml in/out ø 9 mm	P30P-001x01	9/9
Blood pump PU valves 50 ml in/out ø 12 mm	P50P-001	12/12
Blood pump PU valves 60 ml in/out ø 12 mm	P60P-001	12/12

Tab. 12-1 Blood Pump

12.1.2 Overview: Relationship: body weight – pump size




Fig. 12-1 Relationship: body weight - pump size

12.1.3 LV apex cannulae

Designation	Article number	Length of head [mm]	Overall length [mm]	Lumen diameter [mm]
Apex cannula for infants	C14A-040	14	220	5
Apex cannula for small children	C18A-020	18	250	6
Apex cannula for children, staged	C22A-004	28	270	12, 9; head 9
Apex cannula	C27A-001	38	265	12

Tab. 12-2 LV apex cannulae

12.1.4 Atrial cannulae

Designation	Article number	Length of head [mm]	Length of corpus [mm]	Lumen diameter [mm]
Atrial cannula for infants	C15V-040	15	200	5
Atrial cannula for small children	C19V-020	19	250	6
Atrial cannula for children, staged (with mandrin)	C22V-004	22	280	9, 12; head 9
Atrial cannula for children, staged (with mandrin)	C25V-004	25	280	9, 12; head 9
Atrial cannula (with mandrin)	C22V-002	22	330	12
Atrial cannula (with mandrin)	C26V-002	26	330	12

Tab. 12-3 Atrial cannulae

12.1.5 Arterial cannulae

Designation	Article number	Head angle [°]	Length of corpus [mm]	Lumen diameter [mm]
Arterial cannula for infants	C80G-040	80	200	5
Arterial cannula for small children	C80G-021	80	250	6
Arterial cannula for children, staged	C60G-004	60	280	9, 12; head 9
Arterial cannula for children, staged	C85G-004	85	280	9, 12; head 9
Arterial cannula	C60G-002	60	330	12
Arterial cannula	C85G-002	85	330	12

Tab. 12-4 Arterial cannulae

12.1.6 Overview: Which cannulae should be used for which pump?

Pump: connector ø [mm]	Which pump?	Cannula: lumen ø [mm] where cannula joins pump	Which inflow cannula?	Which outflow cannula? (arterial cannula)
6	P10P-001	5/ 6	C15V-040 (AT) C19V-020 (AT) C14A-040 (AP) C18A-020 (AP)	C80G-040 C80G-021
9	P25P-001x01 P30P-001x01	9	C22V-004 (AT;SC) C25V-004 (AT;SC) C22A-004 (AP;SC)	C60G-004 (SC) C85G-004 (SC)
12	P50P-001 P60P-001	12	C22V-004 (AT;SO) C25V-004 (AT;SO) C22V-002 (AT) C26V-002 (AT) C22A-004 (AP;SO) C27A-001 (AP)	C60G-004 (SO) C85G-004 (SO) C60G-002 C85G-002
Explanation:	AT atrial cannula AP apex cannula SO staged (stepped diameter) cannula, original diameter		SC staged (stepped diameter) cannula, diameter after cutting to size	

Tab. 12-5 Which cannulae for which pump?

12.1.7 System accessories

Designation	Article number
Accessory set for blood pumps with PU valves (membrane set, de-airing set and tube connecting set)	T00L-002
Driving tube, red; length: 200 cm	L20H-002
Driving tube, blue; length: 200 cm	L20H-003
Tank unit	1600422

Tab. 12-6 System Accessories

12.1.8 Driving unit

Designation	Article number
EXCOR® Stationary Driving Unit Ikus (115V/ 60Hz) - SW 3.41	D03I-111

Tab. 12-7 Driving unit

12.1.9 Special components


Designation	Article number
Connector set for cannulae ø 6 mm to ø 9 mm	A06-009
Connecting set for cannulae ø 9 mm to ø 12 mm	A09-012
Cannula tunnelling tip	attached to cannula

Tab. 12-8 Special components

12.1.10 Pump-cannula combinations

Cannulation		Blood pums				
∅ inflow cannula	∅ outflow cannula	10 ml	25 ml	30 ml	50 ml	60 ml
5 mm	5 mm	130 bpm				
6 mm	5 mm	130 bpm				
6 mm	6 mm	130 bpm	80 bpm	65 bpm		
9 mm	6 mm		100 bpm	90 bpm		
9 mm	9 mm		130 bpm	130 bpm	130 bpm	105 bpm
12 mm	9 mm				130 bpm	105 bpm
12 mm	12 mm				130 bpm	125 bpm

Tab. 12-9 Pump-cannula combinations



Pump-cannula combinations in which not every parameter combination is recommended (pump rate, % systole, systolic and diastolic pressure,) can lead to incomplete filling and emptying of the blood pump.

Rate value (bpm)

The value indicated is the upper threshold for pump rates. Values that are below the upper threshold are within the acceptable range. Values that are higher than the upper threshold are in a questionable range. The threshold values have been determined (in vitro) taking a mean arterial blood pressure of 120 mmHg as a basis.

Rate value (bpm)

Red marked values displayed on the laptop: These parameter combination (pump rate, % systole, systolic and diastolic pressure) for these pump-cannula combination can lead to incomplete filling and emptying of the blood pump(s). Observe the filling behavior of the blood pump(s)!

in biventricular mode

The lower value of both pump rates (corresponding to the pump sizes used) must also be considered. The higher of the 2 pump rates should be disregarded.

12.1.11 Blood pump combinations in biventricular mode

The following combinations are recommended:

- left pump 10 ml - right pump 10 ml (10 ml/ 10 ml)
- left pump 30 ml - right pump 25 ml (30 ml/ 25 ml)
- left pump 60 ml - right pump 50 ml (60 ml/ 50 ml)

Check whether a blood pump combination that is not recommended has been selected for the patient. The final decision on the combination of blood pumps and cannulae is to be reached by the implanting surgeon, in consultation with *Berlin Heart, Inc* Clinical Affairs.

12.1.12 Relative systolic duration

The relative systolic duration is adjustable in the range of 20% and 70%. The upper and lower threshold (20-30% and 60-70%) are marked in red on the laptop. For these values it cannot be guaranteed that the activated pressure parameters are achievable for each single case.

12.2 Sample copy: EXCOR Implantation log



Berlin Heart®

Implantation Record Form

EXCOR® VAD



This form applies **only** to USA and Canada



Please fill out the form (3 pages), and fax it to Berlin Heart, Inc. *immediately after implantation* (fax: 866.540.5026).

After replacing a blood pump, please fill out the "Pump Replacement" section (page 1), list the supplies used on page 2/3, and fax (3 pages) to Berlin Heart Inc. (fax: 866.540.5026)

Hospital	City/Country
----------	--------------

Patient data (for Berlin Heart registry)

Patient's initials	Sex m <input type="checkbox"/> / f <input type="checkbox"/>	Age	Body size [cm]	Weight [kg]
Patient-No. (BH Site No. followed by the patient No. ie: 004-103)	IABP pre-op n <input type="checkbox"/> y <input type="checkbox"/>		ECMO pre-op n <input type="checkbox"/> y <input type="checkbox"/> since _____ (Date)	
	On transplantation list n <input type="checkbox"/> y <input type="checkbox"/> since _____ (Date)			
Ischemic CMP <input type="checkbox"/>	Idiopathic CMP <input type="checkbox"/>	Acute Myocarditis <input type="checkbox"/>	Postcardiotomy <input type="checkbox"/>	
Acute Myocardial Infarction <input type="checkbox"/>	Congenital <input type="checkbox"/>	Other <input type="checkbox"/>		
PAP mean [mmHg]	CVP [mmHg]	MAP [mmHg]	LVEF % FS %	
CI [l/min/m ²]	NYHA	LVEDP [mmHg]	LVEDD [mm]	
Creatinine [mg/dl]	Total Bilirubin [mg/dl]	Platelet count [/μl]	Leukocytes [/μl]	

Implantation

Date	Surgeon			
Type BVAD <input type="checkbox"/> LVAD <input type="checkbox"/> RVAD <input type="checkbox"/>	Access medial <input type="checkbox"/> lateral <input type="checkbox"/>		Left-sided cannulation atrial <input type="checkbox"/> apical <input type="checkbox"/>	
LVAD	Pump type: PU valve <input type="checkbox"/>	Tilting-disk valve <input type="checkbox"/>		
	Pump size: 10 ml <input type="checkbox"/>	25 ml <input type="checkbox"/>	30 ml <input type="checkbox"/>	50 ml <input type="checkbox"/> 60 ml <input type="checkbox"/> 80 ml <input type="checkbox"/>
RVAD	Pump type: PU valve <input type="checkbox"/>	Tilting-disk valve <input type="checkbox"/>		
	Pump size: 10 ml <input type="checkbox"/>	25 ml <input type="checkbox"/>	30 ml <input type="checkbox"/>	50 ml <input type="checkbox"/> 60 ml <input type="checkbox"/> 80 ml <input type="checkbox"/>

Pump replacement

Left pump <input type="checkbox"/>	Reason for replacement _____
Date	Location of deposit inflow <input type="checkbox"/> outflow <input type="checkbox"/> pump chamber <input type="checkbox"/>
Right pump <input type="checkbox"/>	Reason for replacement _____
Date	Location of deposit inflow <input type="checkbox"/> outflow <input type="checkbox"/> pump chamber <input type="checkbox"/>

 **Berlin Heart®** **Implantation Record Form** **EXCOR® VAD**

 **This form applies *only* to USA and Canada** 

Please record the lot numbers of the used EXCOR® components used and the components to be kept as back-up and fax to Berlin Heart, Inc. *immediately* after implantation (fax: 866.540.5026).

Hospital/City _____ **Date of Implantation** _____

Patient ID (BH Site No. followed by the patient No. ie: 004-103) _____

Ikus-No. _____ **Ikus hours of operation** _____

Replacement Ikus
Ikus-No. _____ **Ikus hours of operation** _____

Item	Lot-No.		Article No.
	used	b/u	

EXCOR Blood Pumps with PU valves			
10 ml in/out Ø 6 mm			P10P-001
25 ml in/out Ø 9 mm			P25P-001x01
30 ml in/out Ø 9 mm			P30P-001x01
50 ml in/out Ø 12 mm			P50P-001
60 ml in/out Ø 12 mm			P60P-001
80 ml in/out Ø 12 mm			P80P-001***
EXCOR Blood Pumps with Tilting-disk valves			
50 ml in/out Ø 12 mm			P50M-001***
60 ml in/out Ø 12 mm			P60M-001***
80 ml in/out Ø 12 mm			P80M-001***
80 ml out/in Ø 12 mm (in/out exchanged)			P80M-005***
80 ml in/out Ø 16 mm			P80M-003***
80 ml out/in Ø 16 mm (in/out exchanged)			P80M-004***
EXCOR Apex Cannulas			
Ø 5 mm, L 22 cm (Apex cannula for infants)			C14A-040
Ø 6 mm, L 25 cm (Apex cannula for small children)			C18A-020
Ø 12/9 mm, L 27 cm (Apex pediatric cannula, staged)			C22A-004
Ø 12 mm, L 26,5 cm (Apex cannula, one-piece)			C27A-001
Ø 16 mm, L 33 cm (Apex cannula)			C41A-050***
EXCOR Atrial Cannulas			
Ø 5 mm, L 20 cm, head 15 mm (Atrial cannula for infants)			C15V-040
Ø 6 mm, L 25 cm, head 19 mm (Atrial cannula for small children)			C19V-020
Ø 12/9 mm, L 28 cm, head 22 mm (Atrial pediatric cannula, staged)			C22V-004
Ø 12/9 mm, L 28 cm, head 25 mm (Atrial pediatric cannula, staged)			C25V-004
Ø 12 mm, L 33 cm, head 22 mm (Atrial cannula)			C22V-002
Ø 12 mm, L 33 cm, head 26 mm (Atrial cannula)			C26V-002
Ø 12 mm, L 33 cm, head 30 mm (Atrial cannula)			C30V-002

*** Not available for general use in the US and Canada

+++ Not available for general use in the US



Berlin Heart®

Implantation Record Form

EXCOR® VAD



This form applies **only** to USA and Canada




Item	Lot-No.		Article No.
	used	b/u	
EXCOR Arterial Cannulas			
Ø 5 mm, L 20 cm (Arterial cannula for infants)			C80G-040
Ø 6 mm, L 25 cm (Arterial cannula for small children)			C80G-021
Ø 12/9 mm, L 26 cm (Graft-adapter pediatric cannula, staged)			C00P-004+++
Ø 12/9 mm, L 28 cm, 85° (Arterial pediatric cannula, staged)			C85G-004
Ø 12/9 mm, L 28 cm, 60° (Arterial pediatric cannula, staged)			C60G-004
Ø 12 mm, L 33 cm, 60° (Arterial cannula)			C60G-002
Ø 12 mm, L 33 cm, 85° (Arterial cannula)			C85G-002
Ø 12 mm, L 26 cm (Graft-adapter cannula)			C00P-001+++
Ø 16/12 mm, L 36 cm, 85° (Arterial cannula, staged)			C85G-050+++
Ø 16 mm, L 26 cm (Graft-adapter cannula)			C00P-050+++
Connecting Set for Cannulas			
Ø 6/9 mm			A06-009
Ø 9/12 mm			A09-012
Ø 12/16 mm			A12-016***
Accessories			
Accessory set Tilting-disk valves			T00L-001***
Accessory set PU-valves			T00L-002
Driving tube, red Ø 6/8 mm, L 2 m			L20H-002
Driving tube, blue Ø 6/8 mm, L 2 m			L20H-003
Tank unit			1600422

*** Not available for general use in the US or Canada
 +++ Not available for general use in the US

Date	Signature
-------------	------------------

12.3 Sample copy: EXCOR pump log

12.3.1 Explanations on the pump log



- 1 transition inflow cannula - inflow connector
- 2 only on pumps with PU valves: inflow stub in front of inflow valve
- 3 inflow valve
- 4 inflow stub behind inflow valve
- 5 area between inflow and outflow stubs
- 6 remaining area of blood chamber
- 7 transition blood chamber - membrane (directly above the reinforcement ring)
- 8 outflow stub in front of outflow valve
- 9 outflow valve
- 10 only on pumps with PU valves: outflow stub behind outflow valve
- 11 transition outflow connector - outflow cannula


Fig. 12-2 Numbering of the checkpoints

ADVICE

To briefly describe the findings, use the following letter codes:

- p - small punctual deposit
- a - small area of deposit
- f - small strand
- t - small thrombus
- ~above the respective letter indicates floating deposits
- P - large punctual deposit
- A - large area of deposit
- F - large strand
- T - large thrombus

Example: Plotting of the deposits



- 1 small laminar deposit
- 2 small suture on the inflow valve
- 3 small specks
- 4 larger laminar deposit


Fig. 12-3 Plotting of the deposits

Example: Notation with letter code

Datum date	Zeit time	Name Sign.	Linke Pumpe/ Left pump		ml: 50ml No.: 0815										
					1	2	3	4	5	6	7	8	9	10	11
01.01.	8:00	z.B.			a		F			p		A			

Fig. 12-4 Notation for letter code

12.3.2 Sample copy: EXCOR pump log



- 1 transition inflow cannula - inflow connector
- 2 only on pumps with PU valves: inflow stub in front of inflow valve
- 3 inflow valve
- 4 inflow stub behind inflow valve
- 5 area between inflow and outflow stubs
- 6 remaining area of blood chamber
- 7 transition blood chamber - membrane (directly above the reinforcement ring)
- 8 outflow stub in front of outflow valve
- 9 outflow valve
- 10 only on pumps with PU valves: outflow stub behind outflow valve
- 11 transition outflow connector - outflow cannula

Fig. 12-5 EXCOR blood pump with checkpoint numbers


Patient:

Lot No.: Pump left

Lot No.: Pump right

Datum date			Zeit time			Name Sign.			Linke Pumpe/ Left pump				Rechte Pumpe/ Right pump																					
									ml:				ml:																					
									No.:				No.:																					
									1	2	3	4	5	6	7	8	9	10	11	1	2	3	4	5	6	8	7	9	10	11				

12.4 Sample copy: EXCOR parameter log

 Berlin Heart, Inc.		Berlin Heart EXCOR [®] Pump Performance Flow Sheet (FM-CA-031, Rev. 1)															
Ikus-No _____-	Date	Time	MODE		RATE		Drive Pressures (mmHg)				% Systole		Membrane Movement				
			L	R	L	R	Left systole	Left diastole	Right systole	Right diastole	L %	R %	LEFT eject	LEFT fill	RIGHT eject	RIGHT fill	
	Sample:	9:15 AM	S		60	60		210	-30	130	-30	40	40	++	+-	++	++

Membrane movement:
 ++ complete or almost complete (max. 1 small wrinkle in membrane) movement of membrane to end position
 +- incomplete ejection or filling movement (membrane has more than 1 small wrinkle at end position of ejection or filling movement)
 .. bad ejection or filling (membrane moves only partially or not at all)

Index

A

accessories	
description - blood pump, cannulae, ~	45
EXCOR components and ~	49
EXCOR Pediatric ~	47
system ~	113
accessory set	
for blood pumps with PU valves	49, 113
for commissioning and operation	47
adverse events	20
ambient temperature	18
anticoagulation	
status	81, 89
therapy	83

B

biventricular	
~ operation, adjusting the operating mode	67
anastomosis	63
assist device, shown in place	45
dysfunction	21
explantation after ~ support	108
overview, assignment	64
blood pump	
combinations in biventricular mode	114
connecting the ~ to the cannulae	63
consists of ...	21
de-airing	51
EXCOR Pediatric ~	46
explantation	95
operation period	11, 22
preparing a replacement ~	90
product range	111
regular checks of ~	79
replacing	89
rinse and fill the ~	52
BTT, bridge to transplant	109

C

cannula tunneling tip	54
available sizes of ~	55
use of ~	55
cannulae	
apex ~	56, 112
arterial ~	112
atrial ~	59, 112
connecting the blood pumps to the ~	63
description	45
EXCOR Pediatric ~	46
outflow ~	60
regular checks of the ~	79
removing the VAD ~	109

shortening	62
types of ~	46
Which ~ should be used for which pump?	113
cardiac transplantation	21
Carmeda® BioActive Surface	12, 46
cautionary measures	
adjusting the parameter values	80
EMC	18
filling behavior of the blood pump(s)	80
visual inspection - deposits	81
checks	
~using the monitor program	81
regular ~ of blood pumps and cannulae	79
Cleaning	16
components	
and accessories	49
permanently active ~	45
preparing the ~ and materials required	49
contraindications	22

D

de-airing	
~ hammer	12, 49
~ needle, figure	51
~ nipple, figure	51
~ set	47, 49, 51, 113
~ step	65
inserting the ~ needle	51
removing the ~ needle	69
deposits	
~ in the pump	89
floating	89, 118
large thrombus	118
large, laminar	118
large, speck-like	118
plotting of the ~	118
small thrombus	118
small, laminar	118
small, speck-like	118
visual inspection	81
dressing	
applying a new ~	78
removing the old ~	76
driving pressure	
possible range	67
systolic	21
driving tube	45
as an active component	45
connector	46
connector, figure	46, 51
function	45

E

echocardiogram, intraoperative trans-esophageal	73
electro-stimulation therapy	19
EMC	

special precautions required	18
events, adverse	20
examination, medical ~ of patient	80
EXCOR Pediatric	49
parameter log, sample copy	120
pump log, sample copy	119
EXCOR Pediatric HDE Study	90
explantation	95
~ after biventricular support	108
~ after univentricular support	108

I

Ikus	
ambient conditions	19
designation	9
device description	21
software version	3
troubleshooting and correcting faults	89
implantation	
anesthesia	73
preparations in the operating room	49
surgical procedure	53
infection	75
risk of ~	75
to prevent ~	13
inspection	
filling behavior of blood pump	80
frequency of ~	79
of pump areas which come in contact with blood	81
visual ~ - deposits	81
instruction and care of the patient	20
interaction with other procedures and therapies	19

K

key switch (main switch)	
switching the Ikus off	108, 109

L

laptop	
software	3
version	3
LVAD	
~ anastomosis	63
~ anesthesia	73
~ assignment	64

M

magnetic resonance imaging	18
manual pump	
driving a blood pump with the ~	92
figure	94
membrane set	47, 49, 113
mobilization	80
monitoring	

~ procedure	73
-------------------	----

N

nuclear	
~ diagnostics	19
~ therapy	19

O

obturator	51
operating manager, obligations of ~	20
operating mode	64

P

packaging	
and sterilization	12
sterile ~	51
parameter	
adjusting the ~ values	80
adjustment, possible range	67
checking the ~s when the pump is started and adjusting them	66
EXCOR Pediatric parameter log	120
platelet concentrates	19, 73
product range and possible combinations	111
power supply	
failure	92
precautions, general	19
pump log	
explanations on the EXCOR Pediatric ~	118
sample copy EXCOR Pediatric ~	119
pump size	
relationship between bodyweight and ~	111
pump, see blood pump	21

R

radiotherapy	18, 19
replacing blood pump	89

S

safety risks	20
sample copy	
Blood pump replacement	120
EXCOR Implantation log	115
EXCOR Pediatric parameter log	120
EXCOR Pediatric pump log	119
single-step mode	
figure	65
removing air from the blood pump(s) in ~	65
sterilization and packaging	12
storage	
and durability	20
suction pressure, diastolic	21
possible range	67

T

technical specifications	114
toggle switch (power switch)	
switch off the Ikus	108, 109
transplantation	
orthotopic heart ~	110
treatment and wound care	75
Trendelenburg position	63, 65, 90, 91
trocar	
as part of the de-airing set	47, 49
figure	51
to de-air the blood pump	51

U

ultrasonic treatment	19
univentricular	
overview, assignment	64
which driving tube?	49
UVAD, univentricular	
explantation after ~ support	108

V

valves	
three-leaflet polyurethane ~	21
visual inspection	
deposits	81
filling and ejecting of the blood pump	16
filling behavior of blood pump	80
regular ~ of blood pumps and cannulae	79

W

weaning	110
wound care	75
and treatment	75

X

X-rays	19
--------------	----

This page was left blank intentionally