

SynCardia (TAH-t)




The Device

The TAH-t replaces the tricuspid valve, right ventricle, pulmonary valve, mitral valve, left ventricle, and aortic valve. Only a small rim of patient atrial tissue remains.

Patient Device Size Selection

- $>1.7M^2$ BSA = 70 cc
- $\leq 1.7M^2$ BSA = 50 cc




Implant Procedure




Scan QR for a medical animation.


Device Components



HOSPITAL CART



DRIVER SYSTEM



*Battery Life: 1 hour

C2 Driver: Typical Settings

Typical Device Settings

- Rate: 125 ± 15 bpm
- Right Drive Pressure: 80–100 mmHg
- Left Drive Pressure: 180–210 mmHg
- Right Vacuum Pressure: 0–(-10) mmHg
- Left Vacuum Pressure: 0–(-13) mmHg
- Left vacuum typically greater than right vacuum
- % Systole: 50 ± 5 %



Anticoagulation



- Start **Heparin** within 72 hours, with goal of PTT 50–70 (should be determined by care team)
- **Aspirin** POD 3–5
- Start **Coumadin** when PO intake is tolerated
 - Maintain INR between 2–3
 - Monitor intake of Vitamin K


Blood Pressure

- Pump is afterload-sensitive, maintain blood pressure goal below 135 mmHg
- Patient will have palpable pulse at the set pump rate


Fluid Status

- Pump is preload-dependent
- Prevent overflow with close attention to fluid status

Pressure



Flow



Eject: Determined by systolic pressure and % systole. Look for full eject waveform.

Fill: Determined by rate, vacuum and % systole. Look for partial fill.

Goal Fill: 70 cc TAH = 50–60 mL 50 cc TAH = 30–40 mL

Avoid Emergencies

- There is **no ECG lead**
- **No chest compressions**
- High afterload, such as code dose epi, could lead to the pump unable to eject.
- **No upper lines** should be placed. If a wire gets into the valve the TAH could stop.
- If concern for **obstruction** to flow consider liver U/S to look at IVC size and for obstruction.
- Biggest concern post-op should be for **bleeding**

Emergency Care

Bleeding

- **Why:** This occurs more frequently in the TAH due to the multiple anastomosis sites
- **Watch for:** Changes in chest tube drainage and signs of low output. This may include lower CO on the device.
- **Treatment:** Alter the anticoagulation strategy of consider surgical intervention.

Tamponade (or Compression)

- **Why:** Fluid accumulation may cause compression around the native atria and the SVC and IVC.
- **Watch for:** Decreased fill volumes and decline in overall CO, from baseline. Watch for unique waveform patterns to determine location of compression.
- **Treatment:** Urgent surgical exploration is needed if CO is inadequate.

Device Malfunction

- **Why:** As with all devices there could be a mechanical malfunction.
- **Watch for:** Driver to read inaccurate numbers or to stop working.
- **Treatment:** If the patient is stable switch drivers. If the driver stops use the handpump while retrieving back up driver.

Device Overfill

- **Why:** The pump should only reach partial fill or the fluid can back into the lungs and cause pulmonary edema.
- **Watch for:** Changes in lung compliance and high fills on the device. Also watch for the full fill waveform. The 50 cc pump does not alarm for overfilling.
- **Treatment:** Increase the rate on device, decrease the vacuums and give diuretics. Support respiratory status.