

# CentriMag™/PediMag™

## Unresponsive Patient

Provide adequate ventilation & Assess perfusion

## Adequate Perfusion?

MAP goal \_\_\_\_\_

YES ↙

NO ↘

### Treat non-VAD causes for AMS

- Stroke
- Hypovolemia
- Bleeding
- Hypoxia/hypoventilation
- Arrhythmias
- Sepsis

### Is the VAD working?

- Assess flows (volume/afterload)
- Assess pump (motor, magnet), console, and circuit
- Verify power is connected

<1 Minute

YES  
↓

NO  
↓

**Start  
CPR**

- Pump/console change

<1 Minute

If no improvement in MAP, **start CPR**

## VAD

☐ CentriMag™

☐ PediMag™

## CANNULA TYPE

☐ Temporary  
(If yes, CPR only at discretion of code leader)

☐ Berlin (long term)

## AORTIC VALVE ANATOMY

☐ Closed (Oversewn)

☐ Open